

# Ned User Manual



## Ned

---

### Robot

Ned is a 6-axis collaborative robot arm designed for Education and Research.



Ned is designed to reproduce all the movements required in the most advanced uses in industry 4.0, with a **precision of 0.5mm** and a **repeatability of 0.5mm**.

Ned's **aluminum structure** makes it exemplary robust, allowing it to accomplish with fluidity the movements required for your robotics projects.

This **cobot** takes advantage of the capacities of the **Raspberry Pi 4**, with a 64-bit ARM V8 high performance processor, 2GB of RAM and an improved connectivity.

Ned is a **collaborative robot based on Ubuntu 18.04 and ROS** (Robot Operating System) Melodic, an open-source solution created for robotics. Through ROS, **Ned has multiple libraries** allowing you to conceive many programs, from the most simple to the most complex ones, responding then, in a flexible way, to your needs.

## Ned tools technical specifications

Mounting a tool on Ned is made easy with our brand new **EasyConnect system**. Simply plug your tool, connect its cable, and it is ready to use.

Ned's package includes a **Custom Gripper**. Its standard jaws can manipulate small objects, and you can 3D print your own custom jaws or buy our [set of jaws](#) ([index.html#overview-standard-jaws](#)) that contains the precision jaws, the flat ones and the XL ones.




Ned can also be used with:

- The [Large Gripper](#) (index.html#large-gripper) which has the ability to grab larger objects while keeping the ability to close entirely;
- The [Adaptive Gripper](#) (index.html#adaptive-gripper) which allow Ned to grab non-standard objects with complex shapes (eg. an egg);
- The [Vacuum Pump](#) (index.html#vacuum-pump) to grab objects with plane and non-porous surfaces;
- The [Electromagnet](#) (index.html#electromagnet) that is useful to manipulate metallic objects, from one to many (eg. screws, bolts...).

**Ned's Ecosystem is designed to let you reproduce advanced use cases of industry 4.0.**

## Vision Set



**The Vision Set gives Ned the ability to see its environment and detect the objects to interact with based on their characteristics** (shape and/or color). This set includes :


- A camera that is mountable on Ned's wrist to give it a dynamic field of view depending on Ned's position;
- A workspace which is designed around a specific repositionnable material;
- Supports for the workspace and the robot, to build a rigid environment for the vision;
- 6 objects (3 rounds, 3 squares) of different colors to use our built-in recognition system based on colors and shapes. The squares are designed to be able to act as containers for the circles. Easily create an industry 4.0 pick & pack application with no particular knowledge of programming;
- Vision jaws for the Custom Gripper, that are specifically designed for the provided objects;
- A calibration tip to set your workspace up.

## Conveyor Belt

**The Conveyor Belt is the key add-on to prototype industry 4.0 production lines.**



The **Standard version** includes the Conveyor Belt, its control box for autonomous use, and an IR sensor to detect the presence of an object on it.



The **Education version** adds the following components:

- A mountable workspace to use vision functions directly on the Conveyor Belt;
- A slope that can help to create complex or multi-robot processes;
- An end-stop to stop the objects when they reach the end of the Conveyor Belt;
- The 6 objects of the Vision Set to extend the quantity of objects you can use on your production line;
- The Vision jaws to manipulate them.

**With Ned's ecosystem, you have the optimal solution to discover, learn and test industry 4.0 processes.**

## Overview of this manual

This manual contains instructions for:

- Safety Information about Ned;
- Mechanical and electrical installation of the robot;
- Software of the robot;
- Maintenance and troubleshooting.

This manual should be read before:

- Installation and electrical connections;
- First usage of the robot;
- Maintenance and repair.

## Safety notice

### Warning

You must read, understand, and follow all safety information in this manual.

It is important that the safety instruction on this manual and the electrical and the mechanical instructions are followed. The robot specifications should be respected to avoid all damage to the machine or to the user.

## Additional content

You can find additional content by clicking on the link[here](https://docs.niryo.com/) (<https://docs.niryo.com/>).

## Installation and commissioning

This section specifies the procedures to follow during unpacking and transportation.

### What the box contains

When you order Ned, you receive a box with the following items:

- 1x Ned Robot
- 1x Custom Gripper
- 1x Power cable
- 1x Universal travel plug adapter
- 1x SD card reader
- 1x Allen Key set

### Unpacking


1. Inspect Ned's packaging to make sure nothing is damaged;
2. Remove the packaging and check for any transport damage or loss;
3. Put Ned in a dedicated operating environment;
4. Make sure Ned is stable.

### Operating conditions

Parameters	Value
Minimum Ambient Temperature	+5°C
Maximum Ambient Temperature	+45°C

## Shipping and transport position

This figure shows Ned in its shipping position, which is a recommended transport position to avoid Ned's damages.



## Safety instructions

### Safety signals and symbols

This section specifies all dangers that can arise when doing the work described in the User Manual. Each danger consists of:

- A caption specifying the danger level (DANGER, WARNING, or CAUTION) and the type of danger;
- A brief description of the danger;
- Instruction about how to eliminate danger.

Designation	Description
DANGER	Warns that an accident will occur if the instructions are not followed, resulting in an injury.
WARNING	Warns that an accident will occur if the instructions are not followed, resulting in a possible injury or product damage.
CAUTION	Warns of possible product damage.
NOTE	Describes an important fact or condition.
TIP	Describes additional information or an easy way to use the function.

Icon	Designation	Description
	ELECTROSTATIC DISCHARGE	Warns of a possible electrostatic discharge that may damage the robot.
	ELECTRICAL SHOCK	Warns of an electrical hazard.

## Safety Overview

This section includes information on general safety risks to be considered when performing installation and service work on Ned.

These safety instructions must be read and followed by any person who deals with Ned.

## Spare parts and special equipment

The installation and the use of external components is possible (such as Arduino, external sensors, devices). But the user must respect the electrical specifications and connections mentioned in the User Manual.

Niryo is not liable for damages or injuries caused by unauthorized modifications to the robot system or the misuse of the robot.

## Complete robot

### Danger

1. Before installing, using, or programming Ned, read the product specifications and Manual;
2. Make sure to install Ned and all electrical equipment according to the specifications and the warnings from chapters [Installation and commissioning](#) ([index.html#installation-and-commissioning](#)) and [Technical specifications](#) ([index.html#technical-specifications](#));
3. Children are not allowed to operate Ned without adult supervision to prevent any possible injury or improper misuse of Ned;
4. Make sure that changing robot software will not cause hazards, damages to the Ned's system or injuries.

### Warning

1. Motors and gearboxes are hot after running Ned for a long time. Avoid touching especially with a higher environment temperature;
2. Removing parts may result in the collapse of Ned. Take the necessary measures to ensure that Ned does not collapse as parts are removed;
3. Make sure to remove the power cable during assembly, wiring, or repair. Doing so will help to prevent equipment damage or accidental short-circuit;
4. The cable packages may be sensitive to mechanical damage;
5. Before using Ned, make sure to firmly fix it on a stable surface to avoid collapse or the fall of Ned;
6. Do not use Ned outside;
7. Do not put Ned in a humid environment or near water;
8. Do not install or operate Ned in dangerous environments (e.g., in the presence a strong magnetic field, dangerous gases, fire or flammables) to avoid dangers which may occur due to external conditions during operation;
9. Respect the technical specifications to avoid damaging the motors. Stepper motors and servo motors may be damaged if excessive force is used;
10. Make sure external change on the motors parameters or the robot software does not cause damage to Ned.
11. Do not open Ned base when the power cable is connected. Do not touch the electronic board, they may be hot.

Niryo is not liable for damages or injuries caused by unauthorized modifications to the robot system or the misuse of the robot.

**Caution**

1. Do not turn off Ned during a sequence or a motion unless absolutely necessary, this may reduce Ned operating age;
2. Try to install Ned arm on a stable surface with enough space to avoid any shock or vibration.

## Voltage related risks, robot

A danger of high voltage is associated with the robot in:

- The power supply of the motors (11.1V, 8V);
- Ned's grippers are powered by 8V. A risk of short circuit and sparks exists. Be sure to make the installations with the power disconnected.

## Safety actions

- Fire extinguishing

**Note**

Use a CARBON DIOXIDE (CO<sub>2</sub>) extinguisher in the event of a fire in Ned's system.

- Manually stopping or overriding the arm: if needed, Ned arm can be stopped manually. This is possible since Ned arm is light, the arm force is limited, and the drivetrain power is limited.

To prevent unnecessary damage and wear of Ned arm, it is recommended to use the normal stopping functions of Ned (by the software or the top button).

**Note**

Ned is not equipped with brakes since its weight and design do not require a holding brake.

## Technical specifications

This section describes the technical specifications of Ned.



Parameters	Value
Weight	6,5 kg
Payload	300 g
Reach	440 mm
Degree of freedom	6 rotating joints
Joints range	$-170^\circ \leq \text{Joint 1} \leq 170^\circ$ $-120^\circ \leq \text{Joint 2} \leq 35^\circ$ $-77^\circ \leq \text{Joint 3} \leq 90^\circ$ $-120^\circ \leq \text{Joint 4} \leq 120^\circ$ $-100^\circ \leq \text{Joint 5} \leq 55^\circ$ $-145^\circ \leq \text{Joint 6} \leq 145^\circ$
Joints range	$-2,97 \text{ rad} \leq \text{Joint 1} \leq 2,97 \text{ rad}$ $-2,09 \text{ rad} \leq \text{Joint 2} \leq 0,61 \text{ rad}$ $-1,34 \text{ rad} \leq \text{Joint 3} \leq 1,57 \text{ rad}$ $-2,09 \text{ rad} \leq \text{Joint 4} \leq 2,09 \text{ rad}$ $-1,75 \text{ rad} \leq \text{Joint 5} \leq 0,96 \text{ rad}$ $-2,53 \text{ rad} \leq \text{Joint 6} \leq -2,53 \text{ rad}$
Joints speed limit	Joint 1 $\leq 150^\circ/\text{s}$ Joint 2 $\leq 115^\circ/\text{s}$ Joint 3 $\leq 140^\circ/\text{s}$ Joint 4 $\leq 180^\circ/\text{s}$ Joint 5 $\leq 180^\circ/\text{s}$ Joint 6 $\leq 180^\circ/\text{s}$
Joints speed limit	Joint 1 $\leq 2.6 \text{ rad/s}$ Joint 2 $\leq 2.0 \text{ rad/s}$ Joint 3 $\leq 2.5 \text{ rad/s}$ Joint 4 $\leq 3.14 \text{ rad/s}$ Joint 5 $\leq 3.14 \text{ rad/s}$ Joint 6 $\leq 3.14 \text{ rad/s}$
TCP max speed	1144 mm/s
Repeatability	0.5 mm

Parameters	Value
Footprint	200x200 mm
Mounting	Table
I/O power supply	5V for digital GPIO 12V for Switches output
I/O ports	Digital inputs/outputs 5V: 6 Digital outputs 12V: 2
I/O interface	2 x USB2.0 2 x USB3.0 1 x ETHERNET GIGABIT
Communication	Modbus TCP (master) TCP/IP CAN (slave, compatible with NiryoStepper motors)
Materials	Aluminum Plastic PLA
Temperature	5 - 45°C
Power supply	DC 11.1V - 6A
Programming Environment	Niryo Studio PyNiryo C++ Python ROS

## Tools technical specifications

### Tools Overview


Ned is compatible with 5 different sorts of tools:

- 3 x Grippers
- 1 x Electromagnet
- 1 x Vacuum Pump

### Specifications

#### Custom Gripper

#### Overview - Standard Jaws




(\_images/GRIPPER\_STANDARD\_SCHEMA\_1.jpg)




(\_images/GRIPPER\_MORS\_STANDARD.png)

### Overview - XL Jaws



(\_images/GRIPPER\_XL\_SCHEMA\_1.jpg)




(\_images/GRIPPER\_XL\_SCHEMA\_2.jpg)




(\_images/GRIPPER\_MORS\_XL.png)

## Overview - Flat Jaws



(\_images/GRIPPER\_FLAT\_SCHEMA\_1.jpg)




(\_images/GRIPPER\_FLAT\_SCHEMA\_2.jpg)



(\_images/GRIPPER\_MORS\_FLAT.png)

## Overview - Precision Jaws




(\_images/GRIPPER\_PRECISION\_SCHEMA\_1.jpg)




(\_images/GRIPPER\_MORS\_PRECISION.png)

## Overview - Vision Jaws



(\_images/GRIpper\_VISION\_SCHEMA\_1.jpg)



(\_images/GRIpper\_VISION\_SCHEMA\_2.jpg)




(\_images/GRIpper\_MORS\_VISION.png)


### Custom Gripper's technical specifications

Parameter	Description
Max operating width	Standard Jaws: 23.82 mm XL Jaws: 73.92 mm Flat Jaws: 35.92 mm Precision Jaws: 21.57 mm Vision Jaws: 48.92 mm
Picking distance from end effector base	Standard Jaws: 85 mm XL Jaws: 85 mm Flat Jaws: 90 mm Precision Jaws: 85 mm Vision Jaws: 85 mm
Motor	XL320 Servo Motor
Weight	Standard Jaws: 109 g XL Jaws: 107 g Flat Jaws: 123 g Precision Jaws: 101 g Vision Jaws: 107 g
Power Supply	7.6 V
Operating temperature	5-45°C

### Large Gripper



(\_images/GRIPPER\_LARGE\_SCHEMA\_1.jpg)




(\_images/GRIPPER\_LARGE\_SCHEMA\_2.jpg)




(\_images/GRIPPER\_LARGE.png)

Parameter	Description
Max operating width	60 mm
Picking distance from end effector base	80 mm
Motor	XL320 Servo Motor
Weight	90 g
Power Supply	7.6 V
Operating temperature	5-45°C

## Adaptive Gripper



(\_images/GRIPPER\_ADAPTATIVE\_SCHEMA\_1.jpg)




(\_images/GRIPPER\_ADAPTATIVE\_SCHEMA\_2.jpg)



(\_images/GRIPPER\_ADAPTATIF.png)

Parameter	Description
Max operating width	90 mm
Picking distance from end effector base	85 mm
Motor	XL320 Servo Motor
Weight	110 g
Power Supply	7.6 V
Operating temperature	5-45°C

## Electromagnet




(\_images/GRIpper\_ELECTROMAGNET\_SCHEMA\_1.jpg)




(\_images/AlMANT.png)

Parameter	Description
Max operating width	20 mm
Picking surface	20 mm
Picking distance from end effector base	22.5 mm
Control interface	Digital output
Weight	42 g
Power Supply	5 V
Operating temperature	5-45°C

## Vacuum Pump



(\_images/GRIpper\_VACUUM\_PUMP\_SCHEMA\_1.jpg)



(\_images/GRIpper\_VACUUM\_PUMP\_SCHEMA\_2.jpg)



(\_images/POMPE.png)


Parameter	Description
Vacuum Pump control	Integrated XL320 servo motor
Weight	10 g
Total weight	170 g
Payload	300 g
Power Supply	7.6 V
Operating temperature	5-45°C
Vacuum Tube dimension	20 mm
Picking distance from end effector base	47 mm

## Mechanical interface


### Mechanical Overview

This chapter introduces the mechanical and the electrical interfaces of Ned.


Ned is a 6-axis collaborative robotic arm. It consists mainly of six robot joints of aluminum with plastic covers. Ned consists of 7 parts :




### Three-view of the robot



## Robot workspace





## Tool mounting

Each tool has the same mechanical connector interface. This option allows a quick-change option of the tool and good stability. See the picture below.



## Vision Set mounting

To mount the Vision Set, all you need to do is attach the camera to Ned's wrist, set up and fix your workspace to Ned's structure to have a stable environment and put the vision jaws on the Custom Gripper.



## Robot calibration

Ned is delivered without a mechanical calibration. Ned can not execute any motion properly without auto-calibration because of a mismatch between the origin of each motor and its corresponding origin stored in the controller.

### Note

Make sure to process the calibration in order to align the origins of every joint.


### Hint

Check [Niryo Studio Manual](https://docs.niryo.com/product/niryo-studio/index.html) (<https://docs.niryo.com/product/niryo-studio/index.html>) to get more information about Ned's calibration.

## Electrical interface

### Electrical interface Overview

Electrical architecture:



## Electrical warnings and cautions

### **⚠ Danger**

- Make sure that all the equipment and the wires are kept dry. If water enters the equipment, disconnect the power;
- Carefully follow the safety instructions of the next section to avoid damaging your robot.


### **⚠ Warning**

Do not unplug or power off the robot during use, you may damage the Raspberry Pi 4 inside the robot, or the micro SD card inside the Raspberry Pi 4

## Panel connectors

### Panel connectors Overview





1	2 x USB port 2.0
2	2 x USB port 3.0
3	Ethernet port
4	Panel LED
5	2 x NiroStepper connector
6	XL320 connector
7	XL430 connector (not used)
8	2 x SW : 12 V digital output
9	2 x GPIO panel
10	Top button

## Power supply

Ned is powered with AC/DC adapter with 11.1V output

Parameters	value
Frequency	50/60 Hz
Input	AC 100-240 V , 1.5 A
Output	11.1 V – 6 A

### ⚠ Warning

- When Ned is turned on, do not disconnect the power supply of Ned. You may damage the Raspberry Pi 4 inside Ned, or the micro-SD card.
- If you modify the original power supply of the robot, you should respect the same specifications.
- Before plugging the power connector and powering on Ned, first make sure that:
  - Ned is on a flat and stable surface;
  - The power switch is turned off (position 0);
  - Ned has enough space to move without hitting something or someone.

## Power Switch

The power switch turns Ned on and off. It enables the power supply.

Always switch it off after shutting down Ned and turn it on to power on Ned.

### Warning

Before plugging or unplugging the power supply, make sure to turn off the switch.

Do not keep the switch on after shutting down Ned, otherwise , the power will be present on the motors and on the shield



## Digital inputs/outputs

- Ned has 6 inputs/outputs that can be configured either as digital input or digital output with Niryo Studio (see How to program the digital I/O on the [Software Manual](#) (<https://docs.niryo.com/product/niryo-studio/index.html>)).

Terminal	Parameters	Min	Max	Unit
<b>Digital Input</b>				
[1A-2A-1B-2B-1C-2B]	Voltage	0	5	V
<b>Digital output</b>				
[1A-2A-1B-2B-1C-2B]	Voltage	0	5	V



- Digital input connection:

This illustration below shows how to connect a button to a digital input.



- Digital output connection:



This example shows how to connect a load in order to be controlled from a digital output (see below).



## Switch connections

Terminal	Parameters	Min	Max	Unit
<b>Digital Output</b>				
SW1, SW2	Voltage	0	11.6	V
	Current	0	0.5	A

Below is an example of how to connect an external load to the switch output.



## Ethernet

An Ethernet connection is provided at the panel connector of Ned.

The Ethernet interface can be used for the following:

- Modbus TCP/IP;
- PyNiryo;
- Connect Ned to your local area network.



## USB Port

The USB can be used to connect the keyboard, mouse, or Vision Set.

Find more information about Vision Set on the [Vision Set Manual](https://docs.niryo.com/product/vision-set/index.html) (<https://docs.niryo.com/product/vision-set/index.html>).

## NiryoStepper connectors

NiryoStepper connectors interface can be used to drive a NiryoStepper motor. It is based on CAN Bus communication protocol.




This interface can be used to drive Niryo's [Conveyor Belt](https://docs.niryo.com/product/conveyor-belt/index.html) (<https://docs.niryo.com/product/conveyor-belt/index.html>).

## DXL connection

### XL320 Connector


This connector can be used to connect Ned's Vacuum Pump tool.



GND    7.2V    DXL-Data1

## XL430 Connector

Not used.



GND    7.2V    DXL-Data2

## Panel LED

Ned's LED panel has several colors which represent different modes and error status. Refer to the following table:

Color / Blinking	Description	Troubleshooting
Solid blue	Hotspot mode, the robot has started and creates its own Wi-Fi network.	N/A
Solid green	Connected mode, the robot is fully started and connected to a Wi-Fi network.	N/A
Alternating between Green/Red	Motor error on connected mode	Check "HARDWARE STATUS" and logs to troubleshoot the error.
Alternating between Blue/Red	Motor error on hotspot mode.	Check "HARDWARE STATUS" and logs to troubleshoot the error.
Solid purple	Robot controller is shutting down	N/A
Alternating white – green	The robot is paused on connected mode	N/A
Alternating white – blue	The robot is paused on hotspot mode	N/A
Solid red	On start-up: the robot is booting; On power off: shutdown completed	N/A

## Top button

The top button of Ned has multiple functionalities depending on how long you press the button and the number of presses.

### Note

**Note**

All functionalities are available for use only after the robot has been fully started (LED is green or blue).

Press count / Delay	Function
<b>Program is executing</b>	
1 press - 0.1 and 3s	Pause the robot program with Learning Mode off.
2 press - 0.1 and 3s	Pause the robot program with Learning Mode on.
<b>Program is paused</b>	
1 press - 0.1 and 3s	Continue the program.
2 press - 0.1 and 3s	Stop the program.
<b>Program is set on autorun</b>	
1 press - 0.1 and 3s	Start an autorun program.
<b>Programming the robot with blockly</b>	
1 press - 0.1 and 3s	Add / save a freedrive joint position.

**Note**

When you pause the robot's program with the top button, you should either press the button to continue or press twice to stop it.

**Hint**

You can test your robot hardware for the first time without using any software tool, you can follow the steps below:

- Plug the power supply;
- Turn on the power switch button;
- Wait until the LED panel turns blue (if it is the first time you use Ned this should take 2 minutes);
- Press the top button: Ned will start by calibrating and then it will execute a pre-saved sequence.

## Tool wiring

### Custom Gripper, Large Gripper and Adaptive Gripper:

There is a cable on the tool end of Ned with 3 pins. This connector can drive Custom Gripper, Large Gripper and Adaptive Gripper.

**Note**

Make sure to respect color order when you are wiring your gripper.

**Warning**


Not respecting the wiring instructions may damage your gripper or Ned.

## Vacuum pump:

Plug the suction cup to Ned's hand and connect the motor wire to the [XL320 Connector](#) placed at the back of Ned.

### **Electromagnet:**

Plug the electromagnet to Ned's hand and connect the cable on the GPIO1 or GPIO 2 connector on the panel connector.



#### **Note**

Information about programming and using the tool are available on the [Software Manual](#) (<https://docs.niryo.com/product/niryo-studio/index.html>)

## **Advanced Programming**

There are many other ways to develop on Ned.

For more information, please refer to the following links:

- [Niryo Studio](https://docs.niryo.com/product/niryo-studio/index.html) (<https://docs.niryo.com/product/niryo-studio/index.html>)
- [Python](https://docs.niryo.com/dev/pyniryo/index.html) (<https://docs.niryo.com/dev/pyniryo/index.html>)
- [ROS](https://docs.niryo.com/dev/ros/index.html) (<https://docs.niryo.com/dev/ros/index.html>)
- [Modbus](https://docs.niryo.com/dev/modbus/index.html) (<https://docs.niryo.com/dev/modbus/index.html>)

[Suggest a modification](#)[Download as PDF](#)

# Niryo Studio User Manual



Niryo Studio is a graphical HMI. It allows a fast and direct control of Ned with an external computer.

Its purpose is to provide users with a complete and simple interface for Ned motion, programming environments and current status of Ned.

Through Niryo Studio application, users can:

- Manage and set Ned parameters
- Control Ned motion
- Program Ned using Blockly
- Execute a program
- Manage, set, control and program Ned tools
- Manage, set, control and program Ned addons

Users must read this manual in order to fully understand all Ned and Niryo Studio functionalities.

## **Note**

If you encounter software problems or if you want to access latest updates for your Ned, please give an internet access to your robot (Ethernet or Wi-Fi) to be able to update it through NiryoStudio. Instructions [here](#) (index.html#ned-software-update)

## **Changelog**

### **September release - New features batch**

## Features

### Direct control

- Jog Pose.

Control the pose of the robot in position and orientation. The arrows show the direction of the translation and rotation command to apply to the robot's pose. You can adjust the command speed.

- Jog Joints.

Control the joints of the robot in orientation. The arrows show the direction of the rotation command to apply to a joint. You can adjust the command speed.

### State and 3D View

- Frames visualization.

Visualize the TCP and the robot base frames. You can select which frames you want to display.

## Programming

- **Try Except** block in *Logic* group.

The function will try to run what is inside of it. A number of attempts can be set to repeat the operation if the robot failed to do it. You can set the function to continue or stop the program if the attempts do not succeed.

- **Linear** and **Try Linear** in *Move Pose* block, in *Movement* group.

**Linear** function makes the robot move linearly. **Try Linear** function makes the robot move linearly but if no linear trajectory can be computed, it will do a *Standard* movement.

- **Linear** and **Try Linear** in *Shift pose* block, in *Movement* group.

**Linear** function makes the robot do a translation movement linearly. **Try linear** makes the robot do a translation movement linearly but if no linear trajectory can be computed, it will do a *Standard* movement.

- **Move Trajectory** function in *Movement* group.

This function builds a trajectory by giving it waypoints. The robot will try to reach way points in the given order. The *distance smoothing* parameter defines the distance from which the trajectory will be smoothed around each way point.

- **Set TCP** block in *Tool* group.

This function changes the TCP (Tool Center Point) coordinates of the robot. TCP coordinates are used with *Pose*-type movements.

- **Activate/Deactivate TCP** block in *Tool* group.

This function activates or deactivates TCP coordinates. When activated, TCP coordinates are defined by the user. When deactivated, the TCP coordinates are equal to the final point of the robot's arm.

## Settings

- Angles unit selection.  
Angles unit selection between degree and radian from the top bar.
- Setup TCP coordinates from *Robot's settings* view.

Change TCP (Tool Center Point) coordinates. TCP coordinates are used as references for *Pose* movements. *Activate* TCP allows the modification of the TCP coordinates. *Scan* will check if a tool is attached to the arm and will update the TCP coordinates according to the connected tool. *Reset* will reset the TCP coordinates to the actual tool (0 if no tool). You can manually modify the TCP coordinates by writing your orientation and translation. Click *Apply* to apply the new coordinates.

## Saved positions

- Play position  
Move to the desired position from position list by clicking on the *Play* button.

## Vision set

- Camera's settings  
Brightness, saturation and contrast adjustments for the camera.

## Improvements

### Programming

- Graphical improvement of **save program** function in the *Blockly program* view.
- Graphical improvement of **Load program from robot** function in the *Blockly program* view.
- Remove **Clear workspace** in the *Blockly program* view.
- Graphical improvement of **Add position** function in the *Blockly program* view.

## Download Niryo Studio

---


You can download the archive containing Niryo Studio from Niryo Website [here](https://niryo.com/en/download/) (<https://niryo.com/en/download/>). Select the version you need depending on your computer operating system.

## How to launch Niryo Studio

---

Once you have downloaded the .zip file you need to extract the archive. Then go into the new created folder and search for an executable named "NiryoStudio". Double click on the executable to start the application.


**You do not need to install anything on your computer.**



## Overview

### Operating interface

When the Niryo Studio launches, the window below appears.



1. The top toolbar allows you to change the angle units, the language, and connect to Ned
2. 3D Ned visualization or camera stream video
3. State section and mode selection
4. The left menu allows you to switch between the different sections of the application
5. The main window of the application

### Navigation interface

The listed icons from top to bottom on the left menu are:

- Direct command: allows direct motion control of Ned and tools
- Blockly Programming: create, edit and execute a program in Blockly.
- Saved positions
- Program settings: run and upload the robot's program
- Vision Set: set and program the Vision Set
- Conveyor Belt: scan and control the Conveyor Belt
- Settings: robot and software settings
- Logs and Status: displays robot, Niryo Studio logs and the hardware status

## Robot connection

There are three ways to connect Ned to Niryo Studio:

- Wired connection (Ethernet cable)
- Wireless connection (Wi-Fi):
  - Hotspot mode (the robot provides its own Wi-Fi network)
  - Connected mode (the robot is connected on your network)



### Using Ned in Hotspot mode

If the LED on the back panel is blue, that means you can connect to Ned via its own Wi-Fi. In this case Ned IP address is 10.10.10.10.

Ned has its own Wi-Fi network, which starts with "NiryoRobot" followed by a series of numbers and letters. To connect, follow the steps below:

- Connect to the Wi-Fi access point using your laptop with the password "**niryorobot**"
- Click the arrow near "Not connected" to access Ned connection



## ROBOT CONNECTION

Select Robot IP ▾ ▾

10.10.10.10




Connect to Niryo Robot

- Select “Niryo Robot Hotspot” and connect



## ROBOT CONNECTION



- To disconnect from Ned, open the connection panel again and click on “Disconnect”

## ROBOT CONNECTION

Select Robot IP

192.168.1.15

192.168.1.15



Disconnect from Niryo Robot

Connected to 192.168.1.15


## Using Ned on your Wi-Fi network

If Ned is connected to the Wi-Fi network (see the [robot settings](#) (`index.html#ned-settings`) section), the LED on the back panel is green.

To connect Ned to Niryo Studio, please follow the following steps:

- Connect your computer to the same Wi-Fi network as Ned,

- If you do not know your robot IP address, open the connection panel and click on the “Wi-Fi + Search” button. This will search for connected robots inside the current Wi-Fi network. Choose the correct IP address and click “Connect Niryo Robot”,
- If you gave a custom name to Ned, select it by its name,
- If you already know your Ned’s IP address, click the arrow, enter your IP address and connect to the robot.



## Using Ned with an Ethernet cable

You can also use an Ethernet cable to physically connect Ned with your computer.

- Connect Ned to your local network (computer or switch) using the LAN port,
- Modify the selected IP address on the connection panel to “Niryo Robot Ethernet (default IP)”,




- Connect to Ned.

## Using Ned in simulation with Niryo Studio

Simulators such as [Rviz](https://docs.niryo.com/dev/ros/source/simulation.html#control-with-ros) (<https://docs.niryo.com/dev/ros/source/simulation.html#control-with-ros>) or [Gazebo](https://docs.niryo.com/dev/ros/source/simulation.html#launch-simulation) (<https://docs.niryo.com/dev/ros/source/simulation.html#launch-simulation>) must be launched on the same computer as the one running Niryo Studio.

- Click the arrow near “Not connected” to access the robot connection,



- Select “Localhost” and connect.

## Robots specificities

### Calibration

Ned needs to be calibrated in three cases:

- At the start-up
- If the steppers motors lose their offset positions and need to re-calibrate
- If the 3D view and Ned physical movements are not similar

In most of the cases, Niryo Studio will alert you by showing “CALIBRATION NEEDED” on the top right.

#### **Note**

In the case of “CALIBRATION NEEDED”, you can’t perform any command on Ned without calibrating it.

To perform a calibration, click on the button “Request new calibration” in the “ROBOT CALIBRATION” section.

### ROBOT CALIBRATION

**Request new calibration**


After clicking, a “Calibration Needed” warning will appear on top of the 3D view.

When the “Calibration Needed” warning appears, you can click on it to perform an auto calibration or you can click on the “Auto calibration” button in the “ROBOT CALIBRATION” section.

**Calibration Needed**

3D View

Camera Streaming

**ROBOT CALIBRATION****Auto Calibration**

By choosing auto-calibration, Ned will execute an automatic sequence by moving its axis that needs to be calibrated until they reach their maximum position, so a correct offset can be applied to each motor. This sequence should take about 20-30 seconds. When the calibration is over, Ned will come back to a resting position.

**Caution**

Before launching the auto-calibration, make sure that there is no obstacle around Ned. You should put Ned in a proper position before performing the calibration. The position in the picture below is recommended. This will avoid damaging Ned, collision, and any hazard motion.



#### **Note**


If the calibration is not successful, the “Calibration Needed” warning will still be displayed. Try the calibration again and if this new calibration does not succeed, it can mean that your robot has a mechanical/hardware issue.

#### **Hint**

You can check if the calibration has succeeded by moving the robot and examine the robot motion in the 3D view.

## Learning mode

Niryo robots have two operating modes: the manual mode (Learning mode) and the automatic mode. The current mode can be determined visually in Niryo Studio, via the Learning Mode button position.



- **Manual mode/learning mode: when Ned is in learning mode, the toggle on Niryo Studio is on.**

In this mode, the user can guide Ned by hand. In this case, there is no torque applied on the joints and motors are in relaxed mode. This mode is recommended when no program is running.

- **Automatic Mode (Learning Mode is deactivated): The Learning Mode does not work under Automatic Mode.**

It is not possible to move Ned by hand. This mode is active when a program is running or when Ned is performing a motion.


#### **Caution**

Do not try to override the motor's torque when learning mode is disabled: you could damage the robot and the motors. You can force to eliminate the risk or you can unplug the power supply ONLY in case of emergency (Ned blocked or application freezed)

### **💡 Hint**

Activate the learning mode when the program is paused or stopped.

## State and 3D View



### 3D View

You can see here the 3D representation of the robot. You can display TCP (Tool Center Point) and Base link frames.

### States and commands

From here, you can:

- View the current joints angle.
- View the current TCP position and orientation.
- Activate or deactivate the Learning mode.
- Stop the current movement.

- Save the current position of the robot.
- Set the arm speed limit.

## Direct control

ARM COMMAND

	Joint (rad)	Pose	JOG TCP	JOG JOINTS (rad)
Joint 1	0	<input type="range"/>	<input type="button"/>	<input type="button" value="BACK TO HOME"/>
Joint 2	0	<input type="range"/>	<input type="button"/>	<input type="button" value="MOVE JOINTS"/>
Joint 3	0	<input type="range"/>	<input type="button"/>	<input type="button" value="STOP COMMAND"/>
Joint 4	0	<input type="range"/>	<input type="button"/>	<input type="button" value="SET TO CURRENT"/>
Joint 5	0	<input type="range"/>	<input type="button"/>	
Joint 6	0	<input type="range"/>	<input type="button"/>	Select Saved Position ▾

### Arm command

On the ARM COMMAND section, you can move Ned's arm directly either by moving the robot joint individually or by translating / rotating Ned pose.

You will find 4 main subsections that this documentation will cover below, but we will start with the main buttons available.

#### Home Position

Click the “BACK TO HOME” button in order to move Ned to its home position.

#### Select Saved position

If you already saved a position, you can select it via the “Select Saved Position” scrolling menu. Once selected, it will upload joints value if your are in joints view, or pose if you are in pose view.

Then you can click on “MOVE JOINTS” or “MOVE POSE” to move the robot according to the selected joints/pose.

#### Stop command

The “STOP COMMAND” button stops Ned's motion.

#### Joints

This method allows you to control each joint individually from its minimum to its maximum value (see the [Robots specificities](#) (`index.html#document-source/specification`) section).

If a joint reaches its joint limit, it can't be driven any further.

**ARM COMMAND**

	Joints (rad)	Pose	JOG TCP	JOG JOINTS (rad)
Joint 1	0	<input type="range"/>		<b>BACK TO HOME</b>
Joint 2	0	<input type="range"/>		<b>MOVE JOINTS</b>
Joint 3	0	<input type="range"/>		<b>STOP COMMAND</b>
Joint 4	0	<input type="range"/>		<b>SET TO CURRENT</b>
Joint 5	0	<input type="range"/>		
Joint 6	0	<input type="range"/>		Select Saved Position ▾

By clicking on the “SET TO CURRENT” button, you will refresh the interface to the actual joints values.

To perform a move joints, you need to:

- Modify the values of the joints either by writing the desired joint value or by moving the cursor,
- Click on the “MOVE JOINTS” button to execute the move command.

Once the command is done (success or not), you will get a notification at the bottom of the screen.

### Important

At any time, you can cancel the current command execution by clicking on the “STOP COMMAND” button.

## Pose

This panel gives you the possibility to move the pose of Ned's TCP (Tool Center Point).

You can edit the values by clicking on them then use the “MOVE POSE” button.

**ARM COMMAND**

	Joints (rad)	Pose	JOG TCP	JOG JOINTS (rad)
Joint 1	0	<input type="range"/>		<b>BACK TO HOME</b>
Joint 2	0	<input type="range"/>		<b>MOVE JOINTS</b>
Joint 3	0	<input type="range"/>		<b>STOP COMMAND</b>
Joint 4	0	<input type="range"/>		<b>SET TO CURRENT</b>
Joint 5	0	<input type="range"/>		
Joint 6	0	<input type="range"/>		Select Saved Position ▾

## JOG TCP


This panel will help you move Ned in a more visual way.

The TCP (Tool Center Point) of the robot is the center of the end effector.

- Holding a translational arrow will make the robot follow its direction,
- Holding a rotational arrow, will make the robot rotate around the TCP.

The robot will move until you release the button.

You can modify the speed of the jogs. A low value will give you more precision while a higher value will bring you faster to the position.



### **💡 Hint**

The performance of the jog TCP can be impacted by your network connection. We advise you to use an Ethernet connection to get the optimal communication and the best experience with the jog TCP feature.

## JOG JOINTS

You can also switch to jog Joints on the ARM COMMAND panel.

The jog Joints allows you to control each joints of the robot by holding down one arrow.

The robot will move until you release the button.

You can also change the movement speed. A low value will give you more precision while a higher value will bring you faster to the position.




### **💡 Hint**

The performance of the jog joints can be impacted by your network connection. We advise you to use an Ethernet connection to get the optimal communication and the best experience with the jog Joints feature.


## Tool command

- The “SCAN” button will detect and add a tool if there is already a connected tool (with motor).
- **If you have a gripper:**


- You can set the opening and closing speed by writing the value or by dragging the cursor (from 100 and 1000),
- You can send open and close commands to the gripper by clicking “OPEN” and “CLOSE” buttons.
- If you have a Vacuum Pump, the “ON” button will pull the air with the Vacuum Pump and the “OFF” button will pull the air out of the pump.



## Digital I/O panel



You can monitor and set I/O signals from or to Ned. The “DIGITAL I/O PANEL” displays the current state of the I/O “SW1”, “SW2”, “GPIO1”, and “GPIO2” connectors on the back of Ned.



### Note


Check [the Electrical Interface](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#electrical-interface) ([https://docs.niryo.com/product/ned/source/hardware/electrical\\_interface.html#electrical-interface](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#electrical-interface)) for electrical and connection details.

- “SW1” and “SW2” connectors can be used to plug a fan, a motor, etc.
- For “GPIO1” and “GPIO2”, you have, from left to right: GND, 5V, and 3 digital pins. The digital pins can be used to communicate with another device (ex: an Arduino board. Check out this [tutorial](https://docs.niryo.com/applications/ned/source/tutorials/control_ned_arduino.html) ([https://docs.niryo.com/applications/ned/source/tutorials/control\\_ned\\_arduino.html](https://docs.niryo.com/applications/ned/source/tutorials/control_ned_arduino.html)) to learn more about that).

## Switches

The switches are already set as OUTPUT mode.

You can change the “SW1” and “SW2” states to “HIGH” or “LOW” by clicking the corresponding button.



## GPIO

### Warning

Make sure to read the [Safety Precautions](https://docs.niryo.com/product/ned/source/hardware/safety_instructions.html) ([https://docs.niryo.com/product/ned/source/hardware/safety\\_instructions.html](https://docs.niryo.com/product/ned/source/hardware/safety_instructions.html)) and the [hardware manual](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#digital-inputs-outputs) ([https://docs.niryo.com/product/ned/source/hardware/electrical\\_interface.html#digital-inputs-outputs](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#digital-inputs-outputs)) before using the digital I/O of Ned.

A digital pin (GPIO) can be set as an INPUT or OUTPUT mode.

You can change the mode (INPUT or OUTPUT) and the state (LOW or HIGH) for each pin, by clicking on the corresponding mode and state on the “DIGITAL I/O PANEL”.

- When the pin is set as an INPUT, you can read the state in this panel.
- In OUTPUT mode, you can set the OUTPUT state of the pin (LOW or HIGH) by clicking on the corresponding button.

#### **💡 Note**

By default, all the GPIO pins are in INPUT mode and HIGH state.

## Electromagnet control

Make sure to read the [Safety Precautions](https://docs.niryo.com/product/ned/source/hardware/safety_instructions.html) ([https://docs.niryo.com/product/ned/source/hardware/safety\\_instructions.html](https://docs.niryo.com/product/ned/source/hardware/safety_instructions.html)) and the [Hardware Manual](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#digital-inputs-outputs) ([https://docs.niryo.com/product/ned/source/hardware/electrical\\_interface.html#digital-inputs-outputs](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#digital-inputs-outputs)) before using the digital I/O of Ned.

The Electromagnet can be set up and controlled on the “DIGITAL I/O PANEL”.

You can connect the electromagnet on the GPIO1 or GPIO2 at the back of Ned (see [Hardware Manual](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#electromagnet) ([https://docs.niryo.com/product/ned/source/hardware/electrical\\_interface.html#electromagnet](https://docs.niryo.com/product/ned/source/hardware/electrical_interface.html#electromagnet))).

Then, you need to set up the digital I/O pin.

- Select the pin where you connected your electromagnet cable.
- Set this pin to OUTPUT.
- Control your electromagnet by changing the state of the selected pin to LOW or HIGH.


#### **💡 Hint**

If you are using the connector provided with the electromagnet, **select 1A pin if you connect to the GPIO1 and 2A if you connect to the GPIO2**.


## Programming

### Project panel

Click the project button on the left menu to start creating and editing Blocky programs.




In the figure below, the top is the project Toolbar, and the left is the Blockly menu.



1. Create a new project: click the “NEW” button to create a new project.
2. Save program: click the “SAVE” button, to save a program. You can save a program on your computer or directly on Ned as an XML file.
3. Load project: click the “LOAD” button to import a Blockly program either from computer or from Ned. The file is a xml.
4. Undo. You can also use CTRL + Z.
5. Redo. You can also use CTRL + MAJ + Z.

6. Play the program displayed on the workspace. Once the program is done (success or not), you will get a notification on the bottom of the screen.
7. Stop the current program execution.
8. Add a position block. This will open a pop-up:



This is the most useful way to add a position command block to your workspace. You can choose between a "Move Joint" block and a "Move Pose" block, and use either the current robot state, or a previously saved position.

9. Select the save mode when using the top button on Ned.

#### **💡 Hint**

While the Niryo Block interface is open, you can add a "Move joint" or "Move pose" block with current joints values and pose values just by pressing and releasing the top button of the robot.

10. Help.
11. Clicking on this icon will center the workspace on your blocks.
12. Workspace zoom control.
13. Workspace zoom control.
14. To delete a block, simply drag it and drop it onto the trash. You can also select it and press the delete key on your keyboard.
15. The "Niryo Robot" functions in the block library cover almost all of the possible commands you can execute on Ned:
  - Arm commands
  - Movement commands
  - Tool commands
  - I/O setup and control
  - Vision commands
  - Conveyor control

#### **💡 Note**

The "NIRYO BLOCKS" panel contains everything you need to create complete programs for Ned. This is probably the panel on which you will spend most of your time.

16. This is the workspace. Your whole program will be there.

#### **💡 Note**


You cannot execute a new sequence while another sequence is running. If you want to start a new sequence, you will have to stop the previous one.

#### **💡 Hint**

Niryo blocks along with the learning mode is the perfect combination to easily create programs in no time. Activate the learning mode, then for each position you want, add a position block corresponding to the current state, and that's it.

## Saved positions

In the “SAVED POSITION” panel, you can find all the saved positions.



When you select a position, you can see its properties:


- Name
- Description
- Values: list of joints values in radian and the position and orientation of the TCP (Tool Center Point).

You can:

- Edit the position’s name and description,
- Move the robot to the desired position,
- Delete the position.

## Program settings

In the “Program Settings”, you can set an autorunning program, upload a Python program and display a program.



## Program Autorun


The “PROGRAM AUTORUN” functionality allows you to execute a set of sequences, directly by pressing the top button of the robot.

Here is how it works:

- Program a sequence with Blockly.
- Save the sequence on Ned.
- Add the sequence to the sequence autorun.
- Now trigger the sequence with the top button, without having to open Niryo Studio.

### **💡 Hint**

This is very useful when you use Ned in an exhibition event or for a demo. You will just need to power on Ned, wait for the LED to turn blue or green, and then press the button to run a given set of programs (that you previously programmed).



1. Display the default program triggered with the button
2. Select the program
3. Choose the “PROGRAM AUTORUN” mode. You have choices:
  - a. ONE SHOT: when you press and release the top button, the given set of sequences will be run once. After that, Ned will go back to a “resting position” and activate the “Learning Mode”. You can start another run by pressing the button again. If you press the button while the set of sequences is running, it will stop the execution.
  - b. LOOP: when you press and release the top button, the selected program will run in a loop, forever. You can press the button again to stop the execution.
4. SELECT button: select the program autorun.

## **Upload program**

You can upload a Python script from your computer to Ned.

## **Programs list**

You can pick a previously saved program from the select box. Click on a sequence to see all the properties and actions available for this program.

You can display the Python code of the program.

## PROGRAMS

```
# !/usr/bin/env python
from niryo_robot_python_ros_wrapper import *
import sys
import rospy
rospy.init_node('niryo_blockly_interpreted_code')
n = NiryoRosWrapper()
n.calibrate_auto()
try:
    n.move_pose(*[0.119, 0, 0.209, -0.012, 0.72, 0.001])
except NiryoRosWrapperException as e:
    sys.stderr.write(str(e))
```

After the program details, you get a series of available actions:

1. Play the selected program (same as the “PLAY” button in Niryo blocks).
2. Stop the current program execution (same as the “STOP” button in Niryo blocks).
3. Open the program in Niryo blocks. This will make you switch to the “Niryo blocks” panel, and the program will be added to the current workspace.

**💡 Hint**

This functionality is very useful when you want to duplicate and create a new program from an existing one.

4. Edit the program. You can modify the name, the description, and the Blockly XML.
5. Delete the program.

## Addons

### Conveyor Belt

The use of the Conveyor Belt is detailed [in its documentation](https://docs.niryo.com/product/conveyor-belt/source/how_to_use.html#on-niryo-studio) ([https://docs.niryo.com/product/conveyor-belt/source/how\\_to\\_use.html#on-niryo-studio](https://docs.niryo.com/product/conveyor-belt/source/how_to_use.html#on-niryo-studio)).


### Vision Set

Discover how to use the vision in the [Vision Set documentation](https://docs.niryo.com/product/vision-set/source/setup.html#software) (<https://docs.niryo.com/product/vision-set/source/setup.html#software>).

## Robot's settings and update

### Ned settings

Ned's related parameters can be set in the “ROBOT SETTINGS” panel.



## Robot Name

You can give a custom name to Ned.

- This name will be used when Ned is in hotspot mode to create the Wi-Fi network name.
- In Niryo Studio, when Ned is in connected mode, you will see that Ned is represented both by its IP address and its custom name.

You can set a custom name under “ROBOT SETTINGS”. After choosing the name, you need to reboot Ned.

## Robot TCP (Tool Center Point)

The TCP (Tool Center Point) of the robot corresponds to the working point of the tool. From here, you can change the frame of the TCP.

**ROBOT TCP (TOOL CENTER POINT)**

No tool selected SCAN

Position (m)		Orientation (rad)	
X	0	Roll	0
Y	0	Pitch	0
Z	0	Yaw	0

**Enable TCP** RESET SET

When the TCP is enabled, every position move is set corresponding to the TCP frame.

When you plug a tool and if the TCP is enabled, the TCP frame will automatically update to the corresponding tool TCP.

## Ned Wi-Fi setting

If Ned is already connected to a Wi-Fi network and you want to change your Wi-Fi network:

- Connect to Ned in “Hotspot mode” or in “Connected mode”



**Note**

If you are not connected to Ned, you will not be able to change its Wi-Fi network.


- Under “ROBOT SETTINGS”, check your current IP address.
- Click the Wi-Fi setup button and fill in the Wi-Fi name (SSID) and password of the Wi-Fi network you want Ned to connect to. Make sure to avoid typos in the name and password. Then, click on “Connect robot to Wi-Fi”.
- Ned will try to connect to the Wi-Fi, depending on the result, the robot's reaction changes:
  - If successfully connected: the robot will reboot itself, the LED will become red then green. Connect to the same wifi network and then connect to Ned (Help: [Using Ned on your Wi-Fi network](#) (index.html#using-ned-on-your-wi-fi-network).)
  - If failed to connect to the SSID: the LED will remain blue (the LED didn't change after ~30 seconds). Verify your password and try to connect again.
  - If the SSID can't be found currently by the Raspberry Pi: you won't be disconnected from Niryo Studio, but an error notification will be shown. You can retry within a few seconds (the wi-fi card may be disturbed by other wi-fi devices) or move your Raspberry Pi closer to the Wi-fi router.

**Warning**

If you changed Ned network configuration from Wi-fi to hotspot, troubles can occur when reconnecting via Niryo Studio, please reboot your robot in this case.

**Language selection**


The default language of Niryo Studio is English but you can change the current language with the language button on the top toolbar.

**Note**


The Blockly editor and the logs remain in English.

**Angle unit selection**

The default angle unit in robotics is radian. But you can change to degree with the corresponding button in the top toolbar.



## Raspberry PI settings



ROS logs:

In this section, you can see:

- The available disk size on the Raspberry Pi. If the number is displayed in orange or red, it means there is no space enough on the micro-SD card.
- The size of the stored ROS logs: after some times, they can take a lot of space.
- You can choose to remove all previous ROS logs when Ned boots. This is very useful for a normal usage of Ned. Basically, you do not need all the ROS logs that are produced and stored, unless you are developing or debugging a specific point (in this case you might want to switch the option to "NO" and then turn it back to "YES" when you're done).

### **Caution**

If you never purge ROS logs (on boot or manually), the micro-SD card will eventually become full and might not work as expected. So only change the option to "NO" if you know what you are doing. In any case, if you do not really understand this functionality, always leave the option to "YES".

- You can also remove all ROS logs immediately by clicking on the "Purge ROS logs" button.
- Shutdown Ned.

### Shutdown Robot

**SHUTDOWN NIRYO ROBOT**

or **REBOOT**

- You can shut down Ned with Niryo Studio. The "SHUTDOWN NIRYO ROBOT" will shut down the Raspberry Pi (same behaviour as if you press the physical top button for 3 seconds). Be sure to wait for the LED to turn red before completely powering off Ned.
- You can reboot Ned by clicking "Reboot".

### **Note**

If you reboot Ned, do not power it off. The LED will turn to purple, then red, then will become blue or green again.

## Software versions

You can check that Ned and relative softwares are up to date by comparing the current and latest version for each of them.

- Niryo Studio: the desktop app to control Ned.
- Robot Software: the current version of the Niryo robot ROS stack running on Ned.
- Niryo Steppers motors: a list of the firmware versions running on each NiryoStepper attached to Ned.  
[Click here to update the firmware version](https://niryo.com/docs/niryo-one/update-your-robot/update-niryo-steppers/) (<https://niryo.com/docs/niryo-one/update-your-robot/update-niryo-steppers/>).

### Note

All the NiryoStepper motors must have the same version.

- Latest available versions (directly updated from our website).

### Note

If an incompatible software version is detected, if the version you have is not supported anymore or an update is available, you will get a pop-up notification when you connect to Ned. Please follow the displayed instructions.



## Ned software Update

Ned software inside the robot can be updated through Niryo Studio in only a few steps.

First, your Ned must have an access to internet, so you have to connect it to an ethernet (you have to edit via ssh the network configuration of Ned) or a [Wi-fi](#) ([index.html#using-ned-on-your-wi-fi-network](#)) network that provides internet.

When your robot is properly configured, at the next boot of the robot, it will check if an update is available. A pop-up will then be displayed at each following connection from Niryo Studio with Ned to inform you that an update is available for your robot.

You then only have to click on the update button of the pop-up or the update button in the [software versions](#) panel.

### Warning

Ned will reboot itself at the end of the process (in case of update success) so it's recommended **not** to do any other operation with Ned before starting an update.

Once the update has started, a pop-up will indicate that Ned is updating and that it's downloading / applying the update.



**Updating...**

Please wait until software update has completed and that this popup closes automatically. (Update time may vary depending of your internet connection speed).

**Note**

Update time is heavily impacted by the internet speed of your connection so if you connect to a wi-fi network, be sure that the raspberry is not too far away from the wi-fi modem.

Once the update has successfully ended, another pop-up will appear (closing the previous one) telling that your robot will reboot and you may be disconnected from Niryo Studio. Now you just have to wait for the robot to initialize properly (the red LED turning into green or blue) before reconnecting through Niryo Studio, the update process is done.

Robot will now reboot in a few seconds, you will be disconnected from Niryo Studio. Please wait until the robot has properly rebooted (LED color change from red to blue or green) in order to reconnect to the robot.

**CLOSE**

**Warning**

In case the update failed, it could be a faulty connection, try to do the update once more and if you still encounter a problem, please contact us at [contact@niryo.com](mailto:contact@niryo.com) (mailto:contact%40niryo.com).

## Logs and Status

### Logs

At the top of the Robot Status page, you can find 2 sections:

- Niryo Studio logs: logs produced by the desktop application.
- ROS logs: logs produced internally by Ned.

#### Hint

- You can save Ned logs or Niryo Studio logs on a text file.
- You can send this file to our technical support to help you in case of problem.

### Hardware status

In the “HARDWARE STATUS” panel, you can find the hardware status of Ned.

In this section, you can find details about the controller temperature, motors voltage, temperature, and error code.

HARDWARE STATUS									
Niryo Robot version : Ned	①	Raspberry Pi Temperature : 45 °C	②	Tool Reboot	⑨	Reboot Motors	⑦	Start Motors Report	⑧
Motor	③	Volt	④	°C	⑤	Error code	⑥		
Joint_1	-	-	31	-	-	-	-		
Joint_2	-	-	27	-	-	-	-		
Joint_3	-	-	30	-	-	-	-		
Joint_4	11.3	-	32	-	-	-	-		
Joint_5	11.2	-	32	-	-	-	-		
Joint_6	7.3	-	29	-	-	-	-		
Tool	7.2	-	29	-	-	-	-		

1. Ned hardware version.
2. Raspberry Pi temperature: the temperature of the Raspberry Pi is displayed in real time. If it's displayed in orange or red, that means that something is not normal and the board is too hot. Make sure that your board temperature does not exceed 70°C.
3. List of connected motors: you should have 6 motors (from Joint\_1 to Joint\_6), if you have a connected tool, it should appear on the list.
4. Displayed measured voltage of servo motors (Dynamixels).
5. Measured temperature of each joint.
6. Error code: displays error code if there is any. Check troubleshooting chapter to resolve the problem.
7. Reboot motors: this button reboots servo motors without rebooting the entire robot, this can help you to resolve motors errors.
8. Start motors report: this option helps you to debug Ned. Click the Motor Report button and Start Motors Report.
9. Tool reboot: this option allows you to reboot the tool's motor. This button is useful when a motor error occurs.



## WARNING

=> Please move the robot in home position and make sure to have space around the robot.

Learning Mode

=> Do not touch the robot until report completion.

**Start Motors Report**




**Stop Motors Report**

**Save logs file**

## Logic

Concept	Explanation	Example
If/Else If/Else block	If the first value is true, then does the first block of statements. Otherwise, if the second value is true, does the second block of statements. If none is true, does the last block of statements.	
Mathematical comparison operator block	Compares two values with an operator [=, ≠, <, ≤, >, ≥]. Returns true if the comparison is true.	
Logical comparison operator block	Compares two values with an operator [and / or]. Returns true if the comparison is true.	
Not block	Returns true if the input is false, and false if the input is true.	
True/False block	Returns either true or false.	
Null block	Returns null.	
Test block	Checks the condition in test block statement. If the condition is true, returns the value of the 'if true' value. Otherwise, returns the 'if false' value.	
Try/Except block	On the "Try" statement you can ask Ned to try to execute any type of action a certain amount of time. You can then choose if Ned should CONTINUE or STOP in case he is "On failure". If not Ned will execute the next block.	

## Loops






Concept	Explanation	Example
Time loop block	Does one or several statements multiple times.	
While loop block	While a value is true, then do some statements.	
For loop block	Executes a loop based on a variable, from a starting number to an ending number with a specified interval, and executes the specified blocks.	
For each item in list loop block	For each item in a list, sets the variable "v" to the item, and then does some statements.	
Break block	Breaks out of the containing loop. Must be in a loop block.	

## Math


Concept	Explanation	Example
Set number block	Sets a value.	
Mathematical operation block	Returns the value of the specified operation with two numbers. [+, -, ×, ÷, ^].	
Check operation block	Checks if a number is even, odd, prime, whole, positive, negative or divisible.	
Round block	Rounds a number up or down.	
Mathematical list operations block	Returns the mathematical operation value of all number in the list [sum, min, max, average, median, modes, standard deviation, random item].	
Division remainder block	Returns the remainder from the division of two numbers.	
Constrain number block	Constrains a number to be between the specified limits (inclusive).	
Random integer block	Returns a random integer between the two specified limits (inclusive).	

## Lists


Concept	Explanation	Example

Concept	Explanation	Example
Create list block	Creates an empty list.	 create empty list
Create list with items block	Creates a list with a specified number of items.	 create list with
Create list with one item block	Creates a list consisting of the given value repeated the specified number of times.	 create list with item  5 times
List length block	Returns the length of a list.	 length of
Is list empty block	Returns true if the list is empty.	 is empty
Item index finder block	Returns the index of the <i>first</i> or the <i>last</i> occurrence of the item in the list. Returns 0 if the item is not found.	 in list  find  occurrence of item 
Get item block	<i>Returns / returns and removes / removes</i> the item at the specified position in a list. #1 is the first item.	 in list  get  # 
Set item index block	<i>Sets / inserts at</i> the item at the specified position in a list. #1 is the first item.	 in list  set  #  as 
Copy list-portion block	Creates a copy of the specified portion of a list.	 in list  get sub-list from  to 
Sort list block	<i>Sorts numeric / alphabetic / alphabetic ignore case by ascending / descending</i> a copy of a list.	 sort  ascending 

## Variables

Concept	Explanation	Example
Create variable block	Creates a variable by its name.	 Create variable...
Set variable block	Sets this variable to be equal to the input.	 set  to 
Change variable block	Changes this variable by the input.	 change  by  1
Get variable block	Returns the variable.	 variable 

## Functions

Concept	Explanation	Example
Create no output function block	Creates a function with no output. You can add input arguments.	 to function

Concept	Explanation	Example
Create function with output block	Creates a function with an output. You can add input arguments.	
Function return block	If a value is true, then ends the function and returns the return value. Must be inside a function.	
Call function block	Runs the function.	

## Utility

Concept	Explanation	Example
Wait time block	Creates a break time in the program.	
Comment block	Adds comments to the code. This block will not be executed. Note: accents are not accepted (é, à, è,...).	
Break Point block	Stops the execution of the program. Press the "Play" button to resume.	

## Arm

Concept	Explanation	Example
Learning mode block	Activates / Deactivates the learning mode.	
Set arm speed block	Sets the arm speed.	

## Movement

Concept	Explanation	Example
Joints block	Creates an object pose according to the robot's joints values.	
Move joints block	Moves the robot according to a Joints block.	
Pose block	Creates an object pose according to the end effector's cartesian coordinates.	
Move pose block	Moves the robot according to a pose block. You can change between "Standard" move (Point to Point), "Linear", the robot will follow a linear trajectory from its position to the desired position, or "Try linear", the robot will try to do a "Linear" trajectory, but if it can't compute the linear trajectory, it will do a "Standard" move.	



Concept	Explanation	Example
Shift block	Shifts the robot pose according to an axe <i>x / y / z / roll / pitch / yaw</i> . You can change between "Standard" move (Point to Point), "Linear", the robot will follow a linear trajectory from its position to the desired position, or "Try linear", the robot will try to do a "Linear" trajectory, but if it can't compute the linear trajectory, it will do a "Standard" move.	
Pick from pose block	Moves the robot's TCP (arm's end point + tool coordinate) to a specified pose and activates the tool to pick an object.	
Place from pose block	Moves the robot's TCP (arm's end point + tool coordinate) to a specified pose and deactivates the tool to place an object.	
Move Trajectory block	The robot will pass through the list of goals with the desired distance smooth. The distance smooth is the radius from the goal where the robot will start to go towards the next point.	

## I/O

Concept	Explanation	Example
Get I/O block	Returns the number of the pin.	
Set I/O mode block	Sets I/O pin mode to <i>input</i> or <i>output</i> .	
Set output state block	Sets output pin state to <i>high</i> or <i>low</i> .	
Get input state block	Returns the input pin state.	
State block	Returns the state value <i>high</i> or <i>low</i> .	
Set Switch state block	Sets the 12V switch state to <i>high</i> or <i>low</i> .	

## Tool

Concept	Explanation	Example
Scan tool block	Scans and updates current tool.	
Grasp block	Activates the <i>gripper / vacuum</i> .	
Release block	Deactivates the <i>gripper / vacuum</i> .	
Open gripper block	Opens the gripper at a certain speed.	
Close gripper block	Closes the gripper at a certain speed.	

Concept	Explanation	Example
Pull air block	Pulls the air in the vacuum pump.	
Push air block	Pushes the air in the vacuum pump.	
Setup electromagnet block	Selects on which pin the electromagnet is connected.	
Activate electromagnet block	Activates the electromagnet power.	
Deactivate electromagnet block	Deactivates the electromagnet power.	
Activate / Deactivate TCP block	Activates or deactivates the TCP.	
Set TCP block	Sets TCP frame	

## Vision

The Vision blocks and Vision templates are detailed in the [Vision Set Documentation](#) ([https://docs.niryo.com/product/vision-set/source/how\\_to\\_use\\_blockly.html#blocks-description](https://docs.niryo.com/product/vision-set/source/how_to_use_blockly.html#blocks-description)).

## Conveyor blocks

The Conveyor blocks are detailed in the [Conveyor Belt documentation](#) ([https://docs.niryo.com/product/conveyor-belt/source/how\\_to\\_use.html#blockly](https://docs.niryo.com/product/conveyor-belt/source/how_to_use.html#blockly)).

## Advanced Programming

There are many other ways to develop on Ned.

For more information, please refer to the following links:


- [Python](#) (<https://docs.niryo.com/dev/pyniryo/index.html>)
- [ROS](#) (<https://docs.niryo.com/dev/ros/index.html>)
- [Modbus](#) (<https://docs.niryo.com/dev/modbus/index.html>)

[Suggest a modification](#)

[Download as PDF](#)

## Ned ROS documentation

This documentation contains everything you need to understand Ned's functioning & how to control it through ROS. It is made as well for users who are using the "physic" robot as those who want to use a virtual version.



### Preamble

Before diving into the software documentation, you can learn more about the robot development in the [Overview](#) ([index.html#document-source/overview](#)) section.

Also, if you want to use a simulated version of Ned, browse the [Simulation](#) ([index.html#document-source/simulation](#)) section.

### Ned Control via ROS

Ned is fully based on ROS.

#### ROS Direct control

##### ! Important

To control the robot directly with ROS, you will need either to be connected in SSH to the physic robot, or to use the simulation.

ROS is the most direct way to control the robot. It allows:

- to send command via the terminal in order to call services, trigger action, ...
- to write a entire Python/C++ node in a script to realize a full process.

See [ROS](#) ([index.html#document-source/ros](#)) section to see all Topic & Services available.

#### Python ROS Wrapper

##### ! Important

To use Python ROS Wrapper, you will need either to be connected in SSH to the physic robot, or to use the simulation.

The Python ROS Wrapper is built on top of ROS to allow a faster development than ROS. Programs are run directly on the robot which allows to trigger them with the robot's button once a computer is no longer needed.

See [Python ROS Wrapper](#) ([index.html#document-source/ros\\_wrapper](#)) to see which functions are accessible and examples on how to use them.

#### More ways


Other methods are available to control the robot allowing the user to code and run programs outside his terminal.

Learn more on this [section](#) ([index.html#document-source/more](#)).

#### ROS Stack overview

Ned is a robot based on Raspberry, Arduino & ROS. It uses ROS to make the interface between Hardware and high-level bindings.

On the following figure, you can see a global overview of the Niryo's robot software in order to understand where are placed each part of the software.




Niryo Robot v3 Software  
Global Overview

Niryo Robot  
ROS Stack

Niryo robot v3 software

### Use Niryo robot through simulation

The simulation allows to control a virtual Ned directly from your computer.



Ned in Gazebo Simulation

In this tutorial, you will learn how to setup simulation on a computer.

#### Note

You can use [Niryo Studio with the simulation](https://docs.niryo.com/product/ned/source/software/niryo_studio.html#connecting-simulation-to-niryo-studio/) ([https://docs.niryo.com/product/ned/source/software/niryo\\_studio.html#connecting-simulation-to-niryo-studio/](https://docs.niryo.com/product/ned/source/software/niryo_studio.html#connecting-simulation-to-niryo-studio/)). To do so, you just have to connect Niryo Studio to "localhost".

### Simulation environment installation

#### Attention

The whole ROS Stack is developed and tested on ROS Melodic which requires Ubuntu 18.04 to run correctly. The using of another ROS version or OS may lead to malfunctions of some packages.

To allow the simulation to run on your computer, you will need to install ROS and some packages.

Installation index:

- Prepare environment
- Install ROS dependencies
- Setup Ned ROS environment

#### Prepare environment

**Note**

All terminal command listed are for Ubuntu users.

Place yourself in the folder of your choice and create a folder **catkin\_ws\_niryo\_ned** as well as a sub-folder **src**:

```
mkdir -p catkin_ws_niryo_ned/src
```

Then go to the folder **catkin\_ws\_niryo\_ned** and clone Ned repository in the folder **src**. For the future operation, be sure to stay in the **catkin\_ws\_niryo\_ned** folder:

```
cd catkin_ws_niryo_ned
git clone https://github.com/NiryoRobotics/ned_ros src
```

**Install ROS dependencies****Install ROS**

You firstly need to install ROS Melodic. To do so, follow the ROS official tutorial[here](http://wiki.ros.org/melodic/Installation/Ubuntu) (<http://wiki.ros.org/melodic/Installation/Ubuntu>) and chose the **Desktop-Full Install**.

**Install additional packages**

To ensure the functioning of all Ned's packages, you need to install several more packages:

**Method 1: Quick installation via ROSDep**

For each packages, we have referenced all the dependencies in their respective *package.xml* file, which allow to install each dependency via *rosdep* command:

```
rosdep update
rosdep install --from-paths src --ignore-src --default-yes --rostdistro melodic --skip-keys "python-rpi.gpio"
```

**Method 2: Full installation**

ROS packages needed are:

- build-essential
- catkin
- python-catkin-pkg
- python-pymodbus
- python-rostdistro
- python-rosPKG
- python-rosdep-modules
- python-rosinstall python-rosinstall-generator
- python-wstool

To install a package on Ubuntu:

```
sudo apt install <package_name>
```

Melodic specific packages needed are:

- moveit
- control
- controllers
- tf2-web-republisher
- rosbridge-server
- joint-state-publisher-gui

To install a ROS Melodic's package on Ubuntu:

```
sudo apt install ros-melodic-<package_name>
```

**Setup Ned ROS environment****Note**

Be sure to be still placed in the **catkin\_ws\_niryo\_ned** folder.

Then perform the **make** of Ned's ROS Stack via the command:

```
catkin_make
```

If no errors occurred during the **make** phase, the setup of your environment is almost complete!

It is necessary to source the configuration file to add all Ned packages to ROS environment. To do so, run the command:

```
source devel/setup.bash
```

It is necessary to run this command each time you launch a new terminal. If you want to make this sourcing appends for all your futur terminals, you can add it to your **bashrc** file:

```
echo "source $(pwd)/devel/setup.bash" >> ~/.bashrc
source ~/.bashrc
```

Installation is now finished!

## Simulation utilization

### Important

- If you haven't follow the step of [Simulation Environment Installation](#), you won't be able to use the simulation.
- Hardware features won't be accessible.

The simulation is a powerful tool which allow to test new programs directly on your computer which prevent to transfer new code on the robot. It also helps for developing purpose → no need to transfer code, compile and restart the robot which is way slower than doing it on a desktop computer.

### Without physics - Visualization


The visualization happens with Rviz which is a powerful tool.

#### Control with trackbar

This visualization allows an easy first control of the robot, and helps to understand joints disposal. You can access it by using the command:

```
roslaunch niryo_robot_description display.launch
```

Rviz should open with a window containing 6 trackbars. Each of these trackbars allows to control the corresponding joint.




Example of trackbars use.

#### Control with ROS

Not only Rviz can display the robot, it can also be linked with ROS controllers to show robot's actions from ROS commands! This method can help you debugging ROS topics, services and also, API scripts.

To run it:

```
roslaunch niryo_robot_bringup desktop_rviz_simulation.launch
```



Rviz opening, with the robot ready to be controlled with ROS!

#### With physics - Simulation

For the simulation, Ned uses Gazebo, a well known tool among the ROS community. It allows:

- Collision.
- World creation → An virtual environment in which the robot can deal with objects.
- Gripper & Camera using.

The Niryo Gripper 1 has been replicated in Gazebo. The Camera is also implemented.

##### **Note**


Gazebo also generates camera distortion, which brings the simulation even closer from the reality!

Launch simulation

A specific world has been created to use Ned in Gazebo with 2 workspaces.

To run it:

```
roslaunch niryo_robot_bringup desktop_gazebo_simulation.launch
```



Gazebo view, with the robot ready to be controlled with ROS!

#### **Note**

You can edit Gazebo world to do your own! It's placed in the folder `worlds` of the package `niryo_robot_gazebo`.

#### Simulation option

The user can disable 3 things by adding the specific string to the command line:

- the Gazebo graphical interface: `gui:=false`.
- the Camera & the Gripper - Vision & Gripper wise functions won't be usable: `gripper_n_camera:=false`.

#### **Hint**

Gazebo can be very slow. If your tests do not need Gripper and Camera, consider using Rviz to alleviate your CPU.

### ROS Stack documentation



ROS (Robot Operating System) is an Open-Source Robotic Framework which allows to ease robot software development. The framework is used in almost each part of Ned software.

The high-level packages (motion planner, vision, ...) are coded in Python to give good readability whereas communication with Hardware is developed in C++ to ensure speed.

#### **Note**


To learn more about ROS, go on [Official ROS Wiki](http://wiki.ros.org/) (<http://wiki.ros.org/>).

In this section, you will have access to all information about each Niryo Robot's ROS packages.

#### **Niryo\_robot\_bringup**

This packages provides config and launch files to start Ned and ROS packages with various parameters.

Launch files are placed in the `launch` folder. Only files with `.launch` extension can be executed.



## On RaspberryPI

### Ned

The file **niryo\_ned\_robot.launch** allows to launch ROS on a Raspberry Pi 4.  
This file is automatically launched when Ned boots (Ned RPi4B image).

Command to launch Ned's ROS Stack:

```
roslaunch niryo_robot Bringup niryo_ned_robot.launch
```

## On Desktop (Simulation)

As the simulation happens on a computer, the hardware-related stuff is not used.

### For both of following launch files, you can set:

- *gui* to "false" in order to disable graphical interface.

### Gazebo simulation

**Run Gazebo simulation. The robot can do everything that is not hardware-related:**

- move, get\_pose.
- use the camera (to disable it, set "camera" parameter to 'false').
- use the Gripper 1 (to disable it, set "simu\_gripper" parameter to 'false').
- save/run programs, go to saved pose, ...

Command to launch the simulation:

```
roslaunch niryo_robot Bringup desktop_gazebo_simulation.launch
```

To disable camera & gripper:

```
roslaunch niryo_robot Bringup desktop_gazebo_simulation.launch gripper_n_camera:=false
```

### Rviz simulation

Run Rviz simulation. You can access same features as Gazebo except Camera & Gripper.

To run it, use the command:

```
roslaunch niryo_robot Bringup desktop_rviz_simulation.launch
```

## Notes - Ned Bringup

*niryo\_robot\_base* files setup many rosparams, these files should be launched before any other package.

### Niryo\_robot\_arm\_commander

This package is the one dealing with all commander related stuff.

It is composed of only one node, which is running separately the arm commander and the tool commander.

#### Commander node

##### The ROS Node is made to interact with:

- The Arm through MoveIt!
- The tools through the tool controller.

##### All commands are firstly received on the actionlib server which:

- Handles concurrent requests.
- Checks if the command can't be processed due to other factors (ex: learning mode).
- Validates parameters.
- Calls required controllers and returns appropriate status and message.

The namespace used is: **/niryo\_robot\_arm\_commander/**

Parameters - Commander

*Commander's Parameters*

Name	Description
reference_frame	Reference frame used by MoveIt! for moveP. Default : 'world'
move_group_commander_name	Name of the group that MoveIt is controlling. By default: "arm"
jog_timer_rate_sec	Publish rate for jog controller
simu_gripper	If you are using the simulated Gripper and want to control the Gripper

## Action Server - Commander¶

*Commander Package Action Servers*

Name	Message Type	Description
robot_action	RobotMove	Command the arm and tools through an action server

## Services - Commander¶

*Commander Package Services*

Name	Message Type	Description
is_active	GetBool	Indicate whereas a command is actually running or not
stop_command	Trigger	Stop the actual command
set_max_velocity_scaling_factor	SetInt	Set a percentage of maximum speed
/niryo_robot/kinematics/forward	GetFK	Compute a Forward Kinematic
/niryo_robot/kinematics/inverse	GetIK	Compute a Inverse Kinematic

## Messages - Commander¶

*Commander Package Messages*

Name	Description
ArmMoveCommand	Message to command the arm
ShiftPose	Message for shifting pose
PausePlanExecution	Pause movement execution

All these services are available as soon as the node is started.

## Dependencies - Commander¶

- actionlib
- actionlib\_msgs
- control\_msgs
- geometry\_msgs
- MoveIt!
- moveit\_msgs
- Niryo\_robot\_msgs
- Niryo robot tools commander package
- python-numpy
- ros\_controllers
- rosbridge\_server
- sensor\_msgs
- std\_msgs
- tf2\_web\_republisher
- trajectory\_msgs

## Action, services &amp; messages files - Commander

## RobotMove (Action)¶

```
# goal
niryo_robot_arm_commander/ArmMoveCommand cmd
---
# result
int32 status
string message
---
# feedback
niryo_robot_msgs/RobotState state
```

## GetFK (Service)¶

```
float32[] joints
---
niryo_robot_msgs/RobotState pose
```

## GetIK (Service)¶

```
niryo_robot_msgs/RobotState pose
---
bool success
float32[] joints
```

## JogShift (Service)¶

```
int32 JOINTS_SHIFT = 1
int32 POSE_SHIFT = 2

int32 cmd

float32[] shift_values

...
int32 status
string message
```

**ArmMoveCommand (Message)¶**

```
int32 JOINTS = 0      # uses joints
int32 POSE = 1        # uses position and rpy
int32 POSITION = 2   # uses position
int32 RPY = 3         # uses rpy
int32 POSE_QUAT = 4  # uses position and orientation
int32 LINEAR_POSE = 5 # uses position and rpy
int32 SHIFT_POSE = 6 # uses shift
int32 SHIFT_LINEAR_POSE = 7 # uses shift
int32 EXECUTE_TRAJ = 8 # uses dist_smoothing, list_poses
int32 DRAW_SPIRAL = 9

int32 cmd_type

float64[] joints
geometry_msgs/Point position
niryo_robot_msgs/RPY rpy
geometry_msgs/Quaternion orientation
niryo_robot_arm_commander/ShiftPose shift

geometry_msgs/Pose[] list_poses
float32 dist_smoothing
```

**PausePlanExecution (Message)¶**

```
int8 STANDBY = 0
int8 PLAY = 1
int8 PAUSE = 2
int8 RESUME = 3
int8 CANCEL = 4

int8 state
```

**ShiftPose (Message)¶**

```
int32 axis_number
float64 value
```

**Niryo\_robot\_description**

This package contains URDF files and meshes (collada + stl) for Ned.

To display Ned on Rviz:

```
roslaunch niryo_robot_description display.launch
```

**Note :** 3D visualization is not available on Ned Raspberry Pi4 image. To use the following commands, you must have setup Ned ros stack on your computer.

**Niryo\_robot\_gazebo****Usage**

This package contains models, materials & Gazebo worlds.

When launching the Gazebo version of the ROS Stack, the file `niryo_robot_gazebo_world.launch.xml` will be called to generate the Gazebo world.

**Create your own world**

Create your world's file and put it on the folder `worlds`. Once it is done, you have to change the parameter `world_name` in the file `niryo_robot_gazebo_world.launch.xml`.

You can take a look at the Gazebo world by launching it without robot by precising the world name in the arg `world_name`:

```
roslaunch niryo_robot_gazebo niryo_gazebo_world.launch world_name:=niryo_cube_world
```

**Niryo\_robot\_msgs**

This package contains standard messages which can be used by all other packages.

**Niryo messages**

## *Ned Messages*

Name	Description
<a href="#">CommandStatus</a>	Enum-wise message for status code
<a href="#">ObjectPose</a>	x, y, z, roll, pitch, yaw
<a href="#">RobotState</a>	position, rpy, quaternion
<a href="#">RPY</a>	roll, pitch, yaw
<a href="#">HardwareStatus</a>	several hardware informations
<a href="#">SoftwareVersion</a>	several software version

## Niryo services

## *Ned Services*

Name	Description
<a href="#">GetBool</a>	Return a bool
<a href="#">GetInt</a>	Return a integer
<a href="#">GetStringList</a>	Return a list of string
<a href="#">SetBool</a>	Set a bool and return status
<a href="#">SetInt</a>	Set a integer and return status
<a href="#">SetString</a>	Set a string and return status
<a href="#">Trigger</a>	Trigger a task

## Niryo message dependencies

- [geometry\\_msgs](#)

## Niryo message files

CommandStatus 

```

int32 val

# overall behavior
int32 SUCCESS = 1
int32 CANCELLED = 2
int32 PREEMPTED = 3

int32 FAILURE = -1
int32 ABORTED = -3
int32 STOPPED = -4

int32 ROS_ERROR = -20

int32 FILE_ALREADY_EXISTS = -30

int32 UNKNOWN_COMMAND = -50
int32 NOT_IMPLEMENTED_COMMAND = -51
int32 INVALID_PARAMETERS = -52

# - Hardware
int32 HARDWARE_FAILURE = -110
int32 HARDWARE_NOT_OK = -111
int32 LEARNING_MODE_ON = -112
int32 CALIBRATION_NOT_DONE = -113
int32 DIGITAL_IO_PANEL_ERROR = -114
int32 LED_MANAGER_ERROR = -115
int32 BUTTON_ERROR = -116
int32 WRONG_MOTOR_TYPE = -117
int32 DXL_WRITE_ERROR = -118
int32 DXL_READ_ERROR = -119
int32 CAN_WRITE_ERROR = -120
int32 CAN_READ_ERROR = -121
int32 NO_CONVEYOR_LEFT = -122
int32 NO_CONVEYOR_FOUND = -123
int32 CONVEYOR_ID_INVALID = -124
int32 CALIBRATION_IN_PROGRESS = -125

# - Vision
int32 VIDEO_STREAM_ON_OFF_FAILURE = -170
int32 VIDEO_STREAM_NOT_RUNNING = -171
int32 OBJECT_NOT_FOUND = -172
int32 MARKERS_NOT_FOUND = -173

# - Commander
# Arm Commander
int32 ARM_COMMANDER_FAILURE = -220
int32 GOAL_STILL_ACTIVE = -221
int32 JOG_CONTROLLER_ENABLED = -222
int32 CONTROLLER_PROBLEMS = -223
int32 SHOULD_RESTART = -224
int32 JOG_CONTROLLER_FAILURE = -225

int32 PLAN_FAILED = -230
int32 NO_PLAN_AVAILABLE = -231
int32 INVERT_KINEMATICS_FAILURE = -232

# Tool Commander
int32 TOOL_FAILURE = -251
int32 TOOL_ID_INVALID = -252
int32 TOOL_NOT_CONNECTED = -253
int32 TOOL_ROS_INTERFACE_ERROR = -254

# - Pose Handlers
int32 POSES_HANDLER_CREATION_FAILED = -300
int32 POSES_HANDLER_REMOVAL_FAILED = -301
int32 POSES_HANDLER_READ_FAILURE = -302
int32 POSES_HANDLER_COMPUTE_FAILURE = -303

int32 WORKSPACE_CREATION_FAILED = -308

# - Programs Manager
int32 PROGRAMS_MANAGER_FAILURE = -320
int32 PROGRAMS_MANAGER_READ_FAILURE = -321
int32 PROGRAMS_MANAGER_UNKNOWN_LANGUAGE = -325
int32 PROGRAMS_MANAGER_NOT_RUNNABLE_LANGUAGE = -326
int32 PROGRAMS_MANAGER_EXECUTION_FAILED = -327
int32 PROGRAMS_MANAGER_STOPPING_FAILED = -328
int32 PROGRAMS_MANAGER_AUTORUN_FAILURE = -329
int32 PROGRAMS_MANAGER_WRITING_FAILURE = -330
int32 PROGRAMS_MANAGER_FILE_ALREADY_EXISTS = -331
int32 PROGRAMS_MANAGER_FILE_DOES_NOT_EXIST = -332

# - Serial
int32 SERIAL_FILE_ERROR = -400
int32 SERIAL_UNKNOWN_ERROR = -401

# - System Api Client
int32 SYSTEM_API_CLIENT_UNKNOWN_ERROR = -440
int32 SYSTEM_API_CLIENT_INVALID_ROBOT_NAME = -441
int32 SYSTEM_API_CLIENT_REQUEST_FAILED = -442

```

**ObjectPose**

```

float64 x
float64 y
float64 z

float64 roll
float64 pitch
float64 yaw

```

**RobotState**

```

geometry_msgs/Point position
niryo_robot_msgs/RPY rpy
geometry_msgs/Quaternion orientation

```

**RPY**

```

# roll, pitch and yaw

float64 roll
float64 pitch
float64 yaw

```

**HardwareStatus**

```

std_msgs/Header header
# Raspberry Pi board
int32 rpi_temperature

# Robot version : 1 Ned (only one until needed modification)
int32 hardware_version

# Motors
bool connection_up
string error_message
bool calibration_needed
bool calibration_in_progress

string[] motor_names
string[] motor_types

int32[] temperatures
float64[] voltages
int32[] hardware_errors
string[] hardware_errors_message

```

**SoftwareVersion**

```

string rpi_image_version
string ros_niryo_robot_version

string[] motor_names
string[] stepper_firmware_versions

```

**Niryo Service files****GetBool**

```

---
bool value

```

**GetInt**

```

---
int32 value

```

**GetNameDescriptionList**

```

---
string[] name_list
string[] description_list

```

**GetStringList**

```

---
string[] string_list

```

**SetBool**

```

bool value
---
int32 status
string message

```

**SetInt**

```

int32 value
---
int32 status
string message

```

**SetString**

```

string value
---
int32 status
string message

```

**Trigger**

```

---
int32 status
string message

```

**Niryo\_robot\_modbus****Niryo\_robot\_poses\_handlers**

This package is in charge of dealing with transforms, workspace, grips and trajectories.

**Poses handlers node****Description - Poses handlers**

The ROS Node is made of several services to deal with transforms, workspace, grips and trajectories.

The namespace used is: `/niryo_robot_poses_handlers/`

**Workspaces**

A workspace is defined by 4 markers that form a rectangle. With the help of the robot's calibration tip, the marker positions are learned. The camera returns poses (x, y, yaw) relative to the workspace. We can then infer the absolute object pose in robot coordinates.

**Grips**

When we know the object pose in robot coordinates, we can't directly send this pose to the robot because we specify the target pose of the tool\_link and not of the actual TCP (tool center point). Therefore we introduce the notion of grip. Each end effector has its own grip that specifies where to place the robot with respect to the object. Currently, the notion of grip is not part of the python/tcp/blockly interface because it would add an extra layer of complexity that is not really necessary for the moment. Therefore we have a default grip for all tools that is selected automatically based on the current tool id. However, everything is ready if you want to define custom grips, e.g. for custom tools or for custom grip positions.

## The vision pick loop

1. Camera detects object relative to markers and sends `x_rel`, `y_rel`, `yaw_rel`.
2. The object is placed on the workspace, revealing the object pose in robot coordinates `x`, `y`, `z`, roll, pitch, yaw.
3. The grip is applied on the absolute object pose and gives the pose the robot should move to.

## Poses & trajectories

List of poses

## Parameters - Poses handlers

*Poses Handlers' Parameters*

Name	Description
<code>workspace_dir</code>	Path to the Workspace storage mother folder
<code>grip_dir</code>	Path to the Grip storage mother folder
<code>poses_dir</code>	Path to the Poses storage mother folder
<code>trajectories_dir</code>	Path to the Trajectory storage mother folder

## Services - Poses handlers

*Poses Handlers' Services*

Name	Message Type	Description
<code>manage_workspace</code>	<code>ManageWorkspace</code>	Save/Delete a workspace
<code>get_workspace_ratio</code>	<code>GetWorkspaceRatio</code>	Get ratio of a workspace
<code>get_workspace_list</code>	<code>GetNameDescriptionList</code>	Get list of workspaces name & description
<code>get_workspace_poses</code>	<code>GetWorkspaceRobotPoses</code>	Get workspace's robot poses
<code>get_target_pose</code>	<code>GetTargetPose</code>	Get saved programs name
<code>manage_pose</code>	<code>ManagePose</code>	Save/Delete a Pose
<code>get_pose</code>	<code>GetPose</code>	Get Pose
<code>get_pose_list</code>	<code>GetNameDescriptionList</code>	Get list of poses name & description
<code>manage_trajectory</code>	<code>ManageTrajectory</code>	Save/Delete a Trajectory
<code>get_trajectory</code>	<code>GetTrajectory</code>	Get Trajectory
<code>get_trajectory_list</code>	<code>GetNameDescriptionList</code>	Get list of trajectories name & description

All these services are available as soon as the node is started.

## Dependencies - Poses handlers

- `geometry_msgs`
- `moveit_msgs`
- `Niryo_robot_msgs`
- `tf`

## Services & messages files - Poses handlers

### GetPose (Service)

```
string name
...
int32 status
string message
niryo_robot_poses_handlers/NiryoPose pose
```

### GetTargetPose (Service)

```
string workspace
float32 height_offset
float32 x_rel
float32 y_rel
float32 yaw_rel
...
int32 status
string message
niryo_robot_msgs/RobotState target_pose
```

### GetTrajectory (Service)

```
string name
...
int32 status
string message
geometry_msgs/Pose[] list_poses
```

### GetWorkspaceRatio (Service)

```
string workspace
---
int32 status
string message
float32 ratio # width/height
```

**GetWorkspaceRobotPoses (Service)¶**

```
string name
---
int32 status
string message
niryo_robot_msgs/RobotState[] poses
```

**ManagePose (Service)¶**

```
int32 cmd
int32 SAVE = 1
int32 DELETE = -1

niryo_robot_poses_handlers/NiryoPose pose
---
int32 status
string message
```

**ManageTrajectory (Service)¶**

```
int32 cmd
int32 SAVE = 1
int32 DELETE = -1

string name
string description

geometry_msgs/Pose[] poses
---
int32 status
string message
```

**ManageWorkspace (Service)¶**

```
int32 SAVE = 1
int32 SAVE_WITH_POINTS = 2
int32 DELETE = -1

int32 cmd

niryo_robot_poses_handlers/Workspace workspace
---

int32 status
string message
```

**NiryoPose (Message)¶**

```
string name
string description

float64[] joints
geometry_msgs/Point position
niryo_robot_msgs/RPY rpy
geometry_msgs/Quaternion orientation
```

**Trajectory (Message)¶**

```
string name
geometry_msgs/Pose[] poses
```

**Workspace (Message)¶**

```
string name # maximum lenght of workspace's name is 30 characters
string description

geometry_msgs/Point[] points
niryo_robot_msgs/RobotState[] poses
```

**Niryo\_robot\_programs\_manager**

This package is in charge of interpreting/running/saving programs. It is used by Niryo Studio.

**Programs manager node**

The ROS Node is made of several services to deal with the storage and running of programs.

Calls are not available from the Python ROS Wrapper, as it made to run its programs with the Python ROS Wrapper.

The namespace used is: `/niryo_robot_programs_manager/`

**Parameters - Programs manager¶**

*Programs Manager's Parameters*

Name	Description
<code>autorun_file_name</code>	Name of the file containing auto run infos
<code>programs_dir</code>	Path to the Program storage mother folder

**Services - Programs manager¶**

*Programs manager Services*

Name	Message Type	Description
execute_program	ExecuteProgram	Executes a program
execute_program_autorun	Trigger	Executes autorun program
get_program	GetProgram	Retrieves saved program
get_program_autorun_infos	GetProgramAutorunInfos	Gets autorun settings
get_program_list	GetProgramList	Gets saved programs' name
manage_program	ManageProgram	Saves and Deletes programs
set_program_autorun	SetProgramAutorun	Sets autorun settings
stop_program	Trigger	Stops the current running program

All these services are available as soon as the node is started whereas on standalone mode or not.

## Dependencies - Programs manager

- Niryo\_robot\_msgs
- python-yaml
- std\_msgs

## Services &amp; messages files - Programs manager

## ExecuteProgram (Service)

```
bool execute_from_string
string name
string code_string
niryo_robot_programs_manager/ProgramLanguage language
---
int16 status
string message
```

## GetProgram (Service)

```
string name
niryo_robot_programs_manager/ProgramLanguage language
---
int32 status
string message
string code
string description
```

## GetProgramAutorunInfos (Service)

```
---
int32 status
string message
niryo_robot_programs_manager/ProgramLanguage language
string name
# Mode
int8 ONE_SHOT = 1
int8 LOOP = 2
int8 mode
```

## GetProgramList (Service)

```
niryo_robot_programs_manager/ProgramLanguage language
---
string[] programs_names
niryo_robot_programs_manager/ProgramLanguageList[] list_of_language_list
string[] programs_description
```

## ManageProgram (Service)

```
# Command
int32 SAVE = 1
int32 DELETE = -1
int8 cmd

# Program Name
string name

# - Creation
niryo_robot_programs_manager/ProgramLanguage language

string code
string description

bool allow_overwrite
---
int16 status
string message
```

## SetProgramAutorun (Service)

```
# Program language
niryo_robot_programs_manager/ProgramLanguage language

# Program Name
string name

# Mode
int8 DISABLE = 0
int8 ONE_SHOT = 1
int8 LOOP = 2

int8 mode

---
int16 status
string message
```

**ProgramLanguage (Message)¶**

```
int8 NONE = -1
int8 ALL = 0

# Runnable
int8 PYTHON2 = 1
int8 PYTHON3 = 2

# Not Runnable
int8 BLOCKLY = 66

int8 used
```

**ProgramLanguageList (Message)¶**

```
niryo_robot_programs_manager/ProgramLanguage[] language_list
```

**ProgramList (Message)¶**

```
string[] programs_names
niryo_robot_programs_manager/ProgramLanguageList[] list_of_language_list
string[] programs_description
```

**Niryo\_robot\_rpi**

This package deals with Raspberry Pi related stuff (Button, fans, I/O, leds, ...).

**Raspberry Pi Node**

The ROS Node manages the following components:

- Physical top button: executes actions when the button is pressed.
- Digital I/O panel: gets commands and sends the current state of digital I/Os. Also controls tools like the Electromagnet.
- Robot fans.
- Led: sets the LED color.
- ROS log: can remove all previous logs on startup to prevent a lack of disk space in the long run (SD cards do not have infinite storage).

The namespace used is: `/niryo_robot_rpi/`

**Note that this package should not be used when you are using Ned ROS stack on your computer in simulation mode. Executes actions when the button is pressed.**

**Publisher - Raspberry Pi¶**
*RPI Package's Publishers*

Name	Message Type	Description
<code>pause_state</code>	<code>PausePlanExecution</code>	Publish the current execution state launched when button is pressed
<code>/niryo_robot/blockly/save_current_point</code>	<code>std_msgs/Int32</code>	Publish current point when user is in Blockly page to save block by pressing button
<code>/niryo_robot/rpi/is_button_pressed</code>	<code>std_msgs/Bool</code>	Publish the button state (true if pressed)
<code>digital_io_state</code>	<code>DigitalIOState</code>	Publish the I/Os state by giving for each its pin / name / mode / state
<code>/niryo_robot/rpi/led_state</code>	<code>std_msgs/Int8</code>	Publish the current led color
<code>ros_log_status</code>	<code>LogStatus</code>	Publish the current log status (log size / available disk / boolean if should delete ros log on startup)

**Services - Raspberry Pi¶**
*RPI Services*

Name	Message Type	Description
<code>shutdown_rpi</code>	<code>SetInt</code>	Shutdown the Raspberry Pi
<code>/niryo_robot/rpi/change_button_mode</code>	<code>SetInt</code>	Change top button mode (autorun program, blockly, nothing, ...)
<code>get_digital_io</code>	<code>GetDigitalIO</code>	Get digital IO state list
<code>set_digital_io_mode</code>	<code>SetDigitalIO</code>	Set a digital IO to the mode given
<code>set_digital_io_state</code>	<code>SetDigitalIO</code>	Set a digital IO to the state given
<code>set_led_state</code>	<code>std_msgs/SetInt</code>	Set led state
<code>set_led_custom_blinker</code>	<code>LedBlinker</code>	Set the led in blink mode with the color given
<code>purge_ros_logs</code>	<code>SetInt</code>	Purge ROS log
<code>set_purge_ros_log_on_startup</code>	<code>SetInt</code>	Modify the permanent settings that tells if robot should purge its ROS log at each boot

Dependencies - Raspberry Pi

- std\_msgs
- actionlib\_msgs
- sensor\_msgs
- Niryo\_robot\_msgs
- Niryo\_robot\_arm\_commander

#### Services & Messages files - Raspberry Pi

##### ChangeMotorConfig (Service)

```
int32[] can_required_motor_id_list
int32[] dxl_required_motor_id_list
...
int32 status
string message
```

##### GetDigitalIO (Service)

```
int32 pin
...
int32 status
string message

int32 pin
string name
int32 mode
int32 state
```

##### LedBlinker (Service)

```
uint8 LED_OFF = 0
uint8 LED_BLUE = 1
uint8 LED_GREEN = 2
uint8 LED_BLUE_GREEN = 3
uint8 LED_RED = 4
uint8 LED_PURPLE = 5
uint8 LED_RED_GREEN = 6
uint8 LED_WHITE = 7

bool activate
uint8 frequency # between 1hz and 100Hz
uint8 color
float32 blinker_duration # 0 for infinite
...
int32 status
string message
```

##### SetDigitalIO (Service)

```
uint8 pin
uint8 value
...
int32 status
string message
```

##### DigitalIOState (Service)

```
# GPIO pin
int32[] pins
# PIN names seen by user to make it simpler
string[] names
# IN/OUT
int32[] modes
# HIGH/LOW
int32[] states
```

##### LogStatus (Service)

```
std_msgs/Header header

# in MB
int32 log_size
int32 available_disk_size
bool purge_log_on_startup
```

#### Niryo robot tools commander package

Provides functionalities to control end-effectors and accessories for Ned.

This package allows to manage the TCP (Tool Center Point) of the robot. If the functionality is activated, all the movements (in Cartesian coordinates [x, y, z, roll, pitch, yaw]) of the robot will be performed according to this TCP. The same program can then work with several tools by adapting the TCP transformation to them. By default this feature is disabled, but can be enabled through the robot services.

##### Tools Commander node

The ROS Node is made of services to equip tool, an action server for tool command and topics for the current tool or the tool state.

The namespace used is: `/niryo_robot_tools_commander/`

##### Action server - tools

*Tools Package Action Server*

Name	Message Type	Description
action_server	ToolAction	Command the tool through an action server

##### Publisher - tools

## Tools Package Publishers

Name	Message Type	Description
current_id	std_msgs/Int32	Publish the current tool ID
tcp	TCP	Publish if the TCP (Tool Center Point) is enabled and transformation between the tool_link and the TCP

## Services - tools¶

## Tools Package Services

Name	Message Type	Description
update_tool	std_srvs/Trigger	Ping/scan for a dxl motor flashed with an ID corresponding to a tool and equip it (if found)
equip_electromagnet	SetInt	Equip the electromagnet with the motor ID given as parameter
enable_tcp	SetBool	Enable or disable the TCP (Tool Center Point) functionality. When we activate it, the transformation will be the last one saved since the robot started. By default it will be the one of the equipped tool.
set_tcp	SetTCP	Activate the TCP (Tool Center Point) functionality and defines a new TCP transformation.
reset_tcp	std_srvs/Trigger	Reset the TCP transformation. By default it will be the one of the equipped tool.

## Dependencies - tools¶

- Niryo\_robot\_msgs
- std\_msgs
- geometry\_msgs

## Action files - tools

## ToolAction (Action)¶

```
# goal
niryo_robot_tools_commander/ToolCommand cmd
...
# result
int32 status
string message
...
# feedback
int32 progression
```

## Messages files - tools

## ToolCommand (Message)¶

```
# Gripper
uint8 OPEN_GRIPPER = 1
uint8 CLOSE_GRIPPER = 2

# Vacuum pump
uint8 PULL_AIR_VACUUM_PUMP = 10
uint8 PUSH_AIR_VACUUM_PUMP = 11

# Tools controlled by digital I/Os
uint8 SETUP_DIGITAL_IO = 20
uint8 ACTIVATE_DIGITAL_IO = 21
uint8 DEACTIVATE_DIGITAL_IO = 22

uint8 cmd_type

# Gripper1= 11, Gripper2=12, Gripper3=13, VacuumPump=31, Electromagnet=30
uint8 tool_id

# if gripper close
uint16 gripper_close_speed

# if gripper open
uint16 gripper_open_speed

# if vacuum pump or electromagnet grove
bool activate

# if tool is set by digital outputs (electromagnet)
# if gpio<0 get value in memory
int8 gpio
```

## TCP (Message)¶

```
bool enabled

geometry_msgs/Point position
niryo_robot_msgs/RPY rpy
geometry_msgs/Quaternion orientation
```

## Services files - tools

## SetTCP (Service)¶

```
geometry_msgs/Point position

#Only one of the two is required.
#if both are filled, the quaternion will be chosen by default
niryo_robot_msgs/RPY rpy
geometry_msgs/Quaternion orientation
...
int32 status
string message
```

## Niryo\_robot\_unit\_tests

This package provides config and launches files to start Ned unit tests.

## Currently, tested packages are:

- [Niryo\\_robot\\_poses\\_handlers](#)
- [Niryo\\_robot\\_programs\\_manager](#)
- [Python ROS Wrapper](#)

## Run tests with launch file

Using a launch file make the task easy as you do not have to bother about your current directory. Nevertheless, it starts a roscore, and so, the process won't end as "/rosout" will be still alive. To overcome this issue, the test node is set as **required**, so do not worry if you see some red lines at the execution's end.

### Launch file with Rviz

It will start the ROS Stack and Rviz, then run the tests:

```
roslaunch niryo_robot_unit_tests simulation_rviz_tests.launch
```

### Launch file with Gazebo

It will start the ROS Stack and Gazebo, then run the tests. Gazebo add Physics, Gripper and Camera to the tests:

```
roslaunch niryo_robot_unit_tests simulation_gazebo_tests.launch
```

### Launch file Headless

It will start the ROS Stack then run the tests in a headless version (no display). It is very useful for CI purposes:

```
roslaunch niryo_robot_unit_tests simulation_headless_tests.launch
```

## Run tests With Python file

Using a python file implies to use the correct relative directory, but it also allows to get status code at the end of the script.

### Python file with Rviz

```
python niryo_robot_unit_tests/scripts/script_test_rviz.py
```

### Python file with Gazebo

```
python niryo_robot_unit_tests/scripts/script_test_gazebo.py
```

### Python file Headless

```
python niryo_robot_unit_tests/scripts/script_test_headless.py
```

## Niryo\_robot\_user\_interface

This packages handle high-level user interface commands coming TCP requests and also system-related features like I/Os, LED and fans.

### You can find their documentations here:

- [TCP Server](#)

### Use Ned's TCP server

Ned is permanently running a TCP Server to acquire requests. This server is built on top of the [Ned Python ROS Wrapper](#) ([index.html#document-source/ros\\_wrapper](#)).

It offers a simple way for developers to create programs for robot to control them via remote communication on a computer, on a mobile or any device with network facilities.

Programs can communicate through network TCP with the robots in any language available.

### Connection

To access the server, you will have to use to robot's IP address and communicate via the **port 40001**.

### Communication

Only one client can communicate with the server (reconnection effective but no multi clients).

The server answers only after the command is done, so it can't deal with multiple commands at the same time.

### Packet convention

#### General format

For easier usage and easier debugging, the communication is based on JSON format.

Every package have this following shape: <json\_packet\_size>{<json\_content>}<payload> .

The JSON packet size is an unsigned short coded on 2 bytes.

The JSON contains command's name & params.

Payload contains *heavy* data like an image.

### Request

#### Format - Request

As no function requests a payload in input, requests have the following.

Format: <json\_packet\_size>{'param\_list': [<param\_1>, <param\_2>, ...], 'command': <command\_str>}

Examples - Request¶

Calibrate auto: {'param\_list': ['AUTO'], 'command': 'CALIBRATE'}

Move joints: {'param\_list': [0.0, 0.0, 0.0, 0.0, 0.0, 0.0], 'command': 'MOVE\_JOINTS'}

Answer¶

Format - Answer¶

Firstly, answers indicate to the user if its command has been well executed. This is indicated in the JSON by the parameter "status".

A successful answer will have the format:

```
{'status': 'OK', 'list_ret_param': [<param_1>, <param_2>, ...], 'payload_size': <payload_size_int>, 'command': <command_str><payload_str>}
```

An unsuccessful answer will have the format: {'status': 'KO', 'message': <message\_str>}

Examples - Answer¶

Calibrate Auto: {'status': 'OK', 'list\_ret\_param': [], 'payload\_size': 0, 'command': 'CALIBRATE'}

Get Pose: {"status": "OK", "list\_ret\_param": [0.2, 0.15, 0.35, 0.5, -0.6, 0.1], "payload\_size": 0, "command": "GET\_POSE"}

Commands enum for TCP server¶

`class CommandEnum`

Enumeration of all commands used

`CALIBRATE= 0`

`SET_LEARNING_MODE= 1`

`GET_LEARNING_MODE= 2`

`SET_ARM_MAX_VELOCITY= 3`

`SET_JOG_CONTROL= 4`

`GET_JOINTS= 10`

`GET_POSE= 11`

`GET_POSE_QUAT= 12`

`MOVE_JOINTS= 20`

`MOVE_POSE= 21`

`SHIFT_POSE= 22`

`MOVE_LINEAR_POSE= 23`

`SHIFT_LINEAR_POSE= 24`

`JOG_JOINTS= 25`

`JOG_POSE= 26`

`FORWARD_KINEMATICS= 27`

`INVERSE_KINEMATICS= 28`

`GET_POSE_SAVED= 50`

`SAVE_POSE= 51`

`DELETE_POSE= 52`

`GET_SAVED_POSE_LIST= 53`

`PICK_FROM_POSE= 60`

`PLACE_FROM_POSE= 61`

`PICK_AND_PLACE= 62`

GET\_TRAJECTORY\_SAVED= 80

EXECUTE\_TRAJECTORY\_FROM\_POSES= 81

EXECUTE\_TRAJECTORY\_SAVED= 82

SAVE\_TRAJECTORY= 83

DELETE\_TRAJECTORY= 84

GET\_SAVED\_TRAJECTORY\_LIST= 85

EXECUTE\_TRAJECTORY\_FROM\_POSES\_AND\_JOINTS= 86

UPDATE\_TOOL= 120

OPEN\_GRIPPER= 121

CLOSE\_GRIPPER= 122

PULL\_AIR\_VACUUM\_PUMP= 123

PUSH\_AIR\_VACUUM\_PUMP= 124

SETUP\_ELECTROMAGNET= 125

ACTIVATE\_ELECTROMAGNET= 126

DEACTIVATE\_ELECTROMAGNET= 127

GET\_CURRENT\_TOOL\_ID= 128

GRASP\_WITH\_TOOL= 129

RELEASE\_WITH\_TOOL= 130

ENABLE\_TCP= 140

SET\_TCP= 141

RESET\_TCP= 142

TOOL\_REBOOT= 145

SET\_PIN\_MODE= 150

DIGITAL\_WRITE= 151

DIGITAL\_READ= 152

GET\_DIGITAL\_IO\_STATE= 153

GET\_HARDWARE\_STATUS= 154

SET\_CONVEYOR= 180

UNSET\_CONVEYOR= 181

CONTROL\_CONVEYOR= 182

GET\_CONNECTED\_CONVEYORS\_ID= 183

GET\_IMAGE\_COMPRESSED= 200

GET\_TARGET\_POSE\_FROM\_REL= 201

GET\_TARGET\_POSE\_FROM\_CAM= 202

VISION\_PICK= 203

MOVE\_TO\_OBJECT= 205

DETECT_OBJECT= 204
GET_CAMERA_INTRINSICS= 210
SAVE_WORKSPACE_FROM_POSES= 220
SAVE_WORKSPACE_FROM_POINTS= 221
DELETE_WORKSPACE= 222
GET_WORKSPACE_RATIO= 223
GET_WORKSPACE_LIST= 224
SET_IMAGE_BRIGHTNESS= 230
SET_IMAGE_CONTRAST= 231
SET_IMAGE_SATURATION= 232
GET_IMAGE_PARAMETERS= 235

## Niryo\_robot\_vision

This package is the one dealing with all vision related stuff.

### Vision Node

The ROS Node is made of several services to deal with video streaming, object detection... The node is working exactly the same way if you chose to use it on simulation or reality.

This node can be launched locally in a standalone mode via the command:

```
roslaunch niryo_robot_vision vision_node_local.launch
```

Configuration (Frame Per Second, Camera Port, Video Resolution) can be edited in the config file:

- For "standard" Node: *niryo\_robot\_vision/config/video\_server\_setup.yaml*
- For local Node: *niryo\_robot\_vision/config/video\_server\_setup\_local.yaml*

The namespace used is: */niryo\_robot\_vision/*

#### Parameters - Vision

*Vision Package's Parameters*

Name	Description
<code>frame_rate</code>	Streams frame rate
<code>simulation_mode</code>	Sets to true if you are using the gazebo simulation. It will adapt how the node get its video stream
<code>debug_compression_quality</code>	Debugs Stream compression quality
<code>stream_compression_quality</code>	Streams compression quality
<code>subsampling</code>	Streams subsampling factor

#### Publisher - Vision

*Vision Package's Publishers*

Name	Message Type	Description
<code>compressed_video_stream</code>	<code>sensor_msgs/CompressedImage</code>	Publishes the last image read as a compressed image
<code>video_stream_parameters</code>	<code>ImageParameters</code>	Publishes the brightness, contrast and saturation settings of the video stream

#### Services - Vision

## Programs manager Services

Name	Message Type	Description
debug_colors	DebugColorDetection	Returns an annotated image to emphasize what happened with color detection
debug_markers	DebugMarkers	Returns an annotated image to emphasize what happened with markers detection
obj_detection_rel	ObjDetection	Object detection service
start_stop_video_streaming	SetBool	Starts or stops video streaming
take_picture	TakePicture	Saves a picture in the specified folder
set_brightness	SetImageParameter	Sets the brightness of the video stream
set_contrast	SetImageParameter	Sets the contrast of the video stream
set_saturation	SetImageParameter	Sets the saturation of the video stream

All these services are available as soon as the node is started.

## Dependencies - Vision

- Niryo\_robot\_msgs
- sensor\_msgs

## Topics files - Vision

## ImageParameters (Topic)

```
float64 brightness_factor
float64 contrast_factor
float64 saturation_factor
```

## Services files - Vision

## DebugColorDetection (Service)

```
string color
...
sensor_msgs/CompressedImage img
```

## DebugMarkers (Service)

```
...
bool markers_detected
sensor_msgs/CompressedImage img
```

## ObjDetection (Service)

```
string obj_type
string obj_color
float32 workspace_ratio
bool ret_image
...
int32 status

niryo_robot_msgs/ObjectPose obj_pose

string obj_type
string obj_color

sensor_msgs/CompressedImage img
```

## TakePicture (Service)

```
string path
...
bool success
```

## SetImageParameter (Service)

```
float64 factor
...
int32 status
string message
```

## Hardware ROS Stack documentation



ROS (Robot Operating System) is an Open-Source Robotic Framework which allows to ease robot software development. The framework is used in almost each part of Ned's software.

The high-level packages (motion planner, vision, ...) are coded in Python to give good readability whereas communication with Hardware is developed in C++ to ensure speed.


**Note**

To learn more about ROS, go on [Official ROS Wiki](http://wiki.ros.org/) (<http://wiki.ros.org/>).

In this section, you will have access to all information about each Niryo robot's ROS hardware stack packages.

## Niryo\_robot\_hardware\_interface

This package handles packages related to the robot's hardware.  
It launches hardware interface nodes, motors communication and driver.



Global overview of hardware stack packages organization.

### Hardware interface node

The ROS Node is made to launch hardware interface and communication:

- Conveyor Belt interface
- Joints interface
- Tools interface
- Fake interface
- Dynamixel driver
- Stepper driver

The namespace used is: `/niryo_robot_hardware_interface/`.

#### Parameters - hardware interface

*Hardware Interface's Parameters*

Name	Description
<code>publish_hw_status_frequency</code>	Publishes rate for hardware status. Default : 2.0'
<code>publish_software_version_frequency</code>	Publishes rate for software status. Default : 2.0'

#### Published topics - hardware interface

*Hardware Interface's Published Topics*

Name	Message Type	Description
<code>hardware_status</code>	<code>niryo_robot_msgs/HardwareStatus</code>	Motors, bus, joints and CPU status
<code>software_version</code>	<code>niryo_robot_msgs/SoftwareVersion</code>	Raspberry and stepper software version

#### Services - hardware interface

*Hardware Interface Package Services*

Name	Message Type	Description
<code>launch_motors_report</code>	Trigger	Starts motors report
<code>reboot_motors</code>	Trigger	Reboots motors
<code>stop_motors_report</code>	Trigger	Stops motors report

#### Dependencies - hardware interface

- `Tools_interface`
- `Joints_interface`
- `Conveyor_interface`
- `CPU_interface`
- `Fake_interface`
- `Niryo_robot_msgs`

## Joints\_interface

This package handles packages related to the robot's joints controller.  
It provides an interface to `ros_control`.

**Joints interface node****The ROS Node is made to:**

- Interface robot's motors to joint trajectory controller, from `ros_control` package.
- Create a controller manager, from `controller_manager` package, provides the infrastructure to load, unload, start and stop controllers.
- Interface with motors calibration.
- Initialize motors parameters.

Parameters - joints interface¶

*Joint Interface's Parameters*

Name	Description
<code>ros_control_loop_frequency</code>	Controls loop frequency. Default: '100.0'
<code>publish_learning_mode_frequency</code>	Publishes rate for learning mode state. Default: '2.0'
<code>calibration_timeout</code>	Waiting time between 2 commands during the calibration process. Default: '30'
<code>calibration_file</code>	File directory where is saved motors calibration value. Default: '/home/niryo/niryo_robot_saved_files/stepper_motor_calibration_offsets.txt'

Published topics - joints interface¶

*Joint Interface's Published Topics*

Name	Message Type	Description
<code>/niryo_robot/learning_mode/state</code>	<code>std_msgs/Bool</code>	Learning mode state

Services - joints interface¶

*Joint Interface Package Services*

Name	Message Type	Description
<code>/niryo_robot/joints_interface/calibrate_motors</code>	<code>SetInt</code>	Starts motors calibration - value can be 1 for auto calibration, 2 for manual
<code>/niryo_robot/joints_interface/request_new_calibration</code>	<code>Trigger</code>	Unsets motors calibration
<code>niryo_robot/learning_mode/activate</code>	<code>Trigger</code>	Either activates or deactivates learning mode

Dependencies - joints interface¶

- `hardware_interface`
- `controller_manager`
- `Dynamixel_driver`
- `Stepper_driver`
- `Niryo_robot_msgs`

**Conveyor\_interface**

This package handles Niryo's Conveyor Belt.  
You can control two Conveyors Belt.

**Conveyor Belt interface node****The ROS Node is made to:**

- Initialize Conveyor Belt motor parameters.
- Set and control Conveyors Belt.
- Publish Conveyors Belt state.

The namespace used is: `/niryo_robot/conveyor/`

Parameters - Conveyor Belt interface¶

*Conveyor Belt Interface's Parameters*

Name	Description
<code>publish_frequency</code>	Publishes rate for conveyors state. Default: '2.0'

Published topics - Conveyor Belt interface¶

*Conveyor Belt Interface's Published Topics*

Name	Message Type	Description
<code>feedback</code>	<code>ConveyorFeedbackArray</code>	Conveyors Belt states

Services - Conveyor Belt interface¶

*Conveyor Belt Interface Package Services*

Name	Message Type	Description
control_conveyor	ControlConveyor	Sends a command to the desired Conveyor Belt
ping_and_set_conveyor	SetConveyor	Scans and sets a new Conveyor Belt

Dependencies - Conveyor Belt interface¶

- std\_msgs
- stepper\_driver

**Services & messages files - Conveyor Belt interface**

ControlConveyor (Service)¶

```
uint8 id
bool control_on
int16 speed
int8 direction
...
int16 status
string message
```

SetConveyor (Service)¶

```
uint8 cmd
uint8 id

uint8 ADD = 1
uint8 REMOVE = 2

...
int16 id
int16 status
string message
```

ConveyorFeedbackArray (Message)¶

```
conveyor_interface/ConveyorFeedback[] conveyors
```

ConveyorFeedback (Message)¶

```
#Conveyor id ( either 12 or 18 )
uint8 conveyor_id
#Conveyor Connection state ( if it is enabled )
bool connection_state
# Conveyor Controls state : ON or OFF
bool running
# Conveyor Speed ( 1-> 100 % )
int16 speed
# Conveyor direction ( backward or forward )
int8 direction
```

**Tools\_interface**

This package handles Niryo's tools.

**Tools interface node****The ROS Node is made to:**

- Set and control tools.
- Publish tool connection state.

Parameters - tools interface¶

*Tools Interface's Parameters*

Name	Description
check_tool_connection_frequency	Publishes and controls rate for tools connection state. Default: '2.0'
id_list	List of tools id Default: '[11,12,13,14,31]'
motor_type_list	List of motor tools type Default: '["xl320","xl320","xl320","xl320","xl320"]'

Published topics - tools interface¶

*Tools Interface's Published Topics*

Name	Message Type	Description
/niryo_robot_hardware/tools/current_id	std_msgs/Int32	Current tool ID

Services - tools interface¶

## Tool/Interface Package Services

Name	Message Type	Description
niryo_robot/tools/ping_and_set_dxl_tool	tools_interface/PingDxlTool	Scans and sets for a tool plugged
niryo_robot/tools/open_gripper	tools_interface/OpenGripper	Opens a gripper tool
niryo_robot/tools/close_gripper	tools_interface/CloseGripper	Closes a gripper tool
niryo_robot/tools/pull_air_vacuum_pump	tools_interface/PullAirVacuumPump	Pulls vacuum pump tool
niryo_robot/tools/push_air_vacuum_pump	tools_interface/PushAirVacuumPump	Pushes vacuum pump tool
niryo_robot/tools/reboot	std_srvs/Trigger	Reboots the motor of the equipped tool

## Dependencies - tools interface¶

- std\_msgs
- std\_srvs
- Dynamixel\_driver

## Services &amp; messages files - tools interface

## PingDxlTool (Service)¶

```
---
```

```
int32 state
uint8 id
```

## OpenGripper (Service)¶

```
uint8 id

int16 open_position
int16 open_speed
int16 open_hold_torque
---
uint8 state
```

## CloseGripper (Service)¶

```
uint8 id

int16 close_position
int16 close_speed
int16 close_hold_torque
int16 close_max_torque
---
uint8 state
```

## PullAirVacuumPump (Service)¶

```
uint8 id

int16 pull_air_position
int16 pull_air_hold_torque
---
uint8 state
```

## PushAirVacuumPump (Service)¶

```
uint8 id

int16 push_air_position
---
uint8 state
```

## CPU\_interface

This package handles CPU states.

## CPU interface node

## The ROS Node is made to:

- monitor CPU temperature.

## Parameters - CPU interface¶

## CPU Interface's Parameters

Name	Description
read_rpi_diagnostics_frequency	Publishes rate for CPU temperature Default: '0.25'
temperature_warn_threshold	CPU temperature [celsius] threshold before a warn message Default: '75'
temperature_shutdown_threshold	CPU temperature [celsius] threshold before shutdown the robot Default: '85'

## Fake\_interface

This package provides fakes hardware interfaces when the robot is used in simulation.

## Fake interface node

**The ROS Node is made to simulate:**

- tools interface
- Conveyor Belt interface
- joints interface

Published topics - fake interface¶

*Fake Interface's Published Topics*

Name	Message Type	Description
/niryo_robot/learning_mode/state	std_msgs/Bool	Learning mode state
/niryo_robot_hardware/tools/current_id	std_msgs/Int32	Current tool ID

Services - fake interface¶

*Fake Interface Package Services*

Name	Message Type	Description
/niryo_robot/joints_interface/calibrate_motors	SetInt	Starts motors calibration - value can be 1 for auto calibration, 2 for manual
/niryo_robot/joints_interface/request_new_calibration	Trigger	Unsets motors calibration
niryo_robot/learning_mode/activate	Trigger	Either activates or deactivates learning mode
niryo_robot/tools/ping_and_set_dxl_tool	tools_interface/PingDxlTool	Scans and sets for a tool plugged
niryo_robot/tools/open_gripper	tools_interface/OpenGripper	Opens a Gripper tool
niryo_robot/tools/close_gripper	tools_interface/CloseGripper	Closes a Gripper tool
niryo_robot/tools/pull_air_vacuum_pump	tools_interface/PullAirVacuumPump	Pulls Vacuum Pump tool
niryo_robot/tools/push_air_vacuum_pump	tools_interface/PushAirVacuumPump	Pushes Vacuum Pump tool
/niryo_robot/conveyor/control_conveyor	ControlConveyor	Sends a command to the desired Conveyor Belt
/niryo_robot/conveyor/ping_and_set_conveyor	SetConveyor	Scans and sets a new Conveyor Belt

Dependencies - fake interface¶

- std\_msgs
- hardware\_interface
- controller\_manager
- Niryo\_robot\_msgs
- Tools\_interface
- Joints\_interface
- Conveyor\_interface

**Dynamixel\_driver**

This package handles dynamixel motors communication through dynamixel sdk.  
It provides an interface to [dynamixel\\_sdk](#).

**Dynamixel Driver Node****The ROS Node is made to:**

- Send commands to dynamixel motors
- Receive dynamixel motors data

Parameters - Dynamixel Driver¶

*Dynamixel Driver's Parameters*

Name	Description
dxl_hardware_control_loop_frequency	Controls loop frequency. Default: '100.0'
dxl_hardware_write_frequency	Writes frequency. Default: '50.0'
dxl_hardware_read_data_frequency	Reads frequency. Default: '15.0'
dxl_hardware_read_status_frequency	Reads dynamixels status frequency. Default: '0.5'
dxl_hardware_check_connection_frequency	Checks dynamixels connection frequency. Default: '2.0'
dxl_motor_id_list	List of dynamixels ID Default: '[2,3,6]'
dxl_motor_type_list	List of dynamixels type Default: "[xl430","xl430","xl320"]'

Services - Dynamixel Driver¶

*Dynamixel Driver Package Services*

Name	Message Type	Description
niryo_robot/dynamixel_driver/set_dx1_leds	SetInt	Controls dynamixel LED
niryo_robot/dynamixel_driver/send_custom_dx1_value	dynamixel_driver/SendCustomDx1Value	Sends a custom dynamixel command

## Dependencies - Dynamixel Driver

- [dynamixel\\_sdk](#)
- [Niryo\\_robot\\_msgs](#)

## Services &amp; Messages files - Dynamixel Driver

## SendCustomDx1Value (Service)

```
# Check XL-320 and XL-430 reference doc for
# the complete register table

int8 motor_type # 3 (XL-320) or 2 (XL-430)
uint8 id
int32 value
int32 reg_address
int32 byte_number
...
int32 status
string message
```

## DxlMotorHardwareStatus (Message)

```
niryo_robot_msgs/MotorHeader motor_identity

uint32 temperature
float64 voltage
uint32 error
string error_msg
```

## DxlMotorCommand (Message)

```
uint8 cmd_type
uint8 CMD_TYPE_POSITION=1
uint8 CMD_TYPE_VELOCITY=2
uint8 CMD_TYPE_EFFORT=3
uint8 CMD_TYPE_TORQUE=4

uint8[] motors_id
uint32[] params
```

## DxlArrayMotorHardwareStatus (Message)

```
std_msgs/Header header
dynamixel_driver/DxlMotorHardwareStatus[] motors_hw_status
```

## Stepper\_driver

This package handles stepper motors communication.

## Stepper driver node

## The ROS Node is made to:

- Send commands to stepper motors.
- Receive stepper motors data.

## Parameters - stepper driver

*Stepper Driver's Parameters*

Name	Description
can.hardware_control_loop_frequency	Controls loop frequency. Default: '1500.0'
can_hw_write_frequency	Writes frequency. Default: '50.0'
can_hw_check_connection_frequency	Checks steppers connection frequency. Default: '2.0'
stepper_motor_id_list	List of steppers ID Default: '[1,2,3]'

## Dependencies - stepper driver

- [mcp\\_can\\_rpi](#)
- [Niryo\\_robot\\_msgs](#)

## Services &amp; messages files - stepper driver

## StepperMotorHardwareStatus (Message)

```
niryo_robot_msgs/MotorHeader motor_identity

string firmware_version
int32 temperature
int32 voltage
int32 error
```

## StepperMotorCommand (Message)¶

```
uint8 cmd_type
uint8 CMD_TYPE_POSITION=1
uint8 CMD_TYPE_VELOCITY=2
uint8 CMD_TYPE_EFFORT=3
uint8 CMD_TYPE_TORQUE=4

uint8[] motors_id
int32[] params
```

## StepperArrayMotorHardwareStatus (Message)¶

```
std_msgs/Header header
stepper_driver/StepperMotorHardwareStatus[] motors_hw_status
```

**Dxl\_debug\_tools**

This package offers scripts to change ping/scan DXL motors and changes register values of these motors.

**Niryo robot - Send DXL custom value**

This script can be launched via:

```
rosrun niryo_robot_debug send_custom_dx1_value.py
```

## Parameters - Send DXL custom value¶

- **-type [Number]:** Motor type (3 for XL-320, 2 for XL-430)
- **-id [Number]:** Motor ID
- **-address [Number]:** Register address to modify
- **-value [Number]:** Value to store at the register address given
- **-size [Number]:** Size in bytes of the value given

**Niryo robot - Dxl debug tools**

This script can be launched via:

```
rosrun niryo_robot_debug dxl_debug_tools
```

## Parameters - Dxl debug tools¶

- **-help / -h:** Print help message
- **-baudrate / -b [Number]:** Baudrate (1000000 by default)
- **-port / -p [Number]:** Set port
- **-id / -i [Number]:** Dxl motor ID (0 by default)
- **-scan:** Scan all Dxl motors on the bus
- **-ping:** Ping specific ID
- **-set-register [Number] [Number] [Number]:** Set a value to a register, parameters are in the order: register address / value / size (in bytes) of the value

**Control with Python ROS Wrapper**

Python Logo

In order to control Ned more easily than calling each topics & services one by one, a Python ROS Wrapper has been built on top of ROS.

For instance, a script realizing a movej via Python ROS Wrapper will look like:

```
niryo_robot = NiryoRosWrapper()
niryo_robot.move_joints(0.1, -0.2, 0.0, 1.1, -0.5, 0.2)
```

What this code is doing in a hidden way:

- It generates a [RobotMove Action Goal](#) and set it as a joint command with the corresponding joints value.
- Sends goal to the Commander Action Server.
- Waits for the Commander Action Server to set Action as finished.
- Checks if action finished with a success.

In this section, we will give some examples on how to use the Python ROS Wrapper to control Ned, as well as a complete documentation of the functions available in the Ned Python ROS Wrapper.

**💡 Hint**

The Python ROS Wrapper forces the user to write his code directly in the robot, or, at least, copy the code on the robot via a terminal command. If you do not want that, and run code directly from your computer you can use the python Package [PyNiryo](#) ([index.html#pyniryo](#)).

**Before running your programs****The variable PYTHONPATH**

The Python interpreter needs to have all used packages in the environment variable **PYTHONPATH**, to do that, you need to have sourced your ROS environment:

- If you are coding directly on your robot, it is made directly in every terminal.
- If your are using simulation, be sure to have followed the setup from [Setup Ned ROS environment](#).

## Required piece of code

To run, your program will need some imports & initialization. We give you below the piece of code you must use to make Python ROS Wrapper work:

```
#!/usr/bin/env python

# Imports
from niroo_robot_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

niryo_robot = NiryoRosWrapper()

# -- YOUR CODE HERE -- #
```

You have now everything you need to control the robot through its Python ROS Wrapper. To run a script, simply use the command `python my_script.py`.

## Examples: Basics

In this file, two short programs are implemented & commented in order to help you understand the philosophy behind the Python ROS Wrapper.

### Danger

If you are using the real robot, make sure the environment around is clear.

## Your first move joint

The following example shows a first use case. It's a simple MoveJ.

```
#!/usr/bin/env python

# Imports
from niroo_robot_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

# Connecting to the ROS Wrapper & calibrating if needed
niryo_robot = NiryoRosWrapper()
niryo_robot.calibrate_auto()

# Moving joint
niryo_robot.move_joints(0.1, -0.2, 0.0, 1.1, -0.5, 0.2)
```

### Code details - First MoveJ

First of all, we indicate to the shell that we are running a Python Script:

```
#!/usr/bin/env python
```

Then, we import the API package to be able to access functions:

```
from niroo_robot_python_ros_wrapper import *
```

Then, we install a ROS Node in order to communicate with ROS master:

```
import rospy
# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')
```

We start a **NiryoRosWrapper** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper](#)) instance:

```
niryo_robot = NiryoRosWrapper()
```

Once the connection is done, we calibrate the robot using its **calibrate\_auto()** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.calibrate\\_auto](#)) function:

```
niryo_robot.calibrate_auto()
```

As the robot is now calibrated, we can do a Move Joints by giving the 6 axis positions in radians! To do so, we use **move\_joints()** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.move\\_joints](#)):

```
niryo_robot.move_joints(0.1, -0.2, 0.0, 1.1, -0.5, 0.2)
```

## Your first pick and place

For our second example, we are going to develop an algorithm of pick and place:

```

#!/usr/bin/env python

# Imports
from niryo_robot_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_robot_example_python_ros_wrapper')

# Connecting to the ROS Wrapper & calibrating if needed
niryo_robot = NiryoRosWrapper()
niryo_robot.calibrate_auto()

# Updating tool
niryo_robot.update_tool()

# Opening Gripper/Pushing Air
niryo_robot.release_with_tool()
# Going to pick pose
niryo_robot.move_pose(0.2, 0.1, 0.14, 0.0, 1.57, 0)
# Picking
niryo_robot.grasp_with_tool()
# Moving to place pose
niryo_robot.move_pose(0.2, -0.1, 0.14, 0.0, 1.57, 0)
# Placing !
niryo_robot.release_with_tool()

```

Code details - first pick and place¶

First of all, we do the imports and start a ROS Node:

```

#!/usr/bin/env python

from niryo_robot_python_ros_wrapper import *
import rospy

rospy.init_node('niryo_robot_example_python_ros_wrapper')

```

Then, create a **NiryoRosWrapper** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper](#)) instance & calibrate the robot:

```

niryo_robot = NiryoRosWrapper()
niryo_robot.calibrate_auto()

```

Then, we equip the tool with **update\_tool()** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.update\\_tool](#))

```
niryo_robot.update_tool()
```

Now that our initialization is done, we can open the Gripper (or push air from the Vacuum pump) with **release\_with\_tool()** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.release\\_with\\_tool](#)), go to the picking pose via **move\_pose()** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.move\\_pose](#)) & then catch an object with **grasp\_with\_tool()** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.grasp\\_with\\_tool](#))!

```

# Opening Gripper/Pushing Air
niryo_robot.release_with_tool()
# Going to pick pose
niryo_robot.move_pose(0.2, 0.1, 0.14, 0.0, 1.57, 0)
# Picking
niryo_robot.grasp_with_tool()

```

We now get to the place pose, and place the object.

```

# Moving to place pose
niryo_robot.move_pose(0.2, -0.1, 0.14, 0.0, 1.57, 0)
# Placing !
niryo_robot.release_with_tool()

```

## Notes - Basics examples

You may not have fully understood how to move the robot and use tools of Ned and that is totally fine because you will find more details on another examples page! The important thing to remember from this page is how to import the library & connect to the robot.

## Examples: Movement

This document shows how to control Ned in order to make Move Joints & Move Pose.

If you want see more, you can look at [API - Joints & Pose](#) ([index.html#joints-pose](#)).

### Danger

If you are using the real robot, make sure the environment around is clear.

## Joints

To do a movej, you should pass 6 floats: (j1, j2, j3, j4, j5, j6) to the method **move\_joints()** ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.move\\_joints](#)):

```
#!/usr/bin/env python

# Imports
from niryo_robot_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

niryo_robot = NiryoRosWrapper()
niryo_robot.move_joints(0.0, 0.0, 0.0, 0.0, 0.0, 0.0)
```

To get joints, we use `get_joints()` ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.get\\_joints](#)):

```
joints = niryo_robot.get_joints()
j1, j2, j3, j4, j5, j6 = joints
```

## Pose

To do a moveP, you should pass 6 floats: (x, y, z, roll, pitch, yaw) to the method `move_pose()` ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.move\\_pose](#)).

See on this example:

```
#!/usr/bin/env python

# Imports
from niryo_robot_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

niryo_robot = NiryoRosWrapper()
niryo_robot.move_pose(0.25, 0.0, 0.25, 0.0, 0.0, 0.0)
```

To get pose, we use `get_pose()` ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.get\\_pose](#)):

```
x, y, z, roll, pitch, yaw = niryo_robot.get_pose()
```

## Examples: tool action

This page shows how to control Ned's tools via the Python ROS Wrapper.

If you want see more, you can look at [API - Tools](#) ([index.html#tools](#)).

### Danger

If you are using the real robot, make sure the environment around it is clear.

## Tool control

### Equip tool

In order to use a tool, it should be mechanically plugged to the robot but also connected software wise.

To do that, we should use the function `update_tool()` ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.update\\_tool](#)) which takes no argument. It will scan motor connections and set the new tool!

The line to equip a new tool is:

```
niryo_robot.update_tool()
```

### Grasping

To grasp with any tool, you can use the function: `grasp_with_tool()` ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.grasp\\_with\\_tool](#)). This action corresponds to:

- Close gripper for Grippers.
- Pull Air for Vacuum pump.
- Activate for Electromagnet.

The code to grasp is:

```
#!/usr/bin/env python

# Imports
from niryo_robot_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

# Connecting to the ROS Wrapper & calibrating if needed
niryo_robot = NiryoRosWrapper()
niryo_robot.calibrate_auto()

# Updating tool
niryo_robot.update_tool()

# Grasping
niryo_robot.grasp_with_tool()
```

To grasp by specifying parameters:

```
#!/usr/bin/env python

# Imports
from niro_ned_example_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

# Connecting to the ROS Wrapper & calibrating if needed
niryo_robot = NiryoRosWrapper()
niryo_robot.calibrate_auto()

# Updating tool
tool_used = ToolID.XXX
niryo_robot.update_tool()

if tool_used in [ToolID.GRIPPER_1, ToolID.GRIPPER_2, ToolID.GRIPPER_3, ToolID.GRIPPER_4]:
    niryo_robot.close_gripper(speed=500)
elif tool_used == ToolID.ELECTROMAGNET_1:
    pin_electromagnet = PinID.XXX
    niryo_robot.setup_electromagnet(pin_electromagnet)
    niryo_robot.activate_electromagnet(pin_electromagnet)
elif tool_used == ToolID.VACUUM_PUMP_1:
    niryo_robot.pull_air_vacuum_pump()
```

#### Releasing

To release with any tool, you can use the function: `release_with_tool()` ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.NiryoRosWrapper.release\\_with\\_tool](#)). This action correspond to:

- Open gripper for Grippers.
- Push Air for Vacuum pump.
- Deactivate for Electromagnet.

The line to release is:

```
niryo_robot.release_with_tool()
```

To release by specifying parameters:

```
#!/usr/bin/env python

# Imports
from niro_ned_example_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

# Connecting to the ROS Wrapper & calibrating if needed
niryo_robot = NiryoRosWrapper()
niryo_robot.calibrate_auto()

# Updating tool
tool_used = ToolID.XXX
niryo_robot.update_tool()

if tool_used in [ToolID.GRIPPER_1, ToolID.GRIPPER_2, ToolID.GRIPPER_3, ToolID.GRIPPER_4]:
    niryo_robot.open_gripper(speed=500)
elif tool_used == ToolID.ELECTROMAGNET_1:
    pin_electromagnet = PinID.XXX
    niryo_robot.setup_electromagnet(pin_electromagnet)
    niryo_robot.deactivate_electromagnet(pin_electromagnet)
elif tool_used == ToolID.VACUUM_PUMP_1:
    niryo_robot.push_air_vacuum_pump(tool_used)
```

#### Pick & place with tools

There are a plenty of ways to realize a pick and place with the ROS Wrapper. Methods will be presented from the lowest to highest level.

Code used will be:

```
# Imports
from niro_ned_example_python_ros_wrapper import *

gripper_used = ToolID.XXX # Tool used for picking

# The pick pose
pick_pose = (0.25, 0., 0.15, 0., 1.57, 0.0)
# The Place pose
place_pose = (0., -0.25, 0.1, 0., 1.57, -1.57)

def pick_n_place_version_x(niryo_ned):
    # -- SOME CODE -- #

if __name__ == '__main__':
    niryo_robot = NiryoRosWrapper()
    niryo_robot.calibrate_auto()
    pick_n_place_version_x(niryo_robot)
```

First solution: the heaviest

Everything is done by hand:

```

def pick_n_place_version_1(niryo_ned):
    height_offset = 0.05 # Offset according to Z-Axis to go over pick & place poses
    gripper_speed = 400

    # Going Over Object
    niryo_ned.move_pose(pick_pose[0], pick_pose[1], pick_pose[2] + height_offset,
                         pick_pose[3], pick_pose[4], pick_pose[5])

    # Opening Gripper
    niryo_ned.open_gripper(gripper_speed)
    # Going to picking place and closing gripper
    niryo_ned.move_pose(pick_pose[0], pick_pose[1], pick_pose[2],
                         pick_pose[3], pick_pose[4], pick_pose[5])
    niryo_ned.close_gripper(gripper_speed)

    # Raising
    niryo_ned.move_pose(pick_pose[0], pick_pose[1], pick_pose[2] + height_offset,
                         pick_pose[3], pick_pose[4], pick_pose[5])

    # Going Over Place pose
    niryo_ned.move_pose(place_pose[0], place_pose[1], place_pose[2] + height_offset,
                         place_pose[3], place_pose[4], place_pose[5])
    # Going to Place pose
    niryo_ned.move_pose(place_pose[0], place_pose[1], place_pose[2],
                         place_pose[3], place_pose[4], place_pose[5])
    # Opening Gripper
    niryo_ned.open_gripper(gripper_speed)
    # Raising
    niryo_ned.move_pose(place_pose[0], place_pose[1], place_pose[2] + height_offset,
                         place_pose[3], place_pose[4], place_pose[5])

```

Second solution: pick from pose & place from pose functions¶

We use predefined functions:

```

def pick_n_place_version_3(niryo_ned):
    # Pick
    niryo_ned.pick_from_pose(*pick_pose)
    # Place
    niryo_ned.place_from_pose(*place_pose)

```

Third solution: all in one¶

We use THE predefined function:

```

def pick_n_place_version_4(niryo_ned):
    # Pick & Place
    niryo_ned.pick_and_place(pick_pose, place_pose)

```

## Examples: Conveyor Belt

This document shows how to use Ned's Conveyor Belt.

If you want see more about Ned's Conveyor Belt functions, you can look at [API - Conveyor](#).

### Note

Imports & initialisation are not mentioned, but you should not forget it!

## Simple Conveyor Belt control

This short example shows how to connect a Conveyor Belt, activate the connection and launch its motor:

```

niryo_robot = NiryoRosWrapper()

# Activating connexion with conveyor and storing ID
conveyor_id = niryo_robot.set_conveyor()

# Running conveyor at 50% of its maximum speed, in Forward direction
niryo_robot.control_conveyor(conveyor_id, True, 100, ConveyorDirection.FORWARD)

# Stopping robot motor
niryo_robot.control_conveyor(conveyor_id, True, 0, ConveyorDirection.FORWARD)

# Deactivating connexion with conveyor
niryo_robot.unset_conveyor(conveyor_id)

```

## Advanced Conveyor Belt control

This example shows how to do a certain amount of pick & place by using the Conveyor Belt with the infrared sensor:

```

def run_conveyor(robot, conveyor):
    robot.control_conveyor(conveyor, bool_control_on=True,
                           speed=50, direction=ConveyorDirection.FORWARD)

# -- Setting variables
sensor_pin_id = PinID.GPIO_1A

catch_nb = 5

# The pick pose
pick_pose = [0.25, 0., 0.15, 0., 1.57, 0.0]
# The Place pose
place_pose = [0.0, -0.25, 0.1, 0., 1.57, -1.57]

# -- MAIN PROGRAM

niryo_robot = NiryoRosWrapper()

# Activating connexion with conveyor
conveyor_id = niryo_robot.set_conveyor()

for i in range(catch_nb):
    run_conveyor(niryo_robot, conveyor_id)
    while niryo_robot.digital_read(sensor_pin_id) == PinState.LOW:
        niryo_robot.wait(0.1)

    # Stopping robot motor
    niryo_robot.control_conveyor(conveyor_id, True, 0, ConveyorDirection.FORWARD)
    # Making a pick & place
    niryo_robot.pick_and_place(pick_pose, place_pose)

# Deactivating connexion with conveyor
niryo_robot.unset_conveyor(conveyor_id)

```

## Examples: Vision

This document shows how to use Ned's Vision Set.

If you want see more about Ned's Vision functions, you can look at [API - Vision](#) ([index.html#vision](#)).

### Beforehand

To realize the following examples, you need to have create a workspace.

As the examples start always the same, there is the code you need to add at the beginning of all of them:

```

#!/usr/bin/env python

# Imports
from niryo_robot_python_ros_wrapper import *
import rospy

# Initializing ROS node
rospy.init_node('niryo_ned_example_python_ros_wrapper')

niryo_robot = NiryoRosWrapper()

# - Constants
workspace_name = "workspace_1" # Robot's Workspace Name

# The observation pose
observation_pose = (0.18, 0., 0.35, 0., 1.57, -0.2)
# The Place pose
place_pose = (0., -0.25, 0.1, 0., 1.57, -1.57)

# - Main Program
# Calibrate robot if robot needs calibration
niryo_robot.calibrate_auto()
# Changing tool
niryo_robot.update_tool()

```



**Simple Vision pick**

This short example shows how to do your first vision pick:

```

niryo_robot.move_pose(*observation_pose)
# Trying to pick target using camera
ret = niryo_robot.vision_pick(workspace_name,
                               height_offset=0.0,
                               shape=ObjectShape.ANY,
                               color=ObjectColor.ANY)
obj_found, shape_ret, color_ret = ret
if obj_found:
    niryo_robot.place_from_pose(*place_pose)

niryo_robot.set_learning_mode(True)

```

## Python ROS Wrapper documentation

This file presents the different Functions, Classes & Enums available with the API.

- [API functions](#)
- [Enums](#)

### API functions

This class allows you to control the robot via internal API.

By controlling, we mean:

- Moving the robot.
- Using Vision.
- Controlling Conveyors Belt.
- Playing with hardware.

List of functions subsections:

- [Main purpose functions](#)
- [Joints & Pose](#)
- [Saved poses](#)
- [Pick & place](#)
- [Trajectories](#)
- [Tools](#)
- [Hardware](#)
- [Conveyor Belt](#)
- [Vision](#)

#### Main purpose functions

##### `class NiryoRosWrapper`

###### `calibrate_auto()`

Call service to calibrate motors then wait for its end. If failed, raise NiryoRosWrapperException

**Returns:** status, message

**Return type:** (int, str)

###### `calibrate_manual()`

Call service to calibrate motors then wait for its end. If failed, raise NiryoRosWrapperException

**Returns:** status, message

**Return type:** (int, str)

###### `get_learning_mode()`

Use /niryo\_robot/learning\_mode/state topic subscriber to get learning mode status

**Returns:** `True` if activate else `False`

**Return type:** bool

###### `set_learning_mode(set_bool)`

Call service to set\_learning\_mode according to set\_bool. If failed, raise NiryoRosWrapperException

**Parameters:** `set_bool (bool)` – `True` to activate, `False` to deactivate

**Returns:** status, message

**Return type:** (int, str)

###### `set_arm_max_velocity(percentage)`

Set relative max velocity (in %)

**Parameters:** `percentage (int)` – Percentage of max velocity

**Returns:** status, message

**Return type:** (int, str)

#### Joints & Pose

##### `class NiryoRosWrapper`

**get\_joints()**

Use /joint\_states topic to get joints status

**Returns:** list of joints value

**Return type:** list[float]

**get\_pose()**

Use /niryo\_robot/robot\_state topic to get pose status

**Returns:** RobotState object (position.x/y/z && rpy.roll/pitch/yaw && orientation.x/y/z/w)

**Return type:** RobotState

**get\_pose\_as\_list()**

Use /niryo\_robot/robot\_state topic to get pose status

**Returns:** list corresponding to [x, y, z, roll, pitch, yaw]

**Return type:** list[float]

**move\_joints(j1, j2, j3, j4, j5, j6)**

Execute Move joints action

**Parameters:**

- j1 (float) –
- j2 (float) –
- j3 (float) –
- j4 (float) –
- j5 (float) –
- j6 (float) –

**Returns:** status, message

**Return type:** (int (<https://docs.python.org/3/library/functions.html#int>), str (<https://docs.python.org/3/library/stdtypes.html#str>))

**move\_to\_sleep\_pose()**

Move to Sleep pose which allows the user to activate the learning mode without the risk of the robot hitting something because of gravity

**Returns:** status, message

**Return type:** (int, str)

**move\_pose(x, y, z, roll, pitch, yaw)**

Move robot end effector pose to a (x, y, z, roll, pitch, yaw) pose.

**Parameters:**

- x (float) –
- y (float) –
- z (float) –
- roll (float) –
- pitch (float) –
- yaw (float) –

**Returns:** status, message

**Return type:** (int (<https://docs.python.org/3/library/functions.html#int>), str (<https://docs.python.org/3/library/stdtypes.html#str>))

**shift\_pose(axis, value)**

Execute Shift pose action

**Parameters:**

- axis (ShiftPose) – Value of RobotAxis enum corresponding to where the shift happens
- value (float) – shift value

**Returns:** status, message

**Return type:** (int (<https://docs.python.org/3/library/functions.html#int>), str (<https://docs.python.org/3/library/stdtypes.html#str>))

**shift\_linear\_pose(axis, value)**

Execute Shift pose action with a linear trajectory

**Parameters:**

- axis (ShiftPose) – Value of RobotAxis enum corresponding to where the shift happens
- value (float) – shift value

**Returns:** status, message

**Return type:** (int (<https://docs.python.org/3/library/functions.html#int>), str (<https://docs.python.org/3/library/stdtypes.html#str>))

**move\_linear\_pose(x, y, z, roll, pitch, yaw)**

Move robot end effector pose to a (x, y, z, roll, pitch, yaw) pose, with a linear trajectory

**Parameters:**

- `x` (`float`) –
- `y` (`float`) –
- `z` (`float`) –
- `roll` (`float`) –
- `pitch` (`float`) –
- `yaw` (`float`) –

**Returns:** status, message

**Return type:** (`int` (<https://docs.python.org/3/library/functions.html#int>), `str` (<https://docs.python.org/3/library/stdtypes.html#str>))

#### `set_jog_use_state(state)`

Turn jog controller On or Off

**Parameters:** `state` (`bool`) – `True` to turn on, else `False`

**Returns:** status, message

**Return type:** (`int, str`)

#### `jog_joints_shift(shift_values)`

Make a Jog on joints position

**Parameters:** `shift_values` (`list[float]`) – list corresponding to the shift to be applied to each joint

**Returns:** status, message

**Return type:** (`int, str`)

#### `jog_pose_shift(shift_values)`

Make a Jog on end-effector position

**Parameters:** `shift_values` (`list[float]`) – list corresponding to the shift to be applied to the position

**Returns:** status, message

**Return type:** (`int, str`)

#### `forward_kinematics(j1, j2, j3, j4, j5, j6)`

Compute forward kinematics

**Parameters:**

- `j1` (`float`) –
- `j2` (`float`) –
- `j3` (`float`) –
- `j4` (`float`) –
- `j5` (`float`) –
- `j6` (`float`) –

**Returns:** list corresponding to [x, y, z, roll, pitch, yaw]

**Return type:** `list` (<https://docs.python.org/3/library/stdtypes.html#list>)[`float` (<https://docs.python.org/3/library/functions.html#float>)]

#### `inverse_kinematics(x, y, z, roll, pitch, yaw)`

Compute inverse kinematics

**Parameters:**

- `x` (`float`) –
- `y` (`float`) –
- `z` (`float`) –
- `roll` (`float`) –
- `pitch` (`float`) –
- `yaw` (`float`) –

**Returns:** list of joints value

**Return type:** `list` (<https://docs.python.org/3/library/stdtypes.html#list>)[`float` (<https://docs.python.org/3/library/functions.html#float>)]

### Saved poses

#### `class NiryoRosWrapper`

##### `move_pose_saved(pose_name)`

Move robot end effector pose to a pose saved

**Parameters:** `pose_name` (`str`) –

**Returns:** status, message

**Return type:** (`int, str`)

##### `get_pose_saved(pose_name)`

Get saved pose from robot intern storage Will raise error if position does not exist

**Parameters:** `pose_name (str)` - Pose Name

**Returns:** x, y, z, roll, pitch, yaw

**Return type:** tuple[float]

#### save\_pose(name, x, y, z, roll, pitch, yaw)

Save pose in robot's memory

**Parameters:**

- `name (str)` -
- `x (float)` -
- `y (float)` -
- `z (float)` -
- `roll (float)` -
- `pitch (float)` -
- `yaw (float)` -

**Returns:** status, message

**Return type:** (int (<https://docs.python.org/3/library/functions.html#int>), str (<https://docs.python.org/3/library/stdtypes.html#str>))

#### delete\_pose(name)

Send delete command to the pose manager service

**Parameters:** `name (str)` -

**Returns:** status, message

**Return type:** (int, str)

#### get\_saved\_pose\_list()

Ask the pose manager service which positions are available

**Returns:** list of positions name

**Return type:** list[str]

### Pick & place

#### class NiryoRosWrapper

##### pick\_from\_pose(x, y, z, roll, pitch, yaw)

Execute a picking from a position. If an error happens during the movement, error will be raised. A picking is described as : - going over the object - going down until height = z - grasping with tool - going back over the object

**Parameters:**

- `x (float)` -
- `y (float)` -
- `z (float)` -
- `roll (float)` -
- `pitch (float)` -
- `yaw (float)` -

**Returns:** status, message

**Return type:** (int (<https://docs.python.org/3/library/functions.html#int>), str (<https://docs.python.org/3/library/stdtypes.html#str>))

##### place\_from\_pose(x, y, z, roll, pitch, yaw)

Execute a placing from a position. If an error happens during the movement, error will be raised. A placing is described as : - going over the place - going down until height = z - releasing the object with tool - going back over the place

**Parameters:**

- `x (float)` -
- `y (float)` -
- `z (float)` -
- `roll (float)` -
- `pitch (float)` -
- `yaw (float)` -

**Returns:** status, message

**Return type:** (int (<https://docs.python.org/3/library/functions.html#int>), str (<https://docs.python.org/3/library/stdtypes.html#str>))

##### pick\_and\_place(pick\_pose, place\_pose, dist\_smoothing=0.0)

Execute a pick and place. If an error happens during the movement, error will be raised. -> Args param is for development purposes

**Parameters:**

- `pick_pose (list[float])` -
- `place_pose (list[float])` -
- `dist_smoothing (float)` - Distance from waypoints before smoothing trajectory

**Returns:** status, message

**Return type:** ([int](https://docs.python.org/3/library/functions.html#int) (<https://docs.python.org/3/library/functions.html#int>), [str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>))

## Trajectories¶

**class** `NiryoRosWrapper`

**get\_trajectory\_saved(trajectory\_name)**

Get saved trajectory from robot intern storage Will raise error if position does not exist

**Parameters:** `trajectory_name (str)` –

**Raises:** `NiryoRosWrapperException` – If trajectory file doesn't exist

**Returns:** list of [x, y, z, qx, qy, qz, qw]

**Return type:** ([list](https://docs.python.org/3/library/functions.html#list)[[float](https://docs.python.org/3/library/functions.html#float)])

**execute\_trajectory\_saved(trajectory\_name)**

Execute trajectory saved in Robot internal storage

**Parameters:** `trajectory_name (str)` –

**Returns:** status, message

**Return type:** ([int](https://docs.python.org/3/library/functions.html#int), [str](https://docs.python.org/3/library/functions.html#str))

**execute\_trajectory\_from\_poses(list\_poses\_raw, dist\_smoothing=0.0)**

Execute trajectory from a list of pose

**Parameters:**

- `list_poses_raw (list[list[float]])` – list of [x, y, z, qx, qy, qz, qw] or list of [x, y, z, roll, pitch, yaw]
- `dist_smoothing (float)` – Distance from waypoints before smoothing trajectory

**Returns:** status, message

**Return type:** ([int](https://docs.python.org/3/library/functions.html#int) (<https://docs.python.org/3/library/functions.html#int>), [str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>))

**execute\_trajectory\_from\_poses\_and\_joints(list\_pose\_joints, list\_type=None, dist\_smoothing=0.0)**

Execute trajectory from list of poses and joints

**Parameters:**

- `list_pose_joints (list[list[float]])` – List of [x,y,z,qx,qy,qz,qw] or list of [x,y,z,roll,pitch,yaw] or a list of [j1,j2,j3,j4,j5,j6]
- `list_type (list[string])` – List of string 'pose' or 'joint', or ['pose'] (if poses only) or ['joint'] (if joints only). If None, it is assumed there are only poses in the list.
- `dist_smoothing (float)` – Distance from waypoints before smoothing trajectory

**Returns:** status, message

**Return type:** ([int](https://docs.python.org/3/library/functions.html#int) (<https://docs.python.org/3/library/functions.html#int>), [str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>))

**save\_trajectory(trajectory\_name, list\_poses\_raw)**

Save trajectory object and send it to the trajectory manager service

**Parameters:**

- `trajectory_name (str)` – name which will have the trajectory
- `list_poses_raw (list[list[float]])` – list of [x, y, z, qx, qy, qz, qw] or list of [x, y, z, roll, pitch, yaw]

**Returns:** status, message

**Return type:** ([int](https://docs.python.org/3/library/functions.html#int) (<https://docs.python.org/3/library/functions.html#int>), [str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>))

**delete\_trajectory(trajectory\_name)**

Send delete command to the trajectory manager service

**Parameters:** `trajectory_name (str)` – name

**Returns:** status, message

**Return type:** ([int](https://docs.python.org/3/library/functions.html#int), [str](https://docs.python.org/3/library/functions.html#str))

**get\_saved\_trajectory\_list()**

Ask the pose trajectory service which trajectories are available

**Returns:** list of trajectory name

**Return type:** ([list](https://docs.python.org/3/library/functions.html#list)[[str](https://docs.python.org/3/library/functions.html#str)])

## Tools¶

**class** `NiryoRosWrapper`

**get\_current\_tool\_id()**

Use /niryo\_robot\_hardware/tools/current\_id topic to get current tool id

**Returns:** Tool Id

**Return type:** ToolID

#### update\_tool()

Call service niryo\_robot\_tools\_commander/update\_tool to update tool

**Returns:** status, message

**Return type:** (int, str)

#### grasp\_with\_tool(pin\_id=-1)

Grasp with the tool linked to tool\_id. This action correspond to - Close gripper for Grippers - Pull Air for Vacuum pump - Activate for Electromagnet

**Parameters:** pin\_id (*PinID*) – [Only required for electromagnet] Pin ID of the electromagnet

**Returns:** status, message

**Return type:** (int, str)

#### release\_with\_tool(pin\_id=-1)

Release with the tool associated to tool\_id. This action correspond to - Open gripper for Grippers - Push Air for Vacuum pump - Deactivate for Electromagnet

**Parameters:** pin\_id (*PinID*) – [Only required for electromagnet] Pin ID of the electromagnet

**Returns:** status, message

**Return type:** (int, str)

#### open\_gripper(speed=500)

Open gripper with a speed 'speed'

**Parameters:** speed (*int*) – Default -> 500

**Returns:** status, message

**Return type:** (int, str)

#### close\_gripper(speed=500)

Close gripper with a speed 'speed'

**Parameters:** speed (*int*) – Default -> 500

**Returns:** status, message

**Return type:** (int, str)

#### pull\_air\_vacuum\_pump()

Pull air

**Returns:** status, message

**Return type:** (int, str)

#### push\_air\_vacuum\_pump()

Pull air

**Returns:** status, message

**Return type:** (int, str)

#### setup\_electromagnet(pin\_id)

Setup electromagnet on pin

**Parameters:** pin\_id (*PinID*) – Pin ID

**Returns:** status, message

**Return type:** (int, str)

#### activate\_electromagnet(pin\_id)

Activate electromagnet associated to electromagnet\_id on pin\_id

**Parameters:** pin\_id (*PinID*) – Pin ID

**Returns:** status, message

**Return type:** (int, str)

#### deactivate\_electromagnet(pin\_id)

Deactivate electromagnet associated to electromagnet\_id on pin\_id

**Parameters:** pin\_id (*PinID*) – Pin ID

**Returns:** status, message

**Return type:** (int, str)

#### enable\_tcp(enable=True)

Enables or disables the TCP function (Tool Center Point). If activation is requested, the last recorded TCP value will be applied. The default value depends on the gripper equipped. If deactivation is requested, the TCP will be coincident with the tool\_link.

**Parameters:** enable (Bool) – True to enable, False otherwise.

**Returns:** status, message

**Return type:** (int, str)

#### set\_tcp(x, y, z, roll, pitch, yaw)

Activates the TCP function (Tool Center Point) and defines the transformation between the tool\_link frame and the TCP frame.

**Parameters:**

- x ([float](#)) –
- y ([float](#)) –
- z ([float](#)) –
- roll ([float](#)) –
- pitch ([float](#)) –
- yaw ([float](#)) –

**Returns:** status, message

**Return type:** (int ([https://docs.python.org/3/library/functions.html#int](#)), str ([https://docs.python.org/3/library/stdtypes.html#str](#)))

#### reset\_tcp()

Reset the TCP (Tool Center Point) transformation. The TCP will be reset according to the tool equipped.

**Returns:** status, message

**Return type:** (int, str)

### Hardware

#### class NiryoRosWrapper

##### set\_pin\_mode(pin\_id, pin\_mode)

Set pin number pin\_id to mode pin\_mode

**Parameters:**

- pin\_id ([PinID](#)) –
- pin\_mode ([PinMode](#)) –

**Returns:** status, message

**Return type:** (int ([https://docs.python.org/3/library/functions.html#int](#)), str ([https://docs.python.org/3/library/stdtypes.html#str](#)))

##### digital\_write(pin\_id, digital\_state)

Set pin\_id state to pin\_state

**Parameters:**

- pin\_id ([PinID](#)) –
- digital\_state ([PinState](#)) –

**Returns:** status, message

**Return type:** (int ([https://docs.python.org/3/library/functions.html#int](#)), str ([https://docs.python.org/3/library/stdtypes.html#str](#)))

##### digital\_read(pin\_id)

Read pin number pin\_id and return its state

**Parameters:** pin\_id ([PinID](#)) –

**Returns:** state

**Return type:** PinState

##### get.hardware\_status()

Get hardware status : Temperature, Hardware version, motors names & types ...

**Returns:** Infos contains in a HardwareStatus object (see niryo\_robot\_msgs)

**Return type:** HardwareStatus

##### get.digital\_io\_state()

Get Digital IO state : Names, modes, states

**Returns:** Infos contains in a DigitalIOState object (see niryo\_robot\_msgs)

**Return type:** DigitalIOState

## Conveyor Belt

`class NiryoRosWrapper``set_conveyor()`

Scan for conveyor on can bus. If conveyor detected, return the conveyor ID

**Raises:** `NiryoRosWrapperException` –

**Returns:** ID

**Return type:** `ConveyorID`

`unset_conveyor(conveyor_id)`

Remove specific conveyor.

**Parameters:** `conveyor_id (ConveyorID)` – Basically, ConveyorID.ONE or ConveyorID.TWO

**Raises:** `NiryoRosWrapperException` –

**Returns:** status, message

**Return type:** `(int, str)`

`control_conveyor(conveyor_id, bool_control_on, speed, direction)`

Control conveyor associated to conveyor\_id. Then stops it if bool\_control\_on is False, else refreshes its speed and direction

**Parameters:**

- `conveyor_id (ConveyorID)` – ConveyorID.ID\_1 or ConveyorID.ID\_2
- `bool_control_on (bool)` – True for activate, False for deactivate
- `speed (int)` – target speed
- `direction (ConveyorDirection)` – Target direction

**Returns:** status, message

**Return type:** `(int (https://docs.python.org/3/library/functions.html#int), str (https://docs.python.org/3/library/stdtypes.html#str))`

## Vision

`class NiryoRosWrapper``get_compressed_image()`

Get last stream image in a compressed format

**Returns:** string containing a JPEG compressed image

**Return type:** `str`

`set_brightness(brightness_factor)`

Modify image brightness

**Parameters:** `brightness_factor (float)` – How much to adjust the brightness. 0.5 will give a darkened image, 1 will give the original image while 2 will enhance the brightness by a factor of 2.

**Returns:** status, message

**Return type:** `(int, str)`

`set_contrast(contrast_factor)`

Modify image contrast

**Parameters:** `contrast_factor (float)` – While a factor of 1 gives original image. Making the factor towards 0 makes the image greyer, while factor>1 increases the contrast of the image.

**Returns:** status, message

**Return type:** `(int, str)`

`set_saturation(saturation_factor)`

Modify image saturation

**Parameters:** `saturation_factor (float)` – How much to adjust the saturation. 0 will give a black and white image, 1 will give the original image while 2 will enhance the saturation by a factor of 2.

**Returns:** status, message

**Return type:** `(int, str)`

`get_target_pose_from_rel(workspace_name, height_offset, x_rel, y_rel, yaw_rel)`

Given a pose (x\_rel, y\_rel, yaw\_rel) relative to a workspace, this function returns the robot pose in which the current tool will be able to pick an object at this pose. The height\_offset argument (in m) defines how high the tool will hover over the workspace. If height\_offset = 0, the tool will nearly touch the workspace.

**Parameters:**

- **workspace\_name** ([str](#)) – name of the workspace
- **height\_offset** ([float](#)) – offset between the workspace and the target height
- **x\_rel** ([float](#)) –
- **y\_rel** ([float](#)) –
- **yaw\_rel** ([float](#)) –

**Returns:** target\_pose

**Return type:** RobotState

#### `get_target_pose_from_cam(workspace_name, height_offset, shape, color)`

First detects the specified object using the camera and then returns the robot pose in which the object can be picked with the current tool

**Parameters:**

- **workspace\_name** ([str](#)) – name of the workspace
- **height\_offset** ([float](#)) – offset between the workspace and the target height
- **shape** ([ObjectShape](#)) – shape of the target
- **color** ([ObjectColor](#)) – color of the target

**Returns:** object\_found, object\_pose, object\_shape, object\_color

**Return type:** ([bool](#) ([https://docs.python.org/3/library/functions.html#bool](#)), [RobotState](#), [str](#) ([https://docs.python.org/3/library/stdtypes.html#str](#)), [str](#) ([https://docs.python.org/3/library/stdtypes.html#str](#)))

#### `vision_pick_w_obs_joints(workspace_name, height_offset, shape, color, observation_joints)`

Move Joints to observation\_joints, then execute a vision pick

#### `vision_pick_w_obs_pose(workspace_name, height_offset, shape, color, observation_pose_list)`

Move Pose to observation\_pose, then execute a vision pick

#### `vision_pick(workspace_name, height_offset, shape, color)`

Picks the specified object from the workspace. This function has multiple phases: 1. detect object using the camera 2. prepare the current tool for picking 3. approach the object 4. move down to the correct picking pose 5. actuate the current tool 6. lift the object

**Parameters:**

- **workspace\_name** ([str](#)) – name of the workspace
- **height\_offset** ([float](#)) – offset between the workspace and the target height
- **shape** ([ObjectShape](#)) – shape of the target
- **color** ([ObjectColor](#)) – color of the target

**Returns:** object\_found, object\_shape, object\_color

**Return type:** ([bool](#) ([https://docs.python.org/3/library/functions.html#bool](#)), [ObjectShape](#) ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.ros\\_wrapper\\_enums.ObjectShape](#)), [ObjectColor](#) ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.ros\\_wrapper\\_enums.ObjectColor](#)))

#### `move_to_object(workspace, height_offset, shape, color)`

Same as `get_target_pose_from_cam` but directly moves to this position

**Parameters:**

- **workspace** ([str](#)) – name of the workspace
- **height\_offset** ([float](#)) – offset between the workspace and the target height
- **shape** ([ObjectShape](#)) – shape of the target
- **color** ([ObjectColor](#)) – color of the target

**Returns:** object\_found, object\_shape, object\_color

**Return type:** ([bool](#) ([https://docs.python.org/3/library/functions.html#bool](#)), [ObjectShape](#) ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.ros\\_wrapper\\_enums.ObjectShape](#)), [ObjectColor](#) ([index.html#niryo\\_robot\\_python\\_ros\\_wrapper.ros\\_wrapper\\_enums.ObjectColor](#)))

#### `detect_object(workspace_name, shape, color)`

**Parameters:**

- **workspace\_name** ([str](#)) – name of the workspace
- **shape** ([ObjectShape](#)) – shape of the target
- **color** ([ObjectColor](#)) – color of the target

**Returns:** object\_found, object\_pose, object\_shape, object\_color

**Return type:** ([bool](#) ([https://docs.python.org/3/library/functions.html#bool](#)), [RobotState](#), [str](#) ([https://docs.python.org/3/library/stdtypes.html#str](#)), [str](#) ([https://docs.python.org/3/library/stdtypes.html#str](#)))

#### `get_camera_intrinsics()`

Get calibration object: camera intrinsics, distortions coefficients

**Returns:** raw camera intrinsics, distortions coefficients

**Return type:** ([list](#), [list](#))

#### `save_workspace_from_poses(name, list_poses_raw)`

Save workspace by giving the poses of the robot to point its 4 corners with the calibration Tip. Corners should be in the good order

**Parameters:**

- **name** ([str](#)) – workspace name, max 30 char.
- **list\_poses\_raw** ([list\[list\]](#)) – list of 4 corners pose

**Returns:** status, message

**Return type:** ([int](#) (<https://docs.python.org/3/library/functions.html#int>), [str](#) (<https://docs.python.org/3/library/stdtypes.html#str>))

#### `save_workspace_from_points(name, list_points_raw)`

Save workspace by giving the poses of its 4 corners in the good order

**Parameters:**

- **name** ([str](#)) – workspace name, max 30 char.
- **list\_points\_raw** ([list\[list\[float\]\]](#)) – list of 4 corners [x, y, z]

**Returns:** status, message

**Return type:** ([int](#) (<https://docs.python.org/3/library/functions.html#int>), [str](#) (<https://docs.python.org/3/library/stdtypes.html#str>))

#### `delete_workspace(name)`

Call workspace manager to delete a certain workspace

**Parameters:** **name** ([str](#)) – workspace name

**Returns:** status, message

**Return type:** ([int, str](#))

#### `get_workspace_poses(name)`

Get the 4 workspace poses of the workspace called 'name'

**Parameters:** **name** ([str](#)) – workspace name

**Returns:** List of the 4 workspace poses

**Return type:** [list\[list\]](#)

#### `get_workspace_ratio(name)`

Give the length over width ratio of a certain workspace

**Parameters:** **name** ([str](#)) – workspace name

**Returns:** ratio

**Return type:** [float](#)

#### `get_workspace_list()`

Ask the workspace manager service names of the available workspace

**Returns:** list of workspaces name

**Return type:** [list\[str\]](#)

## Enums

### `class ShiftPose`

`AXIS_X= 0`

`AXIS_Y= 1`

`AXIS_Z= 2`

`ROT_ROLL= 3`

`ROT_PITCH= 4`

`ROT_YAW= 5`

### `class ToolID`

Tools IDs (need to match tools ids in niro\_niryo\_robot\_tools\_commander package)

`NONE= 0`

`GRIPPER_1= 11`

`GRIPPER_2= 12`

`GRIPPER_3= 13`

**GRIPPER\_4=** 14**ELECTROMAGNET\_1=** 30**VACUUM\_PUMP\_1=** 31**class PinMode**

Pin Mode is either OUTPUT or INPUT

**OUTPUT=** 0**INPUT=** 1**class PinState**

Pin State is either LOW or HIGH

**LOW=** 0**HIGH=** 1**class PinID**

Pins ID

**GPIO\_1A=** 2**GPIO\_1B=** 3**GPIO\_1C=** 16**GPIO\_2A=** 26**GPIO\_2B=** 19**GPIO\_2C=** 6**SW\_1=** 12**SW\_2=** 13**class ConveyorID****NONE=** 0**ID\_1=** 12**ID\_2=** 13**class ConveyorDirection****FORWARD=** 1**BACKWARD=** -1**class ObjectColor****RED=** 'RED'**GREEN=** 'GREEN'**BLUE=** 'BLUE'**ANY=** 'ANY'**class ObjectShape****CIRCLE=** 'CIRCLE'**SQUARE=** 'SQUARE'**ANY=** 'ANY'**Use the Modbus TCP server**

**todo:** vu avec Etienne, supprimer la partie Modbus d'ici et mettre en place une redirection vers doc Modbus.

Ned is permanently running a Modbus TCP Server that enables Ned to communicate with a PLC, or another computer in the same network.

## Connection - Modbus TCP server

The Modbus TCP server is running on port 5020 by default.

## Description - Modbus TCP server

It has been built on top of the [pymodbus library](https://pymodbus.readthedocs.io/en/latest/index.html) (<https://pymodbus.readthedocs.io/en/latest/index.html>).

### All 4 Modbus datastores are implemented:

- Coil.
- Discrete Input.
- Holding Register.
- Input Register.

Each datastore has a different set of functionalities. Note that each datastore contains a completely different set of data.

Address tables start at 0.

### Coil

Each address contains a 1bit value.

READ/WRITE (the stored values correspond to the last given command, not the current robot state).

#### Accepted Modbus functions:

- 0x01: READ\_COILS
- 0x05: WRITE\_SINGLE\_COIL

This datastore can be used to set Digital I/O mode and state.

#### Digital I/O numbers used for Modbus:

- 0/100: 1A
- 1/101: 1B
- 2/102: 1C
- 3/103: 2A
- 4/104: 2B
- 5/105: 2C

Address	Description
0-5	Digital I/O mode (Input = 1, Output = 0)
100-105	Digital I/O state (High = 1, Low = 0)
200-299	Can be used to store your own variables

### Discrete input

Each address contains a 1bit value.

READ-ONLY

#### Accepted Modbus functions:

- 0x02: READ\_DISCRETE\_INPUTS

#### Digital I/O numbers used for Modbus:

- 0/100: 1A
- 1/101: 1B
- 2/102: 1C
- 3/103: 2A
- 4/104: 2B
- 5/105: 2C

Address	Description
0-5	Digital I/O mode (Input = 1, Output = 0)
100-105	Digital I/O state (High = 1, Low = 0)

This datastore can be used to read Digital I/O mode and state of the robot. See [Coil](#) above for digital I/O number mapping.

### Holding register

Each address contains a 16bits value.

READ/WRITE (the stored values correspond to the last given command, not the current robot state)

#### Accepted Modbus functions:

- 0x03: READ\_HOLDING\_REGISTERS
- 0x06: WRITE\_SINGLE\_REGISTER

Address	Description
0-5	Joints (mrad)
10-12	Position x,y,z (mm)
13-15	Orientation roll, pitch, yaw (mrad)
100	Sends Joint Move command with stored joints
101	Sends Pose Move command with stored position and orientation
102	Sends Pose Linear Move command with stored position and orientation
110	Stops current command execution
150	Is executing command flag
151	Last command result (status of the last command)
152	Contains data retrieved from last cmd (depends of the cmd)
153 - 158	Vision - Target pose result
159	Vision - Shape of the object found (-1: ANY, 1: CIRCLE, 2: SQUARE, 3: TRIANGLE, 0: NONE)
160	Vision - Color of the object found (-1: ANY, 1: BLUE, 2: RED, 3: GREEN, 0: NONE)
200-299	Can be used to store your own variables
300	Learning Mode (On = 1, Off = 0)
301	Joystick Enabled (On = 1, Off = 0)
310	Requests new calibration
311	Starts auto calibration
312	Starts manual calibration
401	Gripper open speed (100-1000)
402	Gripper close speed (100-1000)
500	Updates the tool id according to the gripper plugged (gripper 1: 11, gripper 2: 12, gripper 3: 13, vaccum pump: 31)
501	Stores the tool id
510	Opens gripper previously updated
511	Closes gripper with given id
512	Pulls air vacuum pump from given id
513	Pushes air vacuum pump from given id
520	Enables a Conveyor Belt newly connected [on success: store its ID at 152]
521	Detaches / disables Conveyor Belt with the Conveyor Belt ID given at 525
522	Control Conveyor Belt with the Conveyor Belt ID given at 525
523 [related to 522]	Conveyor Belt direction (backward = -1 , forward = 1)
524 [related to 522]	Conveyor Belt speed (0-100)(%)
525 [related to 520/521/522/526]	Stores the Conveyor Belt ID for all related command
526	Stops Conveyor Belt with the Conveyor Belt ID given at 525
600	TCP - Enables or disables the TCP function (Tool Center Point).
601	Activates the TCP function (Tool Center Point) and defines the transformation between the tool_link frame and the TCP frame.
610	Vision - Gets target pose from relative pose, with stored relative pose and height_offset
611	Vision - Gets target pose from camera, with stored workspace name, height offset, shape and color
612	Vision - Vision pick, with stored workspace name, height offset, shape and color
613	Vision - Moves to object, with stored workspace name, height offset, shape and color
614	Vision - Detects object, with stored workspace name, shape and color
620	Vision - Stores workspace's height offset
621	Vision - Stores relative pose x_rel
622	Vision - Stores relative pose y_rel
623	Vision - Stores relative pose yaw_rel
624	Vision - Stores requested shape (-1: ANY, 1: CIRCLE, 2: SQUARE, 3: TRIANGLE)
625	Vision - Stores requested color (-1: ANY, 1: BLUE, 2: RED, 3: GREEN)
626 - max 641	Vision - Stores workspace's name, as a string encoded in 16 bits hex (see examples on how to store a workspace name from a client)

**Input Register**

Each address contains a 16bits value.

READ-ONLY

**Accepted Modbus functions:**

- 0x04: READ\_INPUT\_REGISTERS

Address	Description
0-5	Joints (mrad)
10-12	Position x,y,z (mm)
13-15	Orientation roll, pitch, yaw (mrad)
200	Selected tool ID (0 for no tool)
300	Learning Mode activated
400	Motors connection up (Ok = 1, Not ok = 0)
401	Calibration needed flag
402	Calibration in progress flag
403	Raspberry Pi temperature
404	Raspberry Pi available disk size
405	Raspberry Pi ROS log size
406	RPI software version n.1
407	RPI software version n.2
408	RPI software version n.3
409	Hardware version (1 or 2)
530	Conveyor 1 connection state (Connected = 1 , Not connected = 0)
531	Conveyor 1 control status ( On = 0, Off = 1)
532	Conveyor 1 Speed (0-100 (%))
533	Conveyor 1 direction (Backward = -1, Forward = 1)
540	Conveyor 2 connection state (Connected = 1 , Not connected = 0)
541	Conveyor 2 control status ( On = 0, Off = 1)
542	Conveyor 2 Speed (0-100 (%))
543	Conveyor 2 direction (Backward = -1, Forward = 1)

#### Dependencies - Modbus TCP Server

- [pymodbus library](#)
- [Niryo\\_robot\\_msgs](#)
- [std\\_msgs](#)

#### More ways to control Ned

There is even more ways to control Ned.

If you are a beginner, look at [Blockly](#) section to understand the programming fundamentals.

If you want to go further, maybe experience your own image processing, multi-robot, AI... You can go to [PyNiryo](#) for more information.

#### Blockly

Check out Niryo Studio.

#### PyNiryo

As explained in the page [Use Ned's TCP server](#) (`index.html#use-ned-s-tcp-server`), a TCP Server is running on Ned, which allows it to receive commands from any external device.

PyNiryo is a Python package available on Pip which allows to command the Niryo Robots with easy Python Binding.

[Suggest a modification](#)

[Download as PDF](#)

# Niryo Modbus Documentation



In this document, we will focus on the Modbus/TCP server.

The Modbus/TCP server is running on **port 5020** by default. It has been built on top of the [pymodbus](https://pymodbus.readthedocs.io/en/latest/index.html) (<https://pymodbus.readthedocs.io/en/latest/index.html>) library. This enables you to make Ned communicates with a PLC, or another computer on the same network.

## Modbus Python library installation

To use the Modbus Python library, your workspace must have a Python interpreter with **Python 3** (3.6 or greater) or **Python 2** (2.7 or greater).

### Note

Download Python on the official [Python website](https://www.python.org/) (<https://www.python.org/>) and find more information about the installation on [this website](https://realpython.com/installing-python/) (<https://realpython.com/installing-python/>).

This installation requires the use of [pip](https://pypi.org/project/pip/) (<https://pypi.org/project/pip/>), the package manager included in Python.

Start with the installation of numpy:

```
pip install numpy
```

To use the **Modbus API**, you also need to install Modbus python library [pymodbus](https://pymodbus.readthedocs.io/en/latest/index.html) (<https://pymodbus.readthedocs.io/en/latest/index.html>):

```
pip install pymodbus
```

### Attention

- Pip can require administrator authorizations to install packages. In this case, add

```
sudo
```

before your command lines on Linux.

- If pip is not automatically installed with Python, please visit the following website: [pip installation](https://pypi.org/project/pip/) (<https://pypi.org/project/pip/>).

## Introduction

All 4 Modbus datastores are implemented:[Coils](#), [Discrete inputs](#), [Holding registers](#), [Input registers](#). Each datastore has a different set of functionalities. Note that each datastore contains a completely different set of data.

Discrete Input and Input register are **READ-ONLY** tables. Those have been used to keep the robot state.

Coil and Holding Register are **READ/WRITE** tables. Those have been used to give user commands to the robot. Hence, those 2 tables do not contain the robot state, but the last given command.

Address tables start at 0.

## Coils

---

Each address contains a 1bit value.

**READ/WRITE** (the stored values correspond to the last given command, not the current robot state)

Accepted Modbus functions:

- 0x01: READ\_COILS
- 0x05: WRITE\_SINGLE\_COIL

This datastore can be used to set Digital I/O mode and state. Digital I/O numbers used for Modbus:

- 0: 1A
- 1: 1B
- 2: 1C
- 3: 2A
- 4: 2B
- 5: 2C

Address	Description
0-5	Digital I/O mode (Input = 1, Output = 0)
100-105	Digital I/O state (High = 1, Low = 0)
200-299	Can be used to store your own variables

## Discrete inputs

---

Each address contains a 1bit value.

### READ-ONLY

Accepted Modbus functions:

- 0x02: READ\_DISCRETE\_INPUTS

This datastore can be used to read Digital I/O mode and state. See the [Coils](#) section above for digital I/O number mapping.

Address	Description
0-5	Digital I/O mode (Input = 1, Output = 0)
100-105	Digital I/O state (High = 1, Low = 0)

## Holding registers

Each address contains a 16bit value.

**READ/WRITE** (the stored values correspond to the last given command, not the current robot state)

Accepted Modbus functions:

- 0x03: READ\_HOLDING\_REGISTERS
- 0x06: WRITE\_SINGLE\_REGISTER

Address	Description
0-5	Joints (mrad)
10-12	Position x,y,z (mm)
13-15	Orientation roll, pitch, yaw (mrad)
100	Sends Joint Move command with stored joints
101	Sends Pose Move command with stored position and orientation
102	Sends Linear Pose Move command with stored position and orientation
110	Stops current command execution
150	Is executing command flag
151	Last command result*
152	Last command data result (if not vision related)
153 - 158	Vision - Target pose result
159	Vision - Shape of the object found (-1: ANY, 1: CIRCLE, 2: SQUARE, 3: TRIANGLE, 0: NONE)
160	Vision - Color of the object found (-1: ANY, 1: BLUE, 2: RED, 3: GREEN, 0: NONE)
200-299	Can be used to store your own variables
300	Learning Mode (On = 1, Off = 0)
301	Joystick Enabled (On = 1, Off = 0)
310	Requests new calibration
311	Starts auto calibration
312	Starts manual calibration
401	Gripper open speed (100-1000)
402	Gripper close speed (100-1000)
500	Updates the tool id according to the gripper plugged (gripper 1: 11, gripper 2: 12, gripper 3: 13, vacuum pump: 31)
501	Stores the tool id
510	Opens gripper previously updated

Address	Description
511	Closes gripper previously updated
512	Pulls air vacuum pump with id 31
513	Pushes air vacuum pump with id 31
520	Updates the conveyor id and enable it
521	Detaches or disables the conveyor previously enabled and updated
522	Starts the conveyor previously enabled and updated
523	Sets the conveyor direction (backward = number_to_raw_data(-1), forward = 1)
524	Sets the conveyor speed (0-100)(%)
525	Stores the conveyor id
526	Stops conveyor previously enabled and updated
600	TCP - Enables or disables the TCP function (Tool Center Point).
601	Activates the TCP function (Tool Center Point) and defines the transformation between the tool_link frame and the TCP frame.
610	Vision - Gets target pose from relative pose, with stored relative pose and height_offset
611	Vision - Gets target pose from camera, with stored workspace name, height offset, shape and color
612	Vision - Vision pick, with stored workspace name, height offset, shape and color
613	Vision - Moves to object, with stored workspace name, height offset, shape and color
614	Vision - Detects object, with stored workspace name, shape and color
620	Vision - Stores workspace's height offset
621	Vision - Stores relative pose x_rel
622	Vision - Stores relative pose y_rel
623	Vision - Stores relative pose yaw_rel
624	Vision - Stores requested shape (-1: ANY, 1: CIRCLE, 2: SQUARE, 3: TRIANGLE)
625	Vision - Stores requested color (-1: ANY, 1: BLUE, 2: RED, 3: GREEN)
626 - max 641	Vision - Stores workspace's name, as a string encoded in 16 bits hex (see examples on how to store a workspace name from a client)

\*' The “Last command result” gives you more information about the last executed command:

- 0: no result yet
- 1: success
- 2: command was rejected (invalid params, ...)
- 3: command was aborted
- 4: command was canceled
- 5: command had an unexpected error
- 6: command timeout
- 7: internal error

## Input registers

Each address contains a 16bit value.

**READ-ONLY.**

Accepted Modbus functions:

- 0x04: READ\_INPUT\_REGISTERS

Address	Description
0-5	Joints (mrad)
10-12	Position x,y,z (mm)
13-15	Orientation roll, pitch, yaw (mrad)
200	Selected tool ID (0 for no tool)
300	Learning Mode activated
400	Motors connection up (Ok = 1, Not ok = 0)
401	Calibration needed flag
402	Calibration in progress flag
403	Raspberry Pi temperature
404	Raspberry Pi available disk size
405	Raspberry Pi ROS log size
406	Ned RPI image version n.1
407	Ned RPI image version n.2
408	Ned RPI image version n.3
409	Hardware version (1 or 2)
530	Conveyor 1 connection state (Connected = 1 , Not connected = 0)
531	Conveyor 1 control status ( On = 0, Off = 1)
532	Conveyor 1 Speed (0-100 (%))
533	Conveyor 1 direction (Backward = -1, Forward = 1)
540	Conveyor 2 connection state (Connected = 1 , Not connected = 0)
541	Conveyor 2 control status ( On = 0, Off = 1)
542	Conveyor 2 Speed (0-100 (%))
543	Conveyor 2 direction (Backward = -1, Forward = 1)

## Repository

Examples of Modbus python lib can be found here [Python Modbus examples](https://github.com/NiryoRobotics/ned_ros/tree/master/niryo_robot_modbus/examples/) ([https://github.com/NiryoRobotics/ned\\_ros/tree/master/niryo\\_robot\\_modbus/examples/](https://github.com/NiryoRobotics/ned_ros/tree/master/niryo_robot_modbus/examples/)).

## Examples

In the examples folder, you can find several example scripts that control Ned. These scripts are commented to help you understand every step.

### Client Modbus Test

Calls several functions on the IO of Ned.

## Client Move Command

This script shows the calibration and Ned's moves.

## Client Modbus Conveyor Example

This script shows how to activate the Conveyor Belt through the Modbus Python API, set a direction, a speed, and start and stop the device.

## Client Modbus Vision Example

This script shows how to use the vision pick method from a Modbus Client, through the Modbus Python API. Ned picks a red object seen in its workspace and releases it on its left. Note that we use the **string\_to\_register** method to convert a string into an object storable in registers.

```
#!/usr/bin/env python

from pymodbus.client.sync import ModbusTcpClient
from pymodbus.payload import BinaryPayloadBuilder, BinaryPayloadDecoder
import time
from enum import Enum, unique

# Enums for shape and color. Those enums are the one used by the modbus server to receive requests
@unique
class ColorEnum(Enum):
    ANY = -1
    BLUE = 1
    RED = 2
    GREEN = 3
    NONE = 0

@unique
class ShapeEnum(Enum):
    ANY = -1
    CIRCLE = 1
    SQUARE = 2
    TRIANGLE = 3
    NONE = 0

# Functions to convert variables for/from registers

# Positive number : 0 - 32767
# Negative number : 32768 - 65535
def number_to_raw_data(val):
    if val < 0:
        val = (1 << 15) - val
    return val

def raw_data_to_number(val):
    if (val >> 15) == 1:
        val = - (val & 0x7FFF)
    return val

def string_to_register(string):
    # code a string to 16 bits hex value to store in register
    builder = BinaryPayloadBuilder()
    builder.add_string(string)
    payload = builder.to_registers()
    return payload

# ----- Modbus server related function

def back_to_observation():
    # To change
    # joint_real = [0.057, 0.604, -0.576, -0.078, -1.384, 0.253]
    joint_simu = [0, -0.092, 0, 0, -1.744, 0]

    joint_to_send = list(map(lambda j: int(number_to_raw_data(j * 1000)), joint_simu))
    client.write_registers(0, joint_to_send)
    client.write_register(100, 1)

    while client.read_holding_registers(150, count=1).registers[0] == 1:
```

```

time.sleep(0.01)

def register_workspace_name(ws_name):
    workspace_request_register = string_to_register(ws_name)
    client.write_registers(626, workspace_request_register)

def register_height_offset(height_offset):
    client.write_registers(620, int(number_to_raw_data(height_offset * 1000)))

def auto_calibration():
    print "Calibrate Robot if needed ..."
    client.write_register(311, 1)
    # Wait for end of calibration
    while client.read_input_registers(402, 1).registers[0] == 1:
        time.sleep(0.05)

def get_current_tool_id():
    return client.read_input_registers(200, count=1).registers[0]

def open_tool():
    tool_id = get_current_tool_id()
    if tool_id == 31:
        client.write_register(513, 1)
    else:
        client.write_register(510, 1)
    while client.read_holding_registers(150, count=1).registers[0] == 1:
        time.sleep(0.05)

# Function to call Modbus Server vision pick function
def vision_pick(workspace_str, height_offset, shape_int, color_int):
    register_workspace_name(workspace_str)
    register_height_offset(height_offset)

    client.write_registers(624, number_to_raw_data(shape_int))
    client.write_registers(625, number_to_raw_data(color_int))

    # launch vision pick function
    client.write_registers(612, 1)

    # Wait for end of function
    while client.read_holding_registers(150, count=1).registers[0] == 1:
        time.sleep(0.01)

    # - Check result : SHAPE AND COLOR
    result_shape_int = raw_data_to_number(client.read_holding_registers(159).registers[0])
    result_color_int = raw_data_to_number(client.read_holding_registers(160).registers[0])

    return result_shape_int, result_color_int

# ----- Main programm

if __name__ == '__main__':
    print "... START"

client = ModbusTcpClient('localhost', port=5020)

# ----- Variable definition
# To change
workspace_name = 'gazebo_1'
height_offset = 0.0

# connect to modbus server
client.connect()
print "Connected to modbus server"

# launch auto calibration then go to obs. pose
auto_calibration()
back_to_observation()

# update tool
client.write_registers(500, 1)

print 'VISION PICK - pick a red pawn, lift it and release it'
shape = ShapeEnum.ANY.value
color = ColorEnum.RED.value
shape_picked, color_picked = vision_pick(workspace_name, height_offset, shape, color)

# ---- Go to release pose
joints = [0.866, -0.24, -0.511, 0.249, -0.568, -0.016]
joints_to_send = list(map(lambda j: int(number_to_raw_data(j * 1000)), joints))

```

```
client.write_registers(0, joints_to_send)
client.write_register(100, 1)

# Wait for end of Move command
while client.read_holding_registers(150, count=1).registers[0] == 1:
    time.sleep(0.01)

open_tool()

back_to_observation()

# Activate learning mode and close connexion
client.write_register(300, 1)
client.close()
print "Close connection to modbus server"
print "--- END"
```

[Suggest a modification](#)[Download as PDF](#)

# PyNiryo Documentation



This documentation presents Ned's PyPi package, which is a TCP API made with Python.

It offers a simple way for developers to create programs for robot and to control them via remote communication from their computers. Contrary to the Python ROS Wrapper, the user will not need to be connected on the robot through a terminal.

## Note

This package is able to control Ned in simulation as well as the physical robot.



Ned

## Before getting started

If you haven't already done so, make sure to learn about the ROS robot software by reading [ROS documentation](#).

This documentation also contains everything you need to know if you want to use Ned through simulation.

## Sections organization

This document is organized in 4 main sections.

### Setup

Install & Setup your environment in order to use Ned with PyNiryo.

Firstly, follow [Installation instructions](#) (index.html#document-source/setup/installation), then [find your Robot IP address](#) (index.html#document-source/setup/ip\_address) to be ready.

## Installation

The library uses [Python](#), which must be installed and available in your working environment.

The version should be **equal or above**:

- 2.7 if you are using Python 2
- 3.6 if you are using Python 3

### 💡 Hint

To check your Python version, use the command `python --version` if you are using Python 2 and `python3 --version` if you are using Python3.

The below sections explain how to install the library with [pip](#), the package installer for Python.

### ⚠️ Attention

If you have both Python 2 & Python 3 installed on your computer, the command [pip](#) will install packages in Python 2 version. You should use [pip3](#) instead in order to target Python 3.

## Installation with pip

You need to install Numpy package beforehand:

```
pip install numpy
```

To install Ned's Python package via [pip](#) , simply execute:

```
pip install pyniryo
```

You can find more information about the PyPi package[here](#) (<https://pypi.org/project/pyniryo/>).

If you also want to use Vision functions to do your own Image Processing pipeline install OpenCV via the command:

```
pip install opencv-python
```

## Uninstall

To uninstall the library use:

```
pip uninstall pyniryo
```

## Find your Robot's IP address

In order to use your robot through TCP connection, you will first need to connect to it, which implies that you know its IP address.

The next sections explain how to find your robot IP according to your configuration:

- Hotspot mode
- Simulation or directly on the robot
- Direct ethernet connection
- Computer and robot connected on the same router
- Make IP permanent

## **Hotspot mode**

If you are directly connected to your robot through its wi-fi, the IP address you will need to use is **10.10.10.10**.

## **Simulation or directly on the robot**

In this situation, the robot is running on the same computer as the client, the IP address will be the localhost address **127.0.0.1**.

## **Direct ethernet connection**

If you are directly connected to your robot with an ethernet cable, the static IP of your robot will be **169.254.200.200**.

The reader should note that he may need to change his wired settings to allow the connection. See how to [Connect to Ned via Ethernet on Ubuntu](https://docs.niryo.com/applications/ned/source/tutorials/setup_connect_ned_ethernet.html) ([https://docs.niryo.com/applications/ned/source/tutorials/setup\\_connect\\_ned\\_ethernet.html](https://docs.niryo.com/applications/ned/source/tutorials/setup_connect_ned_ethernet.html)).

## **Computer and robot connected on the same router**

You will need to find the robot's address using **nmap**, or you can also use search button of Niryo Studio to see which robots are available.

You can also [Make IP permanent](#) so that you will not have to search for it next time

## **Make IP permanent**

### **Step 1**

Firstly, you need to be connected to your robot via SSH.

On Ubuntu, use the command line:

```
ssh niryo@<robot_ip_address>
```

The password is **robotics**.

On Windows, you can use [Putty](https://www.putty.org/) (<https://www.putty.org/>). Robot username is **niryo** and password is **robotics**.

### **Step 2**

Find your proxy key:

```
ifconfig
```

Your proxy key is written on the first line and should look something like **eth0**.

### Step 3

Select arbitrarily a number between 50 & 255. It will be your IP address' last number.

Then, edit the file **/etc/network/interfaces**:

```
sudo nano /etc/network/interfaces
```

And add to its end:

```
auto <robot_proxy_key>
iface <robot_proxy_key> inet static
    address 192.168.1.<your_ip_address_last_number>
    broadcast 192.168.1.255
    netmask 255.255.255.0
    gateway 192.168.1.1
    dns-nameservers 192.168.1.1
```

From its next reboot, the robot will appear under the IP **192.168.1.<your\_ip\_address\_last\_number>**.

## Verify your Setup and Get Started

In order to verify your computer's setup, we are going to run a program from it, and see if the robot answers as expected.

### **Note**

Before verifying your setup, be sure that your physical robot (or simulation) is turned on.

Firstly, go in the folder of your choice and create an empty file named "pyniryo\_test.py". This file will contain the checking code.

Edit this file and fill it with the following code:

```
from pyniryo import *
robot_ip_address = "10.10.10.10"
# Connect to robot & calibrate
robot = NiryoRobot(robot_ip_address)
robot.calibrate_auto()
# Move joints
robot.move_joints(0.0, 0.0, 0.0, 0.0, 0.0, 0.0, 0.0)
# Turn learning mode ON
robot.set_learning_mode(True)
# Stop TCP connection
robot.close_connection()
```

### **Attention**

Replace the third line with your Robot IP Address (index.html#document-source/setup/ip\_address) if you are not using Hotspot Mode.

Still on your computer, open a terminal, and place your current directory in the same folder than your file. Then, run the command:

```
python pyniryo_test.py
```

### **Note**

If you are using Python 3, you may need to change `python` to `python3`.

If your robot starts calibrating, then moves, and finally, goes to learning mode, your setup is validated, you can now start coding!

## Examples

Learn how to use the PyNiryo package to implement various tasks.

### Examples: Basics

In this file, two short programs are implemented & commented in order to help you understand the philosophy behind the PyNiryo package.

### **Danger**

If you are using the real robot, make sure the environment around it is clear.

### Your first move joint

The following example shows a first use case. It's a simple MoveJ.

```
from pyniryo import *
robot = NiryoRobot("10.10.10.10")
robot.calibrate_auto()
robot.move_joints(0.2, -0.3, 0.1, 0.0, 0.5, -0.8)
robot.close_connection()
```

### Code Details - First Move J

First of all, we import the library to be able to access functions.

```
from pyniryo import *
```

Then, we instantiate the connection and link the variable `robot` to the robot at the IP Address `10.10.10.10`.

```
robot = NiryoRobot("10.10.10.10")
```

Once the connection is done, we calibrate the robot using its **calibrate\_auto()** (index.html#api.tcp\_client.NiryoRobot.calibrate\_auto) function.

```
robot.calibrate_auto()
```

As the robot is now calibrated, we can do a Move Joints by giving the 6 axis positions in radians! To do so, we use **move\_joints()** (index.html#api.tcp\_client.NiryoRobot.move\_joints).

```
robot.move_joints(0.2, -0.3, 0.1, 0.0, 0.5, -0.8)
```

Our process is now over, we can close the connection with **quit()**.

```
robot.close_connection()
```

## Your first pick and place

In the second example, we are going to develop a pick and place algorithm.

```
from pyniryo import *
robot = NiryoRobot("10.10.10.10")
robot.calibrate_auto()
robot.update_tool()

robot.release_with_tool()
robot.move_pose(0.2, -0.1, 0.25, 0.0, 1.57, 0.0)
robot.grasp_with_tool()

robot.move_pose(0.2, 0.1, 0.25, 0.0, 1.57, 0.0)
robot.release_with_tool()

robot.close_connection()
```

### Code Details - First Pick And Place

First of all, we import the library and start the connection between our computer and the robot. We also calibrate the robot.

```
from pyniryo import *
robot = NiryoRobot("10.10.10.10")
robot.calibrate_auto()
```

Then, we equip the tool with **update\_tool()** (index.html#api.tcp\_client.NiryoRobot.update\_tool).

```
robot.update_tool()
```

Now that our initialization is done, we can open the gripper (or push air from the vacuum) with **release\_with\_tool()** (index.html#api.tcp\_client.NiryoRobot.release\_with\_tool), go to the picking pose via **move\_pose()** (index.html#api.tcp\_client.NiryoRobot.move\_pose) & then catch an object with

**grasp\_with\_tool()** (index.html#api.tcp\_client.NiryoRobot.grasp\_with\_tool)!

```
robot.release_with_tool()
robot.move_pose(0.2, -0.1, 0.25, 0.0, 1.57, 0.0)
robot.grasp_with_tool()
```

We now get to the place pose, and place the object.

```
robot.move_pose(0.2, 0.1, 0.25, 0.0, 1.57, 0.0)
robot.release_with_tool()
```

Our process is now over, we can close the connection.

```
robot.close_connection()
```

**Notes**

You may not have fully understood how to move the robot and use PyNiryo and that is totally fine because you will find more details on the next examples page!

The important thing to remember from this page is how to import the library, connect to the robot & call functions.

**Examples: Movement**

This document shows how to control Ned in order to make Move Joints & Move Pose.

If you want to see more, you can look at [PyNiryo - Joints & Pose](#) (index.html#joints-pose)

**Important**

In the following sections, you are supposed to be already connected to a calibrated robot. The robot's instance is saved in the variable `robot`. To know how to do so, go look at section [Examples: Basics](#) (index.html#document-source/examples/examples\_basics).

**Danger**

If you are using the real robot, make sure the environment around it is clear.

**Joints****Move Joints**

To make a move, you can either provide:

- 6 floats : `j1, j2, j3, j4, j5, j6`
- a list of 6 floats : `[j1, j2, j3, j4, j5, j6]`

It is possible to provide these parameters to the function **move\_joints()** (index.html#api.tcp\_client.NiryoRobot.move\_joints) or via the `joint` setter, at your convenience:

```
# Moving Joints with function & 6 floats
robot.move_joints(0.0, 0.0, 0.0, 0.0, 0.0, 0.0)

# Moving Joints with function & a list of floats
robot.move_joints([-0.5, -0.6, 0.0, 0.3, 0.0, 0.0])

# Moving Joints with setter & 6 floats
robot.joints = 0.2, -0.4, 0.0, 0.0, 0.0, 0.0

# Moving Joints with setter & a list of floats
robot.joints = [-0.2, 0.3, 0.2, 0.3, -0.6, 0.0]
```

You should note that these 4 commands are doing exactly the same thing! In your future scripts, chose the one you prefer, but try to remain consistent to keep a good readability.

## Get Joints

To get actual joint positions, you can use the function **get\_joints()** ([index.html#api.tcp\\_client.NiryoRobot.get\\_joints](#)) or the **joints** getter. Both will return a list of the 6 joints position:

```
# Getting Joints with function
joints_read = robot.get_joints()

# Getting Joints with getter
joints_read = robot.joints
```

### Hint

As we are developing in Python, we can unpack list very easily, which means that we can retrieve joints value in 6 variables by writing **j1, j2, j3, j4, j5, j6 = robot.get\_joints()** .

## Pose

### Move Pose

To perform a moveP, you can provide:

- 6 floats : x, y, z, roll, pitch, yaw
- a list of 6 floats : [x, y, z, roll, pitch, yaw]
- a **PoseObject** ([index.html#api.objects.PoseObject](#))

As for MoveJ, it is possible to provide these parameters to the function **move\_pose()** ([index.html#api.tcp\\_client.NiryoRobot.move\\_pose](#)) or the **pose** setter, at your convenience:

```
pose_target = [0.2, 0.0, 0.2, 0.0, 0.0, 0.0]
pose_target_obj = PoseObject(0.2, 0.0, 0.2, 0.0, 0.0, 0.0)

# Moving Pose with function
robot.move_pose(0.2, 0.0, 0.2, 0.0, 0.0, 0.0)
robot.move_pose(pose_target)
robot.move_pose(pose_target_obj)

# Moving Pose with setter
robot.pose = (0.2, 0.0, 0.2, 0.0, 0.0, 0.0)
robot.pose = pose_target
robot.pose = pose_target_obj
```

Each of these 6 commands are doing the same thing.

### Get Pose

To get end effector actual pose, you can use the function **get\_pose()** (`index.html#api.tcp_client.NiryoRobot.get_pose`) or the **pose** getter. Both will return a **PoseObject** (`index.html#api.objects.PoseObject`):

```
# Getting Joints with function
pose_read = robot.get_pose()

# Getting Joints with getter
pose_read = robot.pose
```

### How to use the PoseObject

The **PoseObject** (`index.html#api.objects.PoseObject`) is a Python object which allows to store all poses' 6 coordinates (x, y, z, roll, pitch, yaw) in one single instance. It can be converted into a list if needed with the method **to\_list()** (`index.html#api.objects.PoseObject.to_list`).

It also allows to create new **PoseObject** (`index.html#api.objects.PoseObject`) with some offset, much easier than copying list and editing only 1 or 2 values. For instance, imagine that we want to shift the place pose by 5 centimeters at each iteration of a loop, you can use the **copy\_with\_offsets()** (`index.html#api.objects.PoseObject.copy_with_offsets`) method:

```
pick_pose = PoseObject(
    x=0.30, y=0.0, z=0.15,
    roll=0, pitch=1.57, yaw=0.0
)
first_place_pose = PoseObject(
    x=0.0, y=0.2, z=0.15,
    roll=0, pitch=1.57, yaw=0.0
)
for i in range(5):
    robot.move_pose(pick_pose)
    new_place_pose = first_place_pose.copy_with_offsets(x_offset=0.05 * i)
    robot.move_pose(new_place_pose)
```

### Examples: Tool action

This page shows how to control Ned's tools.

If you want to see more, you can look at [PyNiryo - Tools](#) (`index.html#tools`)

#### Important

In this section, you are already supposed to be connected to a calibrated robot. The robot instance is saved in the variable **robot**.

#### Danger

If you are using the real robot, make sure the environment around it is clear.

### Tool control

## Equip Tool

In order to use a tool, it should be plugged mechanically to the robot but also connected to the software wise.

To do that, we should use the function **update\_tool()** (`index.html#api.tcp_client.NiryoRobot.update_tool`) which takes no argument. It will scan motor's connections and set the new tool!

The line to equip a new tool is:

```
robot.update_tool()
```

### Note

For the [Grasping](#) and [Releasing](#) sections, this command should be added in your codes! If you want to use a specific tool, you need to store the **ToolID** you are using in a variable named `tool_used`.

## Grasping

To grasp with any tool, you can use the function **grasp\_with\_tool()** (`index.html#api.tcp_client.NiryoRobot.grasp_with_tool`). This action corresponds to:

- Close gripper for Grippers
- Pull Air for Vacuum Pump
- Activate for Electromagnet

The line to grasp is:

```
robot.grasp_with_tool()
```

To grasp an object by specifying the tool:

```
if tool_used in [ToolID.GRIPPER_1, ToolID.GRIPPER_2, ToolID.GRIPPER_3]:
    robot.close_gripper(speed=500)
elif tool_used == ToolID.ELECTROMAGNET_1:
    pin_electromagnet = PinID.XXX
    robot.setup_electromagnet(pin_electromagnet)
    robot.activate_electromagnet(pin_electromagnet)
elif tool_used == ToolID.VACUUM_PUMP_1:
    robot.pull_air_vacuum_pump()
```

## Releasing

To release with any tool, you can use the function **release\_with\_tool()** (`index.html#api.tcp_client.NiryoRobot.release_with_tool`). This action corresponds to:

- Open gripper for Grippers
- Push Air for Vacuum pump
- Deactivate for Electromagnet

To release an object by specifying parameters:

```

if tool_used in [ToolID.GRIPPER_1, ToolID.GRIPPER_2, ToolID.GRIPPER_3]:
    robot.open_gripper(speed=500)
elif tool_used == ToolID.ELECTROMAGNET_1:
    pin_electromagnet = PinID.XXX
    robot.setup_electromagnet(pin_electromagnet)
    robot.deactivate_electromagnet(pin_electromagnet)
elif tool_used == ToolID.VACUUM_PUMP_1:
    robot.push_air_vacuum_pump()

```

## Pick & Place with tools

A Pick & Place is an action which consists in going to a certain pose in order to pick an object and then, going to another pose to place it.

This operation can be proceed as follows:

1. Going over the object with a certain offset to avoid collision;
2. Going down to the object's height;
3. Grasping with tool;
4. Going back to step 1's pose;
5. Going over the place pose with a certain offset to avoid collision;
6. Going down to place's height;
7. Releasing the object with tool;
8. Going back to step 5's pose.

There are plenty of ways to perform a pick and place with PyNiryo. Methods will be presented from the lowest to highest level.

## Code Baseline

For the sake of brevity, every piece of code beside the Pick & Place function will not be rewritten for every method. So that, you will need to use the code and implement the Pick & Place function to it.

```
# Imports
from pyniryo import *

tool_used = ToolID.XXX # Tool used for picking
robot_ip_address = "x.x.x.x" # Robot address

# The pick pose
pick_pose = PoseObject(
    x=0.25, y=0., z=0.15,
    roll=-0.0, pitch=1.57, yaw=0.0,
)
# The Place pose
place_pose = PoseObject(
    x=0.0, y=-0.25, z=0.1,
    roll=0.0, pitch=1.57, yaw=-1.57)

def pick_n_place_version_x(robot):
    # ----- #
    # -- CODE GOES HERE -- #
    # ----- #

    if __name__ == '__main__':
        # Connect to robot
        client = NiryoRobot(robot_ip_address)
        # Calibrate robot if robot needs calibration
        client.calibrate_auto()
        # Changing tool
        client.update_tool()

    pick_n_place_version_x(client)

    # Releasing connection
    client.close_connection()
```

## First Solution: the heaviest

For this first function, every steps are done by hand, as well as poses computing.

### Note

In this example, the tool used is a Gripper. If you want to use another tool than a gripper, do not forget to adapt grasp & release functions!

```
def pick_n_place_version_1(robot):
    height_offset = 0.05 # Offset according to Z-Axis to go over pick & place poses
    gripper_speed = 400

    # Going Over Object
    robot.move_pose(pick_pose.x, pick_pose.y, pick_pose.z + height_offset,
                    pick_pose.roll, pick_pose.pitch, pick_pose.yaw)
    # Opening Gripper
    robot.open_gripper(gripper_speed)
    # Going to picking place and closing gripper
    robot.move_pose(place_pose)
    robot.close_gripper(gripper_speed)

    # Raising
    robot.move_pose(pick_pose.x, pick_pose.y, pick_pose.z + height_offset,
                    pick_pose.roll, pick_pose.pitch, pick_pose.yaw)

    # Going Over Place pose
    robot.move_pose(place_pose.x, place_pose.y, place_pose.z + height_offset,
                    place_pose.roll, place_pose.pitch, place_pose.yaw)
    # Going to Place pose
    robot.move_pose(place_pose)
    # Opening Gripper
    robot.open_gripper(gripper_speed)
    # Raising
    robot.move_pose(place_pose.x, place_pose.y, place_pose.z + height_offset,
                    place_pose.roll, place_pose.pitch, place_pose.yaw)
```

## Second Solution: Use of PoseObject

For the second solution, we use a **PoseObject** (index.html#api.objects.PoseObject) in order to calculate approach poses more easily.

### Note

To see more about **PoseObject** (index.html#api.objects.PoseObject), go look at [PoseObject dedicated section](#) (index.html#how-to-use-the-poseobject)

```
def pick_n_place_version_2(robot):
    height_offset = 0.05 # Offset according to Z-Axis to go over pick & place poses

    pick_pose_high = pick_pose.copy_with_offsets(z_offset=height_offset)
    place_pose_high = place_pose.copy_with_offsets(z_offset=height_offset)

    # Going Over Object
    robot.move_pose(pick_pose_high)
    # Opening Gripper
    robot.release_with_tool()
    # Going to picking place and closing gripper
    robot.move_pose(pick_pose)
    robot.grasp_with_tool()
    # Raising
    robot.move_pose(pick_pose_high)

    # Going Over Place pose
    robot.move_pose(place_pose_high)
    # Going to Place pose
    robot.move_pose(place_pose)
    # Opening Gripper
    robot.release_with_tool(gripper_speed)
    # Raising
    robot.move_pose(place_pose_high)
```

## Third Solution: Pick from pose & Place from pose functions

For those who have already seen the API Documentation, you may have seen pick & place dedicated functions!

In this example, we use **pick\_from\_pose()** (index.html#api.tcp\_client.NiryoRobot.pick\_from\_pose) and **place\_from\_pose()** (index.html#api.tcp\_client.NiryoRobot.place\_from\_pose) in order to split our function in only 2 commands!

```
def pick_n_place_version_3(robot):
    # Pick
    robot.pick_from_pose(pick_pose)
    # Place
    robot.place_from_pose(place_pose)
```

## Fourth Solution: All in one

The example exposed in the previous section could be useful if you want to do an action between the pick & the place phases.

For those who want to do everything in one command, you can use the **pick\_and\_place()** (index.html#api.tcp\_client.NiryoRobot.pick\_and\_place) function!

```
def pick_n_place_version_4(robot):
    # Pick & Place
    robot.pick_and_place(pick_pose, place_pose)
```

## Examples: Conveyor Belt

This document shows how to use Ned's Conveyor Belt.

If you want to see more about Ned's Conveyor Belt functions, you can look at [PyNiryo - Conveyor](#) (index.html#conveyor)

### Danger

If you are using the real robot, make sure the environment around it is clear.

## Simple Conveyor control

This short example shows how to connect a conveyor and launch its motor (control it by setting its speed and direction):

```
from pyniryo import *

# Connecting to robot
robot = NiryoRobot(<robot_ip_address>)

# Activating connexion with Conveyor Belt
conveyor_id = robot.set_conveyor()

# Running the Conveyor Belt at 50% of its maximum speed, in forward direction
robot.run_conveyor(conveyor_id, speed=50, direction=ConveyorDirection.FORWARD)

# Waiting 3 seconds
robot.wait(3)

# Stopping robot's motor
robot.stop_conveyor(conveyor_id)

# Deactivating connexion with the Conveyor Belt
robot.unset_conveyor(conveyor_id)
```

## Advanced Conveyor Belt control

This example shows how to do a certain amount of pick & place by using the Conveyor Belt with the infrared sensor:

```

from pyniryo import *

# -- Setting variables
sensor_pin_id = PinID.GPIO_1A

catch_nb = 5

# The pick pose
pick_pose = PoseObject(
    x=0.25, y=0., z=0.15,
    roll=-0., pitch=1.57, yaw=0.0,
)
# The Place pose
place_pose = PoseObject(
    x=0., y=-0.25, z=0.1,
    roll=0., pitch=1.57, yaw=-1.57)

# -- MAIN PROGRAM

# Connecting to the robot
robot = NiryoRobot(<robot_ip_address>)

# Activating connexion with the Conveyor Belt
conveyor_id = robot.set_conveyor()

for i in range(catch_nb):
    robot.run_conveyor(conveyor_id)
    while robot.digital_read(sensor_pin_id) == PinState.LOW:
        robot.wait(0.1)

    # Stopping robot's motor
    robot.stop_conveyor(conveyor_id)
    # Making a pick & place
    robot.pick_and_place(pick_pose, place_pose)

# Deactivating connexion with the Conveyor Belt
robot.unset_conveyor(conveyor_id)

```

## Examples: Vision

This page shows how to use Ned's Vision Set.

If you want to see more about Ned's Vision functions, you can look at [PyNiryo - Vision](#)  
If you want to see how to do Image Processing, go check out the [Image Processing section](#).

### **Note**

Even if you do not own a Vision Set, you can still realize these examples with the Gazebo simulation version.

### **Danger**

If you are using the real robot, make sure the environment around it is clear.

## Needed piece of code

### **Important**

In order to achieve the following examples, you need to create a vision workspace. In this page, the workspace used is named `workspace_1`. To create it, the user should go on Niryo Studio!

As the examples start always the same, add the following lines at the beginning of codes:

```

# Imports
from pyniryo import *

# - Constants
workspace_name = "workspace_1" # Robot's Workspace Name
robot_ip_address = "x.x.x.x"

# The pose from where the image processing happens
observation_pose = PoseObject(
    x=0.16, y=0.0, z=0.35,
    roll=0.0, pitch=1.57, yaw=0.0,
)
# Place pose
place_pose = PoseObject(
    x=0.0, y=-0.2, z=0.12,
    roll=0.0, pitch=1.57, yaw=-1.57
)

# - Initialization

# Connect to robot
robot = NiryoRobot(robot_ip_address)
# Calibrate robot if the robot needs calibration
robot.calibrate_auto()
# Updating tool
robot.update_tool()

# --- ----- #
# --- CODE GOES HERE --- #
# --- ----- #

robot.close_connection()

```

### Hint

All the following examples are only a part of what can be made with the API in terms of Vision. We advise you to look at [API - Vision](#) (index.html#vision) to understand more deeply

## Simple Vision Pick & Place

The goal of a Vision Pick & Place is the same as a classical Pick & Place, with a close difference: the camera detects where the robot has to go in order to pick!

This short example shows how to do your first Vision pick using the **vision\_pick()** (index.html#api.tcp\_client.NiryoRobot.vision\_pick) function:

```

robot.move_pose(observation_pose)
# Trying to pick target using camera
obj_found, shape_ret, color_ret = robot.vision_pick(workspace_name)
if obj_found:
    robot.place_from_pose(place_pose)

robot.set_learning_mode(True)

```

### Code Details - Simple Vision Pick and Place

To execute a Vision pick, we firstly need to go to a place where the robot will be able to see the workspace:

```
robot.move_pose(observation_pose)
```

Then, we try to perform a Vision pick in the workspace with the **vision\_pick()** (index.html#api.tcp\_client.NiryoRobot.vision\_pick) function:

```
obj_found, shape_ret, color_ret = robot.vision_pick(workspace_name)
```

Variables `shape_ret` and `color_ret` are respectively of type **ObjectShape** and **ObjectColor**, and store the shape and the color of the detected object! We will not use them for this first example.

The `obj_found` variable is a boolean which indicates whereas an object has been found and picked, or not. Thus, if the pick worked, we can place the object at the place pose.

```
if obj_found:  
    robot.place_from_pose(place_pose)
```

Finally, we turn learning mode on:

```
robot.set_learning_mode(True)
```

### **Note**

If your `obj_found` variable indicates `False`, check that:

- Nothing obstructs the camera field of view
- Workspace's 4 markers are visible
- At least 1 object is placed fully inside the workspace

## First conditioning via Vision

In most of use cases, the robot will need to perform more than one Pick & Place. In this example, we will see how to condition multiple objects according to a straight line:

```
# Initializing variables  
offset_size = 0.05  
max_catch_count = 4  
  
# Loop until enough objects have been caught  
catch_count = 0  
while catch_count < max_catch_count:  
    # Moving to observation pose  
    robot.move_pose(observation_pose)  
  
    # Trying to get object via Vision Pick  
    obj_found, shape, color = robot.vision_pick(workspace_name)  
    if not obj_found:  
        robot.wait(0.1)  
        continue  
  
    # Calculate place pose and going to place the object  
    next_place_pose = place_pose.copy_with_offsets(x_offset=catch_count * offset_size)  
    robot.place_from_pose(next_place_pose)  
  
    catch_count += 1  
  
robot.go_to_sleep()
```

### Code Details - First Conditioning via Vision

We want to catch `max_catch_count` objects, and space each of them by `offset_size` meter:

```
offset_size = 0.05
max_catch_count = 4
```

We start a loop until the robot has caught `max_catch_count` objects:

```
catch_count = 0
while catch_count < max_catch_count:
```

For each iteration, we firstly go to the observation pose and then, try to make a Vision pick in the workspace:

```
robot.move_pose(observation_pose)
obj_found, shape, color = robot.vision_pick(workspace_name)
```

If the Vision pick failed, we wait 0.1 second and then, start a new iteration:

```
if not obj_found:
    robot.wait(0.1)
    continue
```

Else, we compute the new place position according to the number of catches, and then, go placing the object at that place:

```
next_place_pose = place_pose.copy_with_offsets(x_offset=catch_count * offset_size)
robot.place_from_pose(next_place_pose)
```

We also increment the `catch_count` variable:

```
catch_count += 1
```

Once the target catch number is achieved, we go to sleep:

```
robot.go_to_sleep()
```

## Multi Reference Conditioning

During a conditioning task, objects may not always be placed as the same place according to their type. In this example, we will see how to align object according to their color, using the color element **ObjectColor** returned by **vision\_pick()** ([index.html#api.tcp\\_client.NiryoRobot.vision\\_pick](#)) function:

```

# Distance between elements
offset_size = 0.05
max_failure_count = 3

# Dict to write catch history
count_dict = {
    ObjectColor.BLUE: 0,
    ObjectColor.RED: 0,
    ObjectColor.GREEN: 0,
}

try_without_success = 0
# Loop until too much failures
while try_without_success < max_failure_count:
    # Moving to observation pose
    robot.move_pose(observation_pose)
    # Trying to get object via Vision Pick
    obj_found, shape, color = robot.vision_pick(workspace_name)
    if not obj_found:
        try_without_success += 1
        robot.wait(0.1)
        continue

    # Choose X position according to how the color line is filled
    offset_x_ind = count_dict[color]

    # Choose Y position according to ObjectColor
    if color == ObjectColor.BLUE:
        offset_y_ind = -1
    elif color == ObjectColor.RED:
        offset_y_ind = 0
    else:
        offset_y_ind = 1

    # Going to place the object
    next_place_pose = place_pose.copy_with_offsets(x_offset=offset_x_ind * offset_size,
                                                    y_offset=offset_y_ind * offset_size)
    robot.place_from_pose(next_place_pose)

    # Increment count
    count_dict[color] += 1
    try_without_success = 0

robot.go_to_sleep()

```

## Code Details - Multi Reference Conditioning

We want to catch objects until Vision Pick failed `max_failure_count` times. Each of the object will be put on a specific column according to its color. The number of catches for each color will be stored on a dictionary `count_dict`.

```

# Distance between elements
offset_size = 0.05
max_failure_count = 3

# Dict to write catch history
count_dict = {
    ObjectColor.BLUE: 0,
    ObjectColor.RED: 0,
    ObjectColor.GREEN: 0,
}

try_without_success = 0
# Loop until too much failures
while try_without_success < max_failure_count:

```

For each iteration, we firstly go to the observation pose and then, try to make a Vision pick in the workspace:

```
robot.move_pose(observation_pose)
obj_found, shape, color = robot.vision_pick(workspace_name)
```

If the Vision pick failed, we wait 0.1 second and then, start a new iteration, without forgetting to increment the failure counter:

```
if not obj_found:
    try_without_success += 1
    robot.wait(0.1)
    continue
```

Else, we compute the new place position according to the number of catches, and then, go place the object at that place:

```
# Choose X position according to how the color line is filled
offset_x_ind = count_dict[color]

# Choose Y position according to ObjectColor
if color == ObjectColor.BLUE:
    offset_y_ind = -1
elif color == ObjectColor.RED:
    offset_y_ind = 0
else:
    offset_y_ind = 1

# Going to place the object
next_place_pose = place_pose.copy_with_offsets(x_offset=offset_x_ind * offset_size,
                                                y_offset=offset_y_ind * offset_size)
robot.place_from_pose(next_place_pose)
```

We increment the `count_dict` dictionary and reset `try_without_success`:

```
count_dict[color] += 1
try_without_success = 0
```

Once the target catch number is achieved, we go to sleep:

```
robot.go_to_sleep()
```

## Sorting Pick with Conveyor

An interesting way to bring objects to the robot, is the use of a Conveyor Belt. In this examples, we will see how to catch only a certain type of object by stopping the conveyor as soon as the object is detected on the workspace.

```

# Initializing variables
offset_size = 0.05
max_catch_count = 4
shape_expected = ObjectShape.CIRCLE
color_expected = ObjectColor.RED

conveyor_id = robot.set_conveyor()

catch_count = 0
while catch_count < max_catch_count:
    # Turning conveyor on
    robot.run_conveyor(conveyor_id)
    # Moving to observation pose
    robot.move_pose(observation_pose)
    # Check if object is in the workspace
    obj_found, pos_array, shape, color = robot.detect_object(workspace_name,
        shape=shape_expected,
        color=color_expected)

    if not obj_found:
        robot.wait(0.5) # Wait to let the conveyor turn a bit
        continue

    # Stopping conveyor
    robot.stop_conveyor(conveyor_id)
    # Making a vision pick
    obj_found, shape, color = robot.vision_pick(workspace_name,
        shape=shape_expected,
        color=color_expected)

    if not obj_found: # If visual pick did not work
        continue

    # Calculate place pose and going to place the object
    next_place_pose = place_pose.copy_with_offsets(x_offset=catch_count * offset_size)
    robot.place_from_pose(next_place_pose)

    catch_count += 1

    # Stopping & unsetting conveyor
    robot.stop_conveyor(conveyor_id)
    robot.unset_conveyor(conveyor_id)

robot.go_to_sleep()

```

## Code Details - Sort Picking

Firstly, we initialize your process: we want the robot to catch 4 red circles. To do so, we set variables `shape_expected` and `color_expected` with `ObjectShape.CIRCLE` and `ObjectColor.RED`.

```

offset_size = 0.05
max_catch_count = 4
shape_expected = ObjectShape.CIRCLE
color_expected = ObjectColor.RED

```

We activate the connection with the Conveyor Belt and start a loop until the robot has caught `max_catch_count` objects:

```

conveyor_id = robot.set_conveyor()

catch_count = 0
while catch_count < max_catch_count:

```

For each iteration, we firstly run the Conveyor Belt (if the latter is already running, nothing will happen), then go to the observation pose:

```

# Turning the Conveyor Belt on
robot.run_conveyor(conveyor_id)
# Moving to observation pose
robot.move_pose(observation_pose)

```

We then check if an object corresponding to our criteria is in the workspace. If not, we wait 0.5 second and then, start a new iteration:

```
obj_found, pos_array, shape, color = robot.detect_object(workspace_name,
                                                       shape=shape_expected,
                                                       color=color_expected)
if not obj_found:
    robot.wait(0.5) # Wait to let the conveyor turn a bit
    continue
```

Else, stop the Conveyor Belt and try to make a Vision pick:

```
# Stopping Conveyor Belt
robot.stop_conveyor(conveyor_id)
# Making a Vision pick
obj_found, shape, color = robot.vision_pick(workspace_name,
                                              shape=shape_expected,
                                              color=color_expected)
if not obj_found: # If visual pick did not work
    continue
```

If Vision Pick succeed, compute new place pose, and place the object:

```
# Calculate place pose and going to place the object
next_place_pose = place_pose.copy_with_offsets(x_offset=catch_count * offset_size)
robot.place_from_pose(next_place_pose)

catch_count += 1
```

Once the target catch number is achieved, we stop the Conveyor Belt and go to sleep:

```
# Stopping & unsetting Conveyor Belt
robot.stop_conveyor(conveyor_id)
robot.unset_conveyor(conveyor_id)

robot.go_to_sleep()
```

## Code templates

As code structures are always the same, we wrote down few templates for you to start your code file with a good form.

### The short template

Very simple, straightforward:

```
from pyniryo import *

# Connect to robot & calibrate
robot = NiryoRobot(<robot_ip_address>)
robot.calibrate_auto()

# --- ----- #
# --- YOUR CODE --- #
# --- ----- #

# Releasing connection
robot.close_connection()
```

### Advanced template

This template let the user define his own process but it handles connection, calibration, tool equipping, and makes the robot go to sleep at the end:

```
from pyniryo import *

local_mode = False # Or True
tool_used = ToolID.GRIPPER_1
# Set robot address
robot_ip_address_rpi = "x.x.x.x"
robot_ip_address_local = "127.0.0.1"

robot_ip_address = robot_ip_address_local if local_mode else robot_ip_address_rpi

def process(niryo_edu):
    # --- ----- #
    # --- YOUR CODE --- #
    # --- ----- #

if __name__ == '__main__':
    # Connect to robot
    robot = NiryoRobot(robot_ip_address)
    # Calibrate robot if robot needs calibration
    robot.calibrate_auto()
    # Equip tool
    robot.update_tool()
    # Launching main process
    process(client)
    # Ending
    robot.go_to_sleep()
    # Releasing connection
    robot.close_connection()
```

## Advanced template for Conveyor Belt

Same as [Advanced template](#) but with a Conveyor Belt

```
from pyniryo import *

# Set robot address
robot_ip_address = "x.x.x.x"

def process(robot, conveyor_id):
    robot.run_conveyor(conveyor_id)

    # --- ----- #
    # --- YOUR CODE --- #
    # --- ----- #

    robot.stop_conveyor()

if __name__ == '__main__':
    # Connect to robot
    robot = NiryoRobot(robot_ip_address)
    # Calibrate robot if robot needs calibration
    robot.calibrate_auto()
    # Equip tool
    robot.update_tool()
    # Activating connexion with conveyor
    conveyor_id = robot.set_conveyor()
    # Launching main process
    process(robot, conveyor_id)
    # Ending
    robot.go_to_sleep()
    # Deactivating connexion with conveyor
    robot.unset_conveyor(conveyor_id)
    # Releasing connection
    robot.close_connection()
```

## Advanced template for Vision

Huge template for Vision users!

```
from pyniryo import *

local_mode = False # Or True
workspace_name = "workspace_1" # Robot's Workspace Name
# Set robot address
robot_ip_address_rpi = "x.x.x.x"
robot_ip_address_local = "127.0.0.1"

robot_ip_address = robot_ip_address_local if local_mode else robot_ip_address_rpi

# The pose from where the image processing happens
observation_pose = PoseObject(
    x=0.18, y=0.0, z=0.35,
    roll=0.0, pitch=1.57, yaw=-0.2,
)

# Center of the conditioning area
place_pose = PoseObject(
    x=0.0, y=-0.23, z=0.12,
    roll=0.0, pitch=1.57, yaw=-1.57
)

def process(robot):
    robot.move_pose(observation_pose)
    catch_count = 0
    while catch_count < 3:
        ret = robot.get_target_pose_from_cam(workspace_name,
                                              height_offset=0.0,
                                              shape=ObjectShape.ANY,
                                              color=ObjectColor.ANY)
        obj_found, obj_pose, shape, color = ret
        if not obj_found:
            continue
        catch_count += 1
        # --- ----- #
        # --- YOUR CODE --- #
        # --- ----- #
        robot.place_from_pose(place_pose)

if __name__ == '__main__':
    # Connect to robot
    robot = NiryoRobot(robot_ip_address)
    # Calibrate robot if robot needs calibration
    robot.calibrate_auto()
    # Equip tool
    robot.update_tool()
    # Launching main process
    process(client)
    # Ending
    robot.go_to_sleep()
    # Releasing connection
    robot.close_connection()
```

## API Documentation

Master controls with PyNiryo with full the detailed functions [here](#) ([index.html#document-source/api\\_doc/api](#)).

Discover also [Vision Functions](#) ([index.html#document-source/vision/image\\_processing\\_overview](#)) to create your own image processing pipelines!

## PyNiryo API Documentation

This file presents the different [Command functions, Enums & Python object classes](#) available with the API.

- [Command functions](#) are used to deal directly with the robot. It could be **`move_joints()`**, **`get_hardware_status()`**, **`vision_pick()`**, or also **`run_conveyor()`**

- Enums are used to pass specific arguments to functions. For instance **PinState**, **ConveyorDirection** , ...
- Python object classes, as **PoseObject** , ease some operations

## Command functions

This section references all existing functions to control your robot, which includes:

- Moving the robot
- Using Vision
- Controlling Conveyor Belts
- Playing with Hardware

All functions to control the robot are accessible via an instance of the class **NiryoRobot**

```
robot = NiryoRobot(<robot_ip_address>)
```

See examples on [Examples: Basics](#) (index.html#examples-basics)

List of functions subsections:

- [TCP Connection](#)
- [Main purpose functions](#)
- [Joints & Pose](#)
- [Saved Poses](#)
- [Pick & Place](#)
- [Trajectories](#)
- [Tools](#)
- [Hardware](#)
- [Conveyor](#)
- [Vision](#)

### TCP Connection

#### **NiryoRobot.connect(ip\_address)**

Connect to the TCP Server

**Parameters:**

**ip\_address** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – IP Address

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None)

#### **NiryoRobot.close\_connection()**

Close connection with robot

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None)

### Main purpose functions

#### **NiryoRobot.calibrate(calibrate\_mode)**

Calibrate (manually or automatically) motors. Automatic calibration will do nothing if motors are already calibrated

**Parameters:**

**calibrate\_mode** (*CalibrateMode*) – Auto or Manual

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.calibrate\_auto()

Start a automatic motors calibration if motors are not calibrated yet

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.need\_calibration()

Return a bool indicating whereas the robot motors need to be calibrate

**Return type:**

[bool](https://docs.python.org/3/library/functions.html#bool) (<https://docs.python.org/3/library/functions.html#bool>)

### NiryoRobot.get\_learning\_mode()

Get learning mode state

**Returns:**

**True** if learning mode is on

**Return type:**

[bool](https://docs.python.org/3/library/functions.html#bool) (<https://docs.python.org/3/library/functions.html#bool>)

### NiryoRobot.set\_learning\_mode(*enabled*)

Set learning mode if param is **True** , else turn it off

**Parameters:**

**enabled** ([bool](https://docs.python.org/3/library/functions.html#bool)) (<https://docs.python.org/3/library/functions.html#bool>) – **True** or **False**

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.set\_arm\_max\_velocity(*percentage\_speed*)

Limit arm max velocity to a percentage of its maximum velocity

**Parameters:**

**percentage\_speed** ([int](https://docs.python.org/3/library/functions.html#int)) (<https://docs.python.org/3/library/functions.html#int>) – Should be between 1 & 100

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.set\_jog\_control(*enabled*)

Set jog control mode if param is True, else turn it off

**Parameters:**

**enabled** ([bool](https://docs.python.org/3/library/functions.html#bool)) (<https://docs.python.org/3/library/functions.html#bool>) – **True** or **False**

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### **static NiryoRobot.wait(duration)**

Wait for a certain time

**Parameters:**

**duration** ([float](https://docs.python.org/3/library/functions.html#float)) (<https://docs.python.org/3/library/functions.html#float>) – duration in seconds

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

## Joins & Pose

### **NiryoRobot.get\_joints()**

Get joints value in radians You can also use a getter

```
joints = robot.get_joints()
joints = robot.joints
```

**Returns:**

List of joints value

**Return type:**

[list](https://docs.python.org/3/library/stdtypes.html#list) (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)]

### **NiryoRobot.get\_pose()**

Get end effector link pose as [x, y, z, roll, pitch, yaw]. x, y & z are expressed in meters / roll, pitch & yaw are expressed in radians You can also use a getter

```
pose = robot.get_pose()
pose = robot.pose
```

**Return type:**

[PoseObject](#) ([index.html#api.objects.PoseObject](#))

### **NiryoRobot.get\_pose\_quat()**

Get end effector link pose in Quaternion coordinates

**Returns:**

Position and quaternion coordinates concatenated in a list : [x, y, z, qx, qy, qz, qw]

**Return type:**

[list](https://docs.python.org/3/library/stdtypes.html#list) (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)]

### **NiryoRobot.move\_joints(\*args)**

Move robot joints. Joints are expressed in radians.

All lines of the next example realize the same operation:

```
robot.joints = [0.2, 0.1, 0.3, 0.0, 0.5, 0.0]
robot.move_joints([0.2, 0.1, 0.3, 0.0, 0.5, 0.0])
robot.move_joints(0.2, 0.1, 0.3, 0.0, 0.5, 0.0)
```

**Parameters:**

**args** (*Union[list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]***float** (<https://docs.python.org/3/library/functions.html#float>)), (<https://docs.python.org/3/library/stdtypes.html#tuple>)**float** (<https://docs.python.org/3/library/functions.html#float>)]]) – either 6 args (1 for each joints) or a list of 6 joints

**Return type:**

**None** (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.move\_pose(\*args)**

Move robot end effector pose to a (x, y, z, roll, pitch, yaw) pose. x, y & z are expressed in meters / roll, pitch & yaw are expressed in radians

All lines of the next example realize the same operation:

```
robot.pose = [0.2, 0.1, 0.3, 0.0, 0.5, 0.0]
robot.move_pose([0.2, 0.1, 0.3, 0.0, 0.5, 0.0])
robot.move_pose(0.2, 0.1, 0.3, 0.0, 0.5, 0.0)
robot.move_pose(PoseObject(0.2, 0.1, 0.3, 0.0, 0.5, 0.0))
```

**Parameters:**

**args** (*Union[tuple* (<https://docs.python.org/3/library/stdtypes.html#tuple>)*]***float** (<https://docs.python.org/3/library/functions.html#float>)), (<https://docs.python.org/3/library/stdtypes.html#list>)**float** (<https://docs.python.org/3/library/functions.html#float>), (<https://index.html#api.objects.PoseObject>)) – either 6 args (1 for each coordinates) or a list of 6 coordinates or a **PoseObject**

**Return type:**

**None** (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.move\_linear\_pose(\*args)**

Move robot end effector pose to a (x, y, z, roll, pitch, yaw) pose with a linear trajectory

**Parameters:**

**args** (*Union[tuple* (<https://docs.python.org/3/library/stdtypes.html#tuple>)*]***float** (<https://docs.python.org/3/library/functions.html#float>)), (<https://docs.python.org/3/library/stdtypes.html#list>)**float** (<https://docs.python.org/3/library/functions.html#float>)), (<https://index.html#api.objects.PoseObject>)) – either 6 args (1 for each coordinates) or a list of 6 coordinates or a PoseObject

**Return type:**

**None** (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.shift\_pose(axis, shift\_value)**

Shift robot end effector pose along one axis

**Parameters:**

- **axis** (*RobotAxis*) – Axis along which the robot is shifted
- **shift\_value** (*float* (<https://docs.python.org/3/library/functions.html#float>)) – In meter for X/Y/Z and radians for roll/pitch/yaw

**Return type:**

**None** (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.shift\_linear\_pose(axis, shift\_value)**

Shift robot end effector pose along one axis, with a linear trajectory

**Parameters:**

- **axis** (*RobotAxis*) – Axis along which the robot is shifted
- **shift\_value** ([float](https://docs.python.org/3/library/functions.html#float)) – In meter for X/Y/Z and radians for roll/pitch/yaw

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None)

**NiryoRobot.jog\_joints(\*args)**

Jog robot joints'. Jog corresponds to a shift without motion planning. Values are expressed in radians.

**Parameters:**

**args** (*Union[list* ([float](https://docs.python.org/3/library/stdtypes.html#list))*, tuple* ([float](https://docs.python.org/3/library/functions.html#float))*, tuple* ([float](https://docs.python.org/3/library/stdtypes.html#tuple))*] – either 6 args (1 for each joints) or a list of 6 joints offset*

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None)

**NiryoRobot.jog\_pose(\*args)**

Jog robot end effector pose. Jog corresponds to a shift without motion planning. Arguments are [dx, dy, dz, d\_roll, d\_pitch, d\_yaw] dx, dy & dz are expressed in meters / d\_roll, d\_pitch & d\_yaw are expressed in radians

**Parameters:**

**args** (*Union[list* ([float](https://docs.python.org/3/library/stdtypes.html#list))*, tuple* ([float](https://docs.python.org/3/library/functions.html#float))*, tuple* ([float](https://docs.python.org/3/library/stdtypes.html#tuple))*] – either 6 args (1 for each coordinates) or a list of 6 offset*

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None)

**NiryoRobot.move\_to\_home\_pose()**

Move to a position where the forearm lays on shoulder

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None)

**NiryoRobot.go\_to\_sleep()**

Go to home pose and activate learning mode

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None)

**NiryoRobot.forward\_kinematics(\*args)**

Compute forward kinematics of a given joints configuration and give the associated spatial pose

**Parameters:**

**args** (*Union[list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]***float** (<https://docs.python.org/3/library/functions.html#float>),  
*tuple* (<https://docs.python.org/3/library/stdtypes.html#tuple>)*]***float** (<https://docs.python.org/3/library/functions.html#float>)]*]*) – either 6 args (1 for each joints) or a list of 6 joints

**Return type:**

**PoseObject** ([index.html#api.objects.PoseObject](#))

**NiryoRobot.inverse\_kinematics(\*args)**

Compute inverse kinematics

**Parameters:**

**args** (*Union[tuple* (<https://docs.python.org/3/library/stdtypes.html#tuple>)*]***float** (<https://docs.python.org/3/library/functions.html#float>),  
*list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]***float** (<https://docs.python.org/3/library/functions.html#float>),  
*PoseObject* ([index.html#api.objects.PoseObject](#))]*]*) – either 6 args (1 for each coordinates) or a list of 6 coordinates or a **PoseObject**

**Returns:**

List of joints value

**Return type:**

**list** (<https://docs.python.org/3/library/stdtypes.html#list>)*]***float** (<https://docs.python.org/3/library/functions.html#float>)]

**Saved Poses****NiryoRobot.get\_pose\_saved(pose\_name)**

Get pose saved in from Ned's memory

**Parameters:**

**pose\_name** (*str* (<https://docs.python.org/3/library/stdtypes.html#str>)) – Pose name in robot's memory

**Returns:**

Pose associated to pose\_name

**Return type:**

**PoseObject** ([index.html#api.objects.PoseObject](#))

**NiryoRobot.save\_pose(pose\_name, \*args)**

Save pose in robot's memory

**Parameters:**

**args** (*Union[list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]***float** (<https://docs.python.org/3/library/functions.html#float>),  
*tuple* (<https://docs.python.org/3/library/stdtypes.html#tuple>)*]***float** (<https://docs.python.org/3/library/functions.html#float>),  
*PoseObject* ([index.html#api.objects.PoseObject](#))]*]*) – either 6 args (1 for each coordinates) or a list of 6 coordinates or a PoseObject

**Return type:**

**None** (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.delete\_pose(pose\_name)**

Delete pose from robot's memory

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.get\_saved\_pose\_list()**

Get list of poses' name saved in robot memory

**Return type:**

[list](https://docs.python.org/3/library/stdtypes.html#list) (<https://docs.python.org/3/library/stdtypes.html#list>) [[str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>)]

**Pick & Place****NiryoRobot.pick\_from\_pose(\*args)**

Execute a picking from a pose.

A picking is described as :

- \* going over the object
- \* going down until height = z
- \* grasping with tool
- \* going back over the object

**Parameters:**

**args** (*Union[*[list](https://docs.python.org/3/library/functions.html#list) (<https://docs.python.org/3/library/functions.html#list>), [tuple](https://docs.python.org/3/library/stdtypes.html#tuple) (<https://docs.python.org/3/library/stdtypes.html#tuple>), [float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>), [float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>), [PoseObject](https://docs.python.org/3/library/api.objects.PoseObject) (<https://docs.python.org/3/library/api.objects.PoseObject>)]*)* – either 6 args (1 for each coordinates) or a list of 6 coordinates or a PoseObject

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.place\_from\_pose(\*args)**

Execute a placing from a position.

A placing is described as :

- \* going over the place
- \* going down until height = z
- \* releasing the object with tool
- \* going back over the place

**Parameters:**

**args** (*Union[*[list](https://docs.python.org/3/library/functions.html#list) (<https://docs.python.org/3/library/functions.html#list>), [tuple](https://docs.python.org/3/library/stdtypes.html#tuple) (<https://docs.python.org/3/library/stdtypes.html#tuple>), [float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>), [float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>), [PoseObject](https://docs.python.org/3/library/api.objects.PoseObject) (<https://docs.python.org/3/library/api.objects.PoseObject>)]*)* – either 6 args (1 for each coordinates) or a list of 6 coordinates or a PoseObject

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.pick\_and\_place(*pick\_pose*, *place\_pos*, *dist\_smoothing*=0.0)

Execute a pick then a place

**Parameters:**

- **pick\_pose** (*Union[list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]float* (<https://docs.python.org/3/library/functions.html#float>),  
[index.html#api.objects.PoseObject](https://index.html#api.objects.PoseObject)]) – Pick Pose : [x, y, z, roll, pitch, yaw] or PoseObject
- **place\_pos** (*Union[list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]float* (<https://docs.python.org/3/library/functions.html#float>),  
[index.html#api.objects.PoseObject](https://index.html#api.objects.PoseObject)]) – Place Pose : [x, y, z, roll, pitch, yaw] or PoseObject
- **dist\_smoothing** (*float* (<https://docs.python.org/3/library/functions.html#float>)) – Distance from waypoints before smoothing trajectory

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

## Trajectories

### NiryoRobot.get\_trajectory\_saved(*trajectory\_name*)

Get trajectory saved in Ned's memory

**Returns:**

Trajectory

**Return type:**

*list* (<https://docs.python.org/3/library/stdtypes.html#list>)*[list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]float* (<https://docs.python.org/3/library/functions.html#float>)]

### NiryoRobot.execute\_trajectory\_from\_poses(*list\_poses*, *dist\_smoothing*=0.0)

Execute trajectory from list of poses

**Parameters:**

- **list\_poses** (*list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]float* (<https://docs.python.org/3/library/functions.html#float>])) – List of [x,y,z,qx,qy,qz,qw] or list of [x,y,z,roll,pitch,yaw]
- **dist\_smoothing** (*float* (<https://docs.python.org/3/library/functions.html#float>)) – Distance from waypoints before smoothing trajectory

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.execute\_trajectory\_from\_poses\_and\_joints(*list\_pose\_joints*, *list\_type*=None, *dist\_smoothing*=0.0)

Execute trajectory from list of poses and joints

**Parameters:**

- **list\_pose\_joints** (*list* (<https://docs.python.org/3/library/stdtypes.html#list>)*]float* (<https://docs.python.org/3/library/functions.html#float>))) – List of [x,y,z,qx,qy,qz,qw] or list of [x,y,z,roll,pitch,yaw] or a list of [j1,j2,j3,j4,j5,j6]

- **list\_type** ([list](https://docs.python.org/3/library/stdtypes.html#list) ([string](https://docs.python.org/3/library/stdtypes.html#string)) – List of string ‘pose’ or ‘joint’, or [‘pose’] (if poses only) or [‘joint’] (if joints only). If None, it is assumed there are only poses in the list.
- **dist\_smoothing** ([float](https://docs.python.org/3/library/functions.html#float) ([float](https://docs.python.org/3/library/functions.html#float)) – Distance from waypoints before smoothing trajectory

**Return type:**[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)**NiryoRobot.execute\_trajectory\_saved(trajectory\_name)**

Execute trajectory from Ned’s memory

**Return type:**[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)**NiryoRobot.save\_trajectory(trajectory\_name, list\_poses)**

Save trajectory in robot memory

**Parameters:**

**list\_poses** ([list](https://docs.python.org/3/library/stdtypes.html#list) (<https://docs.python.org/3/library/stdtypes.html#list>)[\[float\]](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – List of [x,y,z,qx,qy,qz,qw] or list of [x,y,z,roll,pitch,yaw]

**Return type:**[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)**NiryoRobot.delete\_trajectory(trajectory\_name)**

Delete trajectory from robot’s memory

**Return type:**[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)**NiryoRobot.get\_saved\_trajectory\_list()**

Get list of trajectories’ name saved in robot memory

**Return type:**

**list** ([list](https://docs.python.org/3/library/stdtypes.html#list) (<https://docs.python.org/3/library/stdtypes.html#list>)[\[str\]](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>))

**Tools****NiryoRobot.get\_current\_tool\_id()**

Get equipped tool Id

**Return type:**

ToolID

**NiryoRobot.update\_tool()**

Update equipped tool

**Return type:**[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.grasp\_with\_tool()**

Grasp with tool | This action correspond to | - Close gripper for Grippers | - Pull Air for Vacuum pump | - Activate for Electromagnet

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.release\_with\_tool()**

Release with tool | This action correspond to | - Open gripper for Grippers | - Push Air for Vacuum pump | - Deactivate for Electromagnet

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.open\_gripper(*speed=500*)**

Open gripper associated to 'gripper\_id' with a speed 'speed'

**Parameters:**

**speed** ([int](https://docs.python.org/3/library/functions.html#int)) (<https://docs.python.org/3/library/functions.html#int>) – Between 100 & 1000

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.close\_gripper(*speed=500*)**

Close gripper associated to 'gripper\_id' with a speed 'speed'

**Parameters:**

**speed** ([int](https://docs.python.org/3/library/functions.html#int)) (<https://docs.python.org/3/library/functions.html#int>) – Between 100 & 1000

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.pull\_air\_vacuum\_pump()**

Pull air of vacuum pump

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.push\_air\_vacuum\_pump()**

Push air of vacuum pump

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.setup\_electromagnet(*pin\_id*)**

Setup electromagnet on pin

**Parameters:**

**pin\_id** (*PinID*) –

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.activate\_electromagnet(*pin\_id*)**

Activate electromagnet associated to electromagnet\_id on pin\_id

**Parameters:**

**pin\_id** (*PinID*) –

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.deactivate\_electromagnet(pin\_id)

Deactivate electromagnet associated to electromagnet\_id on pin\_id

**Parameters:**

**pin\_id** (*PinID*) –

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.enable\_tcp(enable=True)

Enables or disables the TCP function (Tool Center Point). If activation is requested, the last recorded TCP value will be applied. The default value depends on the gripper equipped. If deactivation is requested, the TCP will be coincident with the tool\_link.

**Parameters:**

**enable** (*Bool*) – True to enable, False otherwise.

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.set\_tcp(\*args)

Activates the TCP function (Tool Center Point) and defines the transformation between the tool\_link frame and the TCP frame.

**Parameters:**

**args** (*Union[list* (<https://docs.python.org/3/library/stdtypes.html#list>)*, float* (<https://docs.python.org/3/library/functions.html#float>), *tuple* (<https://docs.python.org/3/library/stdtypes.html#tuple>), *float* (<https://docs.python.org/3/library/functions.html#float>), *PoseObject* ([index.html#api.objects.PoseObject](#))]) – either 6 args (1 for each coordinates) or a list of 6 coordinates or a PoseObject

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.reset\_tcp()

Reset the TCP (Tool Center Point) transformation. The TCP will be reset according to the tool equipped.

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.tool\_reboot()

Reboot the motor of the tool equipped. Useful when an Overload error occurs. (cf HardwareStatus)

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

## Hardware

### NiryoRobot.set\_pin\_mode(pin\_id, pin\_mode)

Set pin number pin\_id to mode pin\_mode

**Parameters:**

- **pin\_id** (*PinID*) –
- **pin\_mode** (*PinMode*) –

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.digital\_write(pin\_id, digital\_state)

Set pin\_id state to digital\_state

**Parameters:**

- **pin\_id** (*PinID*) –
- **digital\_state** (*PinState*) –

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.digital\_read(pin\_id)

Read pin number pin\_id and return its state

**Parameters:**

**pin\_id** (*PinID*) –

**Return type:**

PinState

### NiryoRobot.get\_hardware\_status()

Get hardware status : Temperature, Hardware version, motors names & types ...

**Returns:**

Infos contains in a [HardwareStatusObject](#)

**Return type:**

[HardwareStatusObject](#) ([index.html#api.objects.HardwareStatusObject](#))

### NiryoRobot.get\_digital\_io\_state()

Get Digital IO state : Names, modes, states

**Returns:**

List of [DigitalPinObject](#) instance

**Return type:**

**list** (<https://docs.python.org/3/library/stdtypes.html#list>) [[DigitalPinObject](#) ([index.html#api.objects.DigitalPinObject](#))]

## Conveyor

### NiryoRobot.set\_conveyor()

Activate a new conveyor and return its ID

**Returns:**

New conveyor ID

**Return type:**

ConveyorID

### NiryoRobot.unset\_conveyor(conveyor\_id)

Remove specific conveyor.

**Parameters:**

**conveyor\_id** (ConveyorID) – Basically, ConveyorID.ONE or ConveyorID.TWO

### NiryoRobot.run\_conveyor(conveyor\_id, speed=50, direction=<ConveyorDirection.FORWARD: 1>)

Run conveyor at id 'conveyor\_id'

**Parameters:**

- **conveyor\_id** (ConveyorID) –
- **speed** ([int](https://docs.python.org/3/library/functions.html#int)) –
- **direction** (ConveyorDirection) –

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.stop\_conveyor(conveyor\_id)

Stop conveyor at id 'conveyor\_id'

**Parameters:**

**conveyor\_id** (ConveyorID) –

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.control\_conveyor(conveyor\_id, control\_on, speed, direction)

Control conveyor at id 'conveyor\_id'

**Parameters:**

- **conveyor\_id** (ConveyorID) –
- **control\_on** ([bool](https://docs.python.org/3/library/functions.html#bool)) (<https://docs.python.org/3/library/functions.html#bool>) –
- **speed** ([int](https://docs.python.org/3/library/functions.html#int)) (<https://docs.python.org/3/library/functions.html#int>) – New speed which is a percentage of maximum speed
- **direction** (ConveyorDirection) – Conveyor direction

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.get\_connected\_conveyors\_id()

**Returns:**

List of the connected conveyors' ID

**Return type:**

[list](https://docs.python.org/3/library/stdtypes.html#list) (<https://docs.python.org/3/library/stdtypes.html#list>)[ConveyorID]

## Vision

### NiryoRobot.get\_img\_compressed()

Get image from video stream in a compressed format. Use **uncompress\_image** from the vision package to uncompress it

**Returns:**

string containing a JPEG compressed image

**Return type:**

[str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>)

### NiryoRobot.set\_brightness(brightness\_factor)

Modify video stream brightness

**Parameters:**

**brightness\_factor** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – How much to adjust the brightness. 0.5 will give a darkened image, 1 will give the original image while 2 will enhance the brightness by a factor of 2.

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.set\_contrast(contrast\_factor)

Modify video stream contrast

**Parameters:**

**contrast\_factor** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – While a factor of 1 gives original image. Making the factor towards 0 makes the image greyer, while factor>1 increases the contrast of the image.

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.set\_saturation(saturation\_factor)

Modify video stream saturation

**Parameters:**

**saturation\_factor** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – How much to adjust the saturation. 0 will give a black and white image, 1 will give the original image while 2 will enhance the saturation by a factor of 2.

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

### NiryoRobot.get\_image\_parameters()

Get last stream image parameters: Brightness factor, Contrast factor, Saturation factor.

Brightness factor: How much to adjust the brightness. 0.5 will give a darkened image, 1 will give the original image while 2 will enhance the brightness by a factor of 2.

Contrast factor: A factor of 1 gives original image. Making the factor towards 0 makes the image greyer, while factor>1 increases the contrast of the image.

Saturation factor: 0 will give a black and white image, 1 will give the original image while 2 will enhance the saturation by a factor of 2.

**Returns:**

Brightness factor, Contrast factor, Saturation factor

**Return type:**

<b>float</b> ( <a href="https://docs.python.org/3/library/functions.html#float">https://docs.python.org/3/library/functions.html#float</a> ), <a href="https://docs.python.org/3/library/functions.html#float">https://docs.python.org/3/library/functions.html#float</a> , <a href="https://docs.python.org/3/library/functions.html#float">https://docs.python.org/3/library/functions.html#float</a>	<b>float</b> <b>float</b>
---	------------------------------

### NiryoRobot.get\_target\_pose\_from\_rel(workspace\_name, height\_offset, x\_rel, y\_rel, yaw\_rel)

Given a pose (x\_rel, y\_rel, yaw\_rel) relative to a workspace, this function returns the robot pose in which the current tool will be able to pick an object at this pose.

The height\_offset argument (in m) defines how high the tool will hover over the workspace. If height\_offset = 0, the tool will nearly touch the workspace.

**Parameters:**

- **workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>)) – name of the workspace
- **height\_offset** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – offset between the workspace and the target height
- **x\_rel** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – x relative pose (between 0 and 1)
- **y\_rel** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – y relative pose (between 0 and 1)
- **yaw\_rel** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – Angle in radians

**Returns:**

target\_pose

**Return type:**

[PoseObject](#) ([index.html#api.objects.PoseObject](#))

### NiryoRobot.get\_target\_pose\_from\_cam(workspace\_name, height\_offset=0.0, shape=<ObjectShape.ANY: 'ANY'>, color=<ObjectColor.ANY: 'ANY'>)

First detects the specified object using the camera and then returns the robot pose in which the object can be picked with the current tool

**Parameters:**

- **workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str) (<https://docs.python.org/3/library/stdtypes.html#str>)) – name of the workspace
- **height\_offset** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – offset between the workspace and the target height
- **shape** ([ObjectShape](#)) – shape of the target
- **color** ([ObjectColor](#)) – color of the target

**Returns:**

object\_found, object\_pose, object\_shape, object\_color

**Return type:**

<b>(bool</b> ( <a href="https://docs.python.org/3/library/functions.html#bool">https://docs.python.org/3/library/functions.html#bool</a> ), <a href="#">index.html#api.objects.PoseObject</a> ), <a href="#">ObjectShape</a> , <a href="#">ObjectColor</a> )	<b>PoseObject</b>
---	-------------------

**NiryoRobot.vision\_pick(workspace\_name, height\_offset=0.0, shape=<ObjectShape.ANY: 'ANY'>, color=<ObjectColor.ANY: 'ANY'>)**

Picks the specified object from the workspace. This function has multiple phases:

1. detect object using the camera
2. prepare the current tool for picking
3. approach the object
4. move down to the correct picking pose
5. actuate the current tool
6. lift the object

Example:

```
robot = NiryoRobot(ip_address="x.x.x.x")
robot.calibrate_auto()
robot.move_pose(<observation_pose>)
obj_found, shape_ret, color_ret = robot.vision_pick(<workspace_name>,
                                                    height_offset=0.0,
                                                    shape=ObjectShape.ANY,
                                                    color=ObjectColor.ANY)
```

**Parameters:**

- **workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – name of the workspace
- **height\_offset** ([float](https://docs.python.org/3/library/functions.html#float)) – offset between the workspace and the target height
- **shape** (*ObjectShape*) – shape of the target
- **color** (*ObjectColor*) – color of the target

**Returns:**

object\_found, object\_shape, object\_color

**Return type:**

([bool](https://docs.python.org/3/library/functions.html#bool)) ([bool](https://docs.python.org/3/library/functions.html#bool)), *ObjectShape*, *ObjectColor*)

**NiryoRobot.move\_to\_object(workspace\_name, height\_offset, shape, color)**

Same as *get\_target\_pose\_from\_cam* but directly moves to this position

**Parameters:**

- **workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – name of the workspace
- **height\_offset** ([float](https://docs.python.org/3/library/functions.html#float)) – offset between the workspace and the target height
- **shape** (*ObjectShape*) – shape of the target
- **color** (*ObjectColor*) – color of the target

**Returns:**

object\_found, object\_shape, object\_color

**Return type:**

([bool](https://docs.python.org/3/library/functions.html#bool)) ([bool](https://docs.python.org/3/library/functions.html#bool)), *ObjectShape*, *ObjectColor*)

**NiryoRobot.detect\_object(workspace\_name, shape=<ObjectShape.ANY: 'ANY'>, color=<ObjectColor.ANY: 'ANY'>)**

Detect object in workspace and return its pose and characteristics

**Parameters:**

- **workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – name of the workspace
- **shape** ([ObjectShape](#)) – shape of the target
- **color** ([ObjectColor](#)) – color of the target

**Returns:**

object\_found, object\_pose, object\_shape, object\_color

**Return type:**

[\(bool\)](#) (<https://docs.python.org/3/library/functions.html#bool>), [PoseObject](#)  
 ([index.html#api.objects.PoseObject](#)), [str](#) (<https://docs.python.org/3/library/stdtypes.html#str>), [str](#)  
 (<https://docs.python.org/3/library/stdtypes.html#str>)

**NiryoRobot.get\_camera\_intrinsics()**

Get calibration object: camera intrinsics, distortions coefficients

**Returns:**

camera intrinsics, distortions coefficients

**Return type:**

[\(list\)](#) (<https://docs.python.org/3/library/stdtypes.html#list>) [[list](#)  
 (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](#)  
 (<https://docs.python.org/3/library/functions.html#float>)] ] , [list](#)  
 (<https://docs.python.org/3/library/stdtypes.html#list>) [[list](#)  
 (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](#)  
 (<https://docs.python.org/3/library/functions.html#float>)]])

**NiryoRobot.save\_workspace\_from\_robot\_poses(workspace\_name, pose\_origin, pose\_2, pose\_3, pose\_4)**

Save workspace by giving the poses of the robot to point its 4 corners with the calibration Tip. Corners should be in the good order. Markers' pose will be deduced from these poses

Poses should be either a list [x, y, z, roll, pitch, yaw] or a PoseObject

**Parameters:**

- **workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – workspace name, maximum lenght 30 char.
- **pose\_origin** ([Union\[list\]](#) (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](#)  
 (<https://docs.python.org/3/library/functions.html#float>)], [PoseObject](#)  
 ([index.html#api.objects.PoseObject](#))) –
- **pose\_2** ([Union\[list\]](#) (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](#)  
 (<https://docs.python.org/3/library/functions.html#float>)], [PoseObject](#)  
 ([index.html#api.objects.PoseObject](#))) –
- **pose\_3** ([Union\[list\]](#) (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](#)  
 (<https://docs.python.org/3/library/functions.html#float>)], [PoseObject](#)  
 ([index.html#api.objects.PoseObject](#))) –
- **pose\_4** ([Union\[list\]](#) (<https://docs.python.org/3/library/stdtypes.html#list>) [[float](#)  
 (<https://docs.python.org/3/library/functions.html#float>)], [PoseObject](#)  
 ([index.html#api.objects.PoseObject](#))) –

**Return type:**

[None](#) (<https://docs.python.org/3/library/constants.html#None>)

**NiryoRobot.save\_workspace\_from\_points(workspace\_name, point\_origin, point\_2, point\_3, point\_4)**

Save workspace by giving the points of workspace's 4 corners. Points are written as [x, y, z] Corners should be in the good order.

#### Parameters:

- **workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – workspace name, maximum lenght 30 char.
- **point\_origin** ([list](https://docs.python.org/3/library/stdtypes.html#list) ([float](https://docs.python.org/3/library/functions.html#float) ([float](https://docs.python.org/3/library/functions.html#float)))) –
- **point\_2** ([list](https://docs.python.org/3/library/stdtypes.html#list) ([float](https://docs.python.org/3/library/functions.html#float) ([float](https://docs.python.org/3/library/functions.html#float)))) –
- **point\_3** ([list](https://docs.python.org/3/library/stdtypes.html#list) ([float](https://docs.python.org/3/library/functions.html#float) ([float](https://docs.python.org/3/library/functions.html#float)))) –
- **point\_4** ([list](https://docs.python.org/3/library/stdtypes.html#list) ([float](https://docs.python.org/3/library/functions.html#float) ([float](https://docs.python.org/3/library/functions.html#float)))) –

#### Return type:

[None](https://docs.python.org/3/library/constants.html#None) ([None](https://docs.python.org/3/library/constants.html#None))

### NiryoRobot.delete\_workspace(workspace\_name)

Delete workspace from robot's memory

#### Parameters:

**workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) –

#### Return type:

[None](https://docs.python.org/3/library/constants.html#None) ([None](https://docs.python.org/3/library/constants.html#None))

### NiryoRobot.get\_workspace\_ratio(workspace\_name)

Get workspace ratio from robot's memory

#### Parameters:

**workspace\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) –

#### Return type:

[float](https://docs.python.org/3/library/functions.html#float) ([float](https://docs.python.org/3/library/functions.html#float))

### NiryoRobot.get\_workspace\_list()

Get list of workspaces' name store in robot's memory

#### Return type:

[list](https://docs.python.org/3/library/stdtypes.html#list) ([str](https://docs.python.org/3/library/stdtypes.html#str) ([str](https://docs.python.org/3/library/stdtypes.html#str)))

## Enums

Enums are used to pass specific parameters to functions.

For instance, **shift\_pose()** will need a parameter from **RobotAxis** enum

```
robot.shift_pose(RobotAxis.Y, 0.15)
```

List of enums:

- **CalibrateMode**

- **RobotAxis**
- **ToolID**
- **PinMode**
- **PinState**
- **PinID**
- **ConveyorID**
- **ConveyorDirection**
- **ObjectColor**
- **ObjectShape**

**`class CalibrateMode(value)`**

Enumeration of Calibration Modes

**AUTO= 0**

**MANUAL= 1**

**`class RobotAxis(value)`**

Enumeration of Robot Axis : it used for Shift command

**X= 0**

**Y= 1**

**Z= 2**

**ROLL= 3**

**PITCH= 4**

**YAW= 5**

**`class ToolID(value)`**

Enumeration of Tools IDs

**NONE= 0**

**GRIPPER\_1= 11**

**GRIPPER\_2= 12**

**GRIPPER\_3= 13**

**ELECTROMAGNET\_1= 30**

**VACUUM\_PUMP\_1= 31**

**`class PinMode(value)`**

## Enumeration of Pin Modes

**INPUT= 0****OUTPUT= 1*****class PinState(value)***

Pin State is either LOW or HIGH

**LOW= 0****HIGH= 1*****class PinID(value)***

Enumeration of Robot Pins

**GPIO\_1A= 0****GPIO\_1B= 1****GPIO\_1C= 2****GPIO\_2A= 3****GPIO\_2B= 4****GPIO\_2C= 5*****class ConveyorID(value)***

Enumeration of Conveyor IDs used for Conveyor control

**NONE= 0****ID\_1= 12****ID\_2= 13*****class ConveyorDirection(value)***

Enumeration of Conveyor Directions used for Conveyor control

**FORWARD= 1****BACKWARD= -1*****class ObjectColor(value)***

Enumeration of Colors available for image processing

**RED= 'RED'**

**BLUE= 'BLUE'****GREEN= 'GREEN'****ANY= 'ANY'****`class ObjectShape(value)`**

Enumeration of Shapes available for image processing

**SQUARE= 'SQUARE'****CIRCLE= 'CIRCLE'****ANY= 'ANY'****Python object classes****Special objects****`class PoseObject(x, y, z, roll, pitch, yaw)`**

Pose object which stores x, y, z, roll, pitch &amp; yaw parameters

**`copy_with_offsets(x_offset=0.0, y_offset=0.0, z_offset=0.0, roll_offset=0.0, pitch_offset=0.0, yaw_offset=0.0)`**

Create a new pose from copying from copying actual pose with offsets

**Return type:**`PoseObject` ([index.html#api.objects.PoseObject](#))**`to_list()`**

Return a list [x, y, z, roll, pitch, yaw] corresponding to the pose's parameters

**Return type:**

`list` (<https://docs.python.org/3/library/stdtypes.html#list>) [`float` (<https://docs.python.org/3/library/functions.html#float>)]

**`class HardwareStatusObject(rpi_temperature, hardware_version, connection_up, error_message, calibration_needed, calibration_in_progress, motor_names, motor_types, motors_temperature, motors_voltage, hardware_errors)`**

Object used to store every hardware information

**`class DigitalPinObject(pin_id, name, mode, state)`**

Object used to store information on digital pins

**Start with Image Processing**

Discover how to create your own image processing pipelines!

**Overview & examples**

This file illustrates few image processing pipeline using vision module from niryo\_edu package. This module is based in [OpenCV](https://opencv.org/) (<https://opencv.org/>) and its functions are detailed in [Functions documentation](#) (index.html#document-source/vision/image\_processing\_api).

The package niryo\_edu comes up with the **vision** module which contains image processing functions including thresholding, blob detection, ...

To use it, add to your imports `from pyniryo.vision import *`.

#### Note

It is also possible to merge both import lines by using `from pyniryo import *`.

## Play with Robot video stream

We are firstly going to take a look at robot's functions which can be find at [API - Vision](#) (index.html#vision)

### Get & display image from stream

Ned can share its video stream through TCP. As sending raw images will lead to heavy packets which can saturate the network, it sends compressed images. You access it through the robot's function: **get\_img\_compressed()** (index.html#api.tcp\_client.NiryoRobot.get\_img\_compressed). Once your image is received, you firstly need to uncompress via **uncompress\_image()** (index.html#vision.image\_functions.uncompress\_image) and you can then display it with **show\_img\_and\_wait\_close()** (index.html#vision.image\_functions.show\_img\_and\_wait\_close).

```
from pyniryo import *

# Connecting to robot
robot = NiryoRobot("10.10.10.10")

# Getting image
img_compressed = robot.get_img_compressed()
# Uncompressing image
img = uncompress_image(img_compressed)

# Displaying
show_img_and_wait_close("img_stream", img)
```

#### Note

**show\_img\_and\_wait\_close()** (index.html#vision.image\_functions.show\_img\_and\_wait\_close) will wait for the user to press either Q or Esc key, before closing the window.

## Undistort and display video stream

In this section, we are going to display the raw video stream & the undistorted video stream.

As Ned's camera is passing raw images to the robot, these images are distorted due to the camera lens. In order to undistort them, we need to use Ned's camera intrinsics.

To undistort the raw image, we use **undistort\_image()** (index.html#vision.image\_functions.undistort\_image) which needs to be called with the parameters given by Ned through **get\_camera\_intrinsics()** (index.html#api.tcp\_client.NiryoRobot.get\_camera\_intrinsics).

Once, we have both raw & undistorted images, we can concatenate them in order to display them in once with **concat\_imgs()** (index.html#vision.image\_functions.concat\_imgs). Finally, we display the image **show\_img()** (index.html#vision.image\_functions.show\_img).

```
from pyniryo import *

observation_pose = PoseObject(
    x=0.18, y=0.0, z=0.35,
    roll=0.0, pitch=1.57, yaw=-0.2,
)

# Connecting to robot
robot = NiryoRobot("10.10.10.10")
robot.calibrate_auto()

# Getting calibration param
mtx, dist = robot.get_camera_intrinsics()
# Moving to observation pose
robot.move_pose(observation_pose)

while "User do not press Escape neither Q":
    # Getting image
    img_compressed = robot.get_img_compressed()
    # Uncompressing image
    img_raw = uncompress_image(img_compressed)
    # Undistorting
    img_undistort = undistort_image(img_raw, mtx, dist)

    # - Display
    # Concatenating raw image and undistorted image
    concat_ims = concat_imgs((img_raw, img_undistort))

    # Showing images
    key = show_img("Images raw & undistorted", concat_ims, wait_ms=30)
    if key in [27, ord("q")]: # Will break loop if the user press Escape or Q
        break
```


### **Note**

To see more about camera distortion/undistortion, go on [OpenCV Documentation about Camera Calibration](#)

([https://docs.opencv.org/2.4/modules/calib3d/doc/camera\\_calibration\\_and\\_3d\\_reconstruction.html](https://docs.opencv.org/2.4/modules/calib3d/doc/camera_calibration_and_3d_reconstruction.html)).

## Pure image processing functions

In order to illustrate functions, we are going to use the following image.



### **Attention**

In this section it is supposed that:

- You have imported `pyniryo.vision`
- The variable `img` is containing the image on which image processing is applied

## Color thresholding

Color thresholding is very useful in order to detect object with an uniform color. The implemented function to realize this operation is `threshold_hsv()` ([index.html#vision.image\\_functions.threshold\\_hsv](#)).

The following code is using parameters from `ColorHSV` ([index.html#vision.enums.ColorHSV](#)) enum in order to threshold Red features & *hand made* parameters to extract Blue:

```
img_threshold_red = threshold_hsv(img_test, *ColorHSV.RED.value)


blue_min_hsv = [90, 85, 70]
blue_max_hsv = [125, 255, 255]


img_threshold_blue = threshold_hsv(img_test, list_min_hsv=blue_min_hsv,
                                   list_max_hsv=blue_max_hsv, reverse_hue=False)

show_img("img_threshold_red", img_threshold_red)

show_img_and_wait_close("img_threshold_blue", img_threshold_blue)
```

*Images result*

Thresh color	Image result
Blue	

Thresh color	Image result
Red	

## Morphological transformations

Morphological transformations are some simple operations based on the image shape. It is normally performed on binary images. It needs two inputs, one is our original image, second one is called structuring element or kernel which decides the nature of operation. Two basic morphological operators are **Erosion** ([https://en.wikipedia.org/wiki/Mathematical\\_morphology#Erosion](https://en.wikipedia.org/wiki/Mathematical_morphology#Erosion)) and **Dilation** ([https://en.wikipedia.org/wiki/Mathematical\\_morphology#Dilation](https://en.wikipedia.org/wiki/Mathematical_morphology#Dilation)).

Then its variant forms like **Opening** ([https://en.wikipedia.org/wiki/Mathematical\\_morphology#Opening](https://en.wikipedia.org/wiki/Mathematical_morphology#Opening)), **Closing** ([https://en.wikipedia.org/wiki/Mathematical\\_morphology#Closing](https://en.wikipedia.org/wiki/Mathematical_morphology#Closing)) also comes into play. Learn more on [Wikipedia page](https://en.wikipedia.org/wiki/Mathematical_morphology) ([https://en.wikipedia.org/wiki/Mathematical\\_morphology](https://en.wikipedia.org/wiki/Mathematical_morphology)).

The implemented function to realize these operations is **morphological\_transformations()** ([index.html#vision.image\\_functions.morphological\\_transformations](#)). It uses **MorphoType** ([index.html#vision.enums.MorphoType](#)) and **KernelType** ([index.html#vision.enums.KernelType](#)) to determine which operation should be applied on the image.

The code shows how to do a Closing & an Erosion:

```
img_threshold = threshold_hsv(img_test, *ColorHSV.ANY.value)



img_close = morphological_transformations(img_threshold, morpho_type=MorphoType.CLOSE,
                                         kernel_shape=(11, 11), kernel_type=KernelType.ELLIPSE)

img_erode = morphological_transformations(img_threshold, morpho_type=MorphoType.ERODE,
                                         kernel_shape=(9, 9), kernel_type=KernelType.RECT)

show_img("img_threshold", img_threshold)
show_img("img_erode", img_erode)
show_img_and_wait_close("img_close", img_close)
```

*Images result*

Morpho type	Image result
-------------	--------------

Morpho type	Image result
None	
Erode	
Close	

### Contours finder

Contours can be explained simply as a curve joining all the continuous points (along the boundary), having same color or intensity. The contours are a useful tool for shape analysis and object detection and recognition.

See more

on [OpenCV](#)

[Documentation](#)

([https://docs.opencv.org/3.4/d3/d05/tutorial\\_py\\_table\\_of\\_contents\\_contours.html](https://docs.opencv.org/3.4/d3/d05/tutorial_py_table_of_contents_contours.html)).

The implemented function to realize these operations is **biggest\_contours\_finder()** ([index.html#vision.image\\_functions.biggest\\_contours\\_finder](#)) which takes a Black & White image, and extracts the biggest (in term of area) contours from it.

The code to extract and draw the 3 biggest contours from an image is the following:


```
img_threshold = threshold_hsv(img_test, *ColorHSV.ANY.value)
img_threshold = morphological_transformations(img_threshold, morpho_type=MorphoType.OPEN,
                                              kernel_shape=(11, 11), kernel_type=KernelType.ELLIPSE)

cnts = biggest_contours_finder(img_threshold, 3)

img_contours = draw_contours(img_threshold, cnts)

show_img("init", img_threshold)
show_img_and_wait_close("img with contours", img_contours)
```

*Images result*



## Find object center position


In order to catch an object, we need to find a pose from where the end effector can grasp the object. The following method uses contours which have been found in the previous section and finds their barycenter and orientation via the functions **get\_contour\_barycenter()** (index.html#vision.image\_functions.get\_contour\_barycenter) & **get\_contour\_angle()** (index.html#vision.image\_functions.get\_contour\_angle).

```
img_threshold = threshold_hsv(img_test, *ColorHSV.ANY.value)
img_threshold = morphological_transformations(img_threshold, morpho_type=MorphoType.OPEN,
                                              kernel_shape=(11, 11), kernel_type=KernelType.ELLIPSE)

cnt = biggest_contour_finder(img_threshold)

cnt_barycenter = get_contour_barycenter(cnt)
cnt_angle = get_contour_angle(cnt)
```

*Images result*



#### **Note**


The drawn vector is normal to the contour's length because we want Ned to catch the object by the width rather than the length. Indeed, it leads to least cases where the gripper cannot open enough.

#### **Markers extraction**

As image processing happens in a workspace, it is important to extract the workspace beforehand! To do so, you can use the function **extract\_img\_workspace()** (index.html#vision.image\_functions.extract\_img\_workspace).

```
status, im_work = extract_img_workspace(img, workspace_ratio=1.0)
show_img("init", img_test)
show_img_and_wait_close("img_workspace", img_workspace)
```

*Images result*



### Debug mode




If Ned's functions are failing, you can use Debug functions which are **debug\_threshold\_color()** (index.html#vision.image\_functions.debug\_threshold\_color) & **debug\_markers()** (index.html#vision.image\_functions.debug\_markers) in order to display what the robot sees.

You can use the functions as follow:

```
debug_color = debug_threshold_color(img_test, ColorHSV.RED)
_status, debug_markers_im = debug_markers(img_test, workspace_ratio=1.0)

show_img("init", img_test)
show_img("debug_color", debug_color)
show_img_and_wait_close("debug_markers", debug_markers_im)
```

*Images result*

Original	
Debug red	
Debug Markers	

## Do your own image processing!

Now that you are a master in image processing, let's look at full examples.

### Display video stream with extracted workspace

In the current state, the following code will display the video stream and the extracted workspace image. You can add your own image processing functions maybe to apply additional process.

```

from pyniryo import *

# Connecting to robot
robot = NiryoRobot("10.10.10.10")
robot.calibrate_auto()

# Getting calibration param
mtx, dist = robot.get_camera_intrinsics()
# Moving to observation pose
robot.move_pose(*observation_pose.to_list())

while "User do not press Escape neither Q":
    # Getting image
    img_compressed = robot.get_img_compressed()
    # Uncompressing image
    img_raw = uncompress_image(img_compressed)
    # Undistorting
    img_undistort = undistort_image(img_raw, mtx, dist)
    # Trying to find markers
    workspace_found, res_img_markers = debug_markers(img_undistort)
    # Trying to extract workspace if possible
    if workspace_found:
        img_workspace = extract_img_workspace(img_undistort, workspace_ratio=1.0)
    else:
        img_workspace = None

    # --- ----- #
    # --- YOUR CODE --- #
    # --- ----- #

    # - Display
    # Concatenating raw image and undistorted image
    concat_ims = concat_imgs((img_raw, img_undistort))
    # Concatenating extracted workspace with markers annotation
    if img_workspace is not None:
        resized_img_workspace = resize_img(img_workspace, height=res_img_markers.shape[0])
        res_img_markers = concat_imgs((res_img_markers, resized_img_workspace))
    # Showing images
    show_img("Images raw & undistorted", concat_ims)
    key = show_img("Markers", res_img_markers, wait_ms=30)
    if key in [27, ord("q")]: # Will break loop if the user press Escape or Q
        break

```

## Vision pick via your image processing pipeline

You may want to send coordinate to Ned in order to pick the object of your choice! To do that, use the function **get\_target\_pose\_from\_rel()** ([index.html#api.tcp\\_client.NiryoRobot.get\\_target\\_pose\\_from\\_rel](#)) which converts a relative pose in the workspace into a pose in the robot's world!

```

# Initializing variables
obj_pose = None
try_without_success = 0
count = 0
color_hsv_setting = ColorHSV.ANY.value

mtx, dist = robot.get_camera_intrinsics()
# Loop
while try_without_success < 5:
    # Moving to observation pose
    robot.move_pose(observation_pose)

    img_compressed = robot.get_img_compressed()
    img = uncompress_image(img_compressed)
    img = undistort_image(img, mtx, dist)
    # extracting working area
    im_work = extract_img_workspace(img, workspace_ratio=1.0)
    if im_work is None:
        print("Unable to find markers")
        try_without_success += 1
        if display_stream:
            cv2.imshow("Last image saw", img)
            cv2.waitKey(25)
        continue

# Applying Threshold on ObjectColor
img_thresh = threshold_hsv(im_work, *color_hsv_setting)

if display_stream:
    show_img("Last image saw", img, wait_ms=0)
    show_img("Image thresh", img_thresh, wait_ms=30)
# Getting biggest contour/blob from threshold image
contour = biggest_contour_finder(img_thresh)
if contour is None or len(contour) == 0:
    print("No blob found")
    obj_found = False
else:
    img_thresh_rgb = cv2.cvtColor(img_thresh, cv2.COLOR_GRAY2BGR)
    drawContours(img_thresh_rgb, [contour])
    show_img("Image thresh", img_thresh, wait_ms=30)

    # Getting contour/blob center and angle
cx, cy = get_contour_barycenter(contour)
cx_rel, cy_rel = relative_pos_from_pixels(im_work, cx, cy)
angle = get_contour_angle(contour)

# Getting object world pose from relative pose
obj_pose = robot.get_target_pose_from_rel(workspace_name,
                                             height_offset=0.0,
                                             x_rel=cx_rel, y_rel=cy_rel,
                                             yaw_rel=angle)
obj_found = True

if not obj_found:
    try_without_success += 1
    continue

# Everything is good, so we going to object
robot.pick_from_pose(obj_pose)

# Computing new place pose
offset_x = count % grid_dimension[0] - grid_dimension[0] // 2
offset_y = (count // grid_dimension[1]) % 3 - grid_dimension[1] // 2
offset_z = count // (grid_dimension[0] * grid_dimension[1])
place_pose = center_conditioning_pose.copy_with_offsets(0.05 * offset_x, 0.05 * offset_y, 0.025 * offset_z)

# Placing
robot.place_from_pose(place_pose)

try_without_success = 0
count += 1

```

## Functions documentation

This file presents the different functions and [Enums Image Processing](#) available for image processing.

These functions are divided in subsections:

- [Pure image processing](#) are used to deal with the thresholding, contours detection, ..
- [Workspaces wise](#) section contains functions to extract workspace and deal with the relative position in the workspace
- The section [Show](#) allows to display images
- [Image Editing](#) contains lot of function which can compress images, add text to image, ...

## Pure image processing

### `image_functions.threshold_hsv(list_min_hsv, list_max_hsv, reverse_hue=False, use_s_prime=False)`

Take BGR image (OpenCV imread result) and return thresholded image according to values on HSV (Hue, Saturation, Value) Pixel will worth 1 if a pixel has a value between min\_v and max\_v for all channels

#### Parameters:

- **img** (`numpy.array`) – image BGR if `rgb_space = False`
- **list\_min\_hsv** (`list` (<https://docs.python.org/3/library/stdtypes.html#list>)) [`int` (<https://docs.python.org/3/library/functions.html#int>)] – list corresponding to [min\_value\_H,min\_value\_S,min\_value\_V]
- **list\_max\_hsv** (`list` (<https://docs.python.org/3/library/stdtypes.html#list>)) [`int` (<https://docs.python.org/3/library/functions.html#int>)] – list corresponding to [max\_value\_H,max\_value\_S,max\_value\_V]
- **use\_s\_prime** (`bool` (<https://docs.python.org/3/library/functions.html#bool>)) – True if you want to use S channel as  $S' = S \times V$  else classic
- **reverse\_hue** (`bool` (<https://docs.python.org/3/library/functions.html#bool>)) – Useful for Red color cause it is at both extremum

#### Returns:

threshold image

#### Return type:

`numpy.array`

### `image_functions.debug_threshold_color(color_hsv)`

Return masked image to see the effect of color threshold

#### Parameters:

- **img** (`numpy.array`) – OpenCV image
- **color\_hsv** (`ColorHSV` (<index.html#vision.enums.ColorHSV>)) – Color used for debug

#### Returns:

Masked image

#### Return type:

`numpy.array`

### `image_functions.morphological_transformations(morpho_type=<MorphaType.CLOSE: 3>, kernel_shape=(5, 5), kernel_type=<KernelType.ELLIPSE: 2>)`

Take black & white image and apply morphological transformation

#### Parameters:

- **im\_thresh** (`numpy.array`) – Black & White Image

- **morpho\_type** (*MorphoType*) (index.html#vision.enums.MorphoType)) – CLOSE/OPEN/ERODE/DILATE => See on OpenCV/Google if you do not know these words
- **kernel\_shape** (*tuple*) (<https://docs.python.org/3/library/stdtypes.html#tuple>)*[float]* (<https://docs.python.org/3/library/functions.html#float>) – tuple corresponding to the size of the kernel
- **kernel\_type** (*KernelType*) (index.html#vision.enums.KernelType)) – RECT/ELLIPSE/CROSS => see on OpenCV

**Returns:**

image after processing

**Return type:**

numpy.array

**image\_functions.get\_contour\_barycenter()**

Return barycenter of an OpenCV Contour

**Parameters:**

**contour** (OpenCV Contour) –

**Returns:**

Barycenter X, Barycenter Y

**Return type:**

**int** (<https://docs.python.org/3/library/functions.html#int>), **int** (<https://docs.python.org/3/library/functions.html#int>)

**image\_functions.get\_contour\_angle()**

Return orientation of a contour according to the smallest side in order to be well oriented for gripper

**Parameters:**

**contour** (OpenCV Contour) – contour

**Returns:**

Angle in radians

**Return type:**

**float** (<https://docs.python.org/3/library/functions.html#float>)

**image\_functions.biggest\_contour\_finder()****image\_functions.biggest\_contours\_finder(nb\_contours\_max=3)**

Function to find the biggest contour in an binary image

**Parameters:**

- **img** (numpy.array) – Binary Image
- **nb\_contours\_max** (**int** (<https://docs.python.org/3/library/functions.html#int>)) – maximal number of contours which will be returned

**Returns:**

biggest contours found

**Return type:**

**list** (<https://docs.python.org/3/library/stdtypes.html#list>)[OpenCV Contour]

**image\_functions.draw\_contours(contours)**

Draw a list of contour on an image and return the drawing image

**Parameters:**

- **img** (`numpy.array`) – Image
- **contours** ([list](https://docs.python.org/3/library/stdtypes.html#list) (<https://docs.python.org/3/library/stdtypes.html#list>)) [*OpenCV Contour*] – contours list

**Returns:**

Image with drawing

**Return type:**

`numpy.array`

## Workspaces wise

### `image_functions.extract_img_workspace(workspace_ratio=1.0)`

Extract working area from an image thanks to 4 Niryo's markers

**Parameters:**

- **img** (`numpy.array`) – OpenCV image which contain 4 Niryo's markers
- **workspace\_ratio** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – Ratio between the width and the height of the area represented by the markers

**Returns:**

extracted and warped working area image

**Return type:**

`numpy.array`

### `image_functions.debug_markers(workspace_ratio=1.0)`

Display detected markers on an image

**Parameters:**

- **img** (`numpy.array`) – OpenCV image which contain Niryo's markers
- **workspace\_ratio** ([float](https://docs.python.org/3/library/functions.html#float) (<https://docs.python.org/3/library/functions.html#float>)) – Ratio between the width and the height of the area represented by the markers

**Returns:**

(status, annotated image)

**Return type:**

`numpy.array`

### `image_functions.relative_pos_from_pixels(x_pixels, y_pixels)`

Transform a pixels position to a relative position

**Parameters:**

- **img** (`numpy.array`) – Image where the object is detected
- **x\_pixels** ([int](https://docs.python.org/3/library/functions.html#int) (<https://docs.python.org/3/library/functions.html#int>)) – coordinate X
- **y\_pixels** ([int](https://docs.python.org/3/library/functions.html#int) (<https://docs.python.org/3/library/functions.html#int>)) – coordinate Y

**Returns:**

X relative, Y relative

**Return type:**

**float** (<https://docs.python.org/3/library/functions.html#float>), **float** (<https://docs.python.org/3/library/functions.html#float>)

## Show

### `image_functions.show_img_and_check_close(img)`

Display an image and check whether the user want to close

**Parameters:**

- **window\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – window's name
- **img** ([numpy.array](https://numpy.org/doc/stable/reference/generated/numpy.array.html)) – Image

**Returns:**

boolean indicating if the user wanted to leave

**Return type:**

[bool](https://docs.python.org/3/library/functions.html#bool) (<https://docs.python.org/3/library/functions.html#bool>)

### `image_functions.show_img(img, wait_ms=1)`

Display an image during a certain time

**Parameters:**

- **window\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – window's name
- **img** ([numpy.array](https://numpy.org/doc/stable/reference/generated/numpy.array.html)) – Image
- **wait\_ms** ([int](https://docs.python.org/3/library/functions.html#int)) – Wait time in milliseconds

**Returns:**

value of the key pressed during the display

**Return type:**

[int](https://docs.python.org/3/library/functions.html#int) (<https://docs.python.org/3/library/functions.html#int>)

### `image_functions.show_img_and_wait_close(img)`

Display an image and wait that the user close it

**Parameters:**

- **window\_name** ([str](https://docs.python.org/3/library/stdtypes.html#str)) – window's name
- **img** ([numpy.array](https://numpy.org/doc/stable/reference/generated/numpy.array.html)) – Image

**Returns:**

None

**Return type:**

[None](https://docs.python.org/3/library/constants.html#None) (<https://docs.python.org/3/library/constants.html#None>)

## Image Editing

### `image_functions.compress_image(quality=90)`

Compress OpenCV image

**Parameters:**

- **img** ([numpy.array](https://numpy.org/doc/stable/reference/generated/numpy.array.html)) – OpenCV Image
- **quality** ([int](https://docs.python.org/3/library/functions.html#int)) – integer between 1 - 100. The higher it is, the less information will be lost, but the heavier the compressed image will be

**Returns:**

status & string representing compressed image

**Return type:**

`bool` (<https://docs.python.org/3/library/functions.html#bool>), `str` (<https://docs.python.org/3/library/stdtypes.html#str>)

**image\_functions.uncompress\_image()**

Take a compressed img and return an OpenCV image

**Parameters:**

`compressed_image` (`str` (<https://docs.python.org/3/library/stdtypes.html#str>)) – compressed image

**Returns:**

OpenCV image

**Return type:**

`numpy.array`

**image\_functions.add\_annotation\_to\_image(text, write\_on\_top=True)**

Add Annotation to an image

**Parameters:**

- `img` (`numpy.array`) – Image
- `text` (`str` (<https://docs.python.org/3/library/stdtypes.html#str>)) – text string
- `write_on_top` (`bool` (<https://docs.python.org/3/library/functions.html#bool>)) – if you write the text on top

**Returns:**

img with text written on it

**Return type:**

`numpy.array`

**image\_functions.undistort\_image(mtx, dist)**

Use camera intrinsics to undistort raw image

**Parameters:**

- `img` (`numpy.array`) – Raw Image
- `mtx` (`list` (<https://docs.python.org/3/library/stdtypes.html#list>)) [`float` (<https://docs.python.org/3/library/stdtypes.html#list>)] – Camera Intrinsics matrix
- `dist` (`list` (<https://docs.python.org/3/library/stdtypes.html#list>)) [`float` (<https://docs.python.org/3/library/stdtypes.html#list>)] – Distortion Coefficient

**Returns:**

Undistorted image

**Return type:**

`numpy.array`

**image\_functions.resize\_img(width=None, height=None, inter=3)**

Resize an image. The user should precise only width or height if he wants to keep image's ratio

**Parameters:**

- **img** (`numpy.array`) – OpenCV Image
- **width** ([int](https://docs.python.org/3/library/functions.html#int)) – Target Width
- **height** ([int](https://docs.python.org/3/library/functions.html#int)) – Target Height
- **inter** (`long`) – OpenCV interpolation flag

**Returns:**

resized image

**Return type:**

`numpy.array`

**image\_functions.concat\_imgs(`axis=1`)**

Concat multiple images along 1 axis

**Parameters:**

- **tuple\_imgs** ([tuple](https://docs.python.org/3/library/stdtypes.html#tuple)) [`numpy.array`] – tuple of images
- **axis** ([int](https://docs.python.org/3/library/functions.html#int)) – 0 means vertically and 1 means horizontally

**Returns:**

Concat image

**Return type:**

`numpy.array`

## Enums Image Processing

Enums are used to pass specific parameters to functions.

List of enums:

- **ColorHSV**
- **ColorHSVPrime**
- **ObjectType**
- **MorphoType**
- **KernelType**

**class ColorHSV(`value`)**

MIN HSV, MAX HSV, Invert Hue (bool)

**BLUE= ([90, 50, 85], [125, 255, 255], False)**

**RED= ([15, 80, 75], [170, 255, 255], True)**

**GREEN= ([40, 60, 75], [85, 255, 255], False)**

**ANY= ([0, 50, 100], [179, 255, 255], False)**

**class ColorHSVPrime(`value`)**

MIN HSV, MAX HSV, Invert Hue (bool)

**BLUE= ([90, 70, 100], [115, 255, 255], False)****RED= ([15, 70, 100], [170, 255, 255], True)****GREEN= ([40, 70, 100], [85, 255, 255], False)****ANY= ([0, 70, 140], [179, 255, 255], False)*****class ObjectType(value)***

An enumeration.

**SQUARE= 4****TRIANGLE= 3****CIRCLE= -1****ANY= 0*****class MorphoType(value)***

An enumeration.

**ERODE= 0****DILATE= 1****OPEN= 2****CLOSE= 3*****class KernelType(value)***

An enumeration.

**RECT= 0****ELLIPSE= 2****CROSS= 1**

## Indices and tables

- [Index](#) (genindex.html)
- [Module Index](#) (py-modindex.html)
- [Search Page](#) (search.html)

[Suggest a modification](#)

[Download as PDF](#)